

E-spares cat Now Available Online  
at [www.invacare.co.uk](http://www.invacare.co.uk)

Customer Name :	Account No :
Dealer :	Order No :
Address :	Order Qty :
Tel No:	Order Date :
Delivery address :	



**Standard Features**



- Max user weight 180kg (28 stone)
- Compact, manoeuvrable, **CWD TDX SP2**
- Modulite seating system
- Tension Adjustable backrest
- Modulite Telescopic adjustable Seat Widths from 380mm - 530mm (15" - 21") additional 50mm, (2") by arms
- Modulite Seat Depths 410mm - 530mm Stepless Adjustment (16" - 20")
- Modulite Detachable T-Armrests, height, width and depth adjustable with standard armpad
- Manually settable seat tilt and backrest recline
- Standard 80° Hangers with fixed footplate
- ACS 2 Control System with REM A Joystick
- Centre (Driving Wheels) 3.00 by 8.00 Grey Puncture Proof (Solid Inlay)
- Choice of 7 modern colours for Chassis
- Suspension and Sure Step Technology
- 6km/h 4 Pole Motors and 73Ah Batteries
- 4 Point tie down system
- Lapbelt with autolock (adjustable both sides)
- Wheel lock / parking brake
- 8A (200 - 250V) Charger

**Standard Base Price £4,800**

Invacare® **Modulite™**  
An evolution in seating designed around you.



\* Configurations may vary from the image shown

Please make your selection with a cross in the box ◆ Standard feature ○ Free of Charge NC = Incompatibility with option XXX99

Code	Description	DEALER
ARJ0000	<b>TDX SP2 Basic Model</b>	◆
ARZ0000	TDX SP2 Special Model - PLEASE COMPLETE DETAILS AT THE END OF THE FORM	POA
<b>SEATING TYPE : MODULITE MULTI ADJUSTABLE SEATING SYSTEM</b>		
<b>(FOR 'MAX' SEATING PLEASE MOVE FURTHER DOWN THE FORM)</b>		
ARJ2100	Modulite Seating Unit	◆
<b>SEAT FRAME AND SEAT PLATE</b>		
ARJ2111	Telescopic Seat Frame with metal seat plate (with stepless adjustment on armrests, legrests and backrest)	£4,800

Customer Ref:

Order Number:

MODULITE SEAT SIZE SETTING CHART:

			SEAT WIDTHS				
			330mm* (13")	ARJ2121 380mm (15")	ARJ2122 430mm (17")	ARJ2123 480mm (19")	ARJ2124 530mm (21")
Seat Depth ARJ2129	Standard Seat Rail 410 - 510mm	ARJ2132 (16") 410mm	○	○	○	○	○
		ARJ2134 (18") 460mm	○	○	○	○	○
		ARJ2136 (20") 510mm	○	○	○	○	○
Seat Depth ARJ2130	Extended Seat Rail 510 - 560mm	ARJ2136 (20") 510mm	○	○	○	○	○
		ARJ2138 (22") 560mm	○	○	○	○	○

Invacare® **Modulite™**  
An evolution in seating designed around you.



\*330mm (13") seat width achieved on a 380mm seat frame (2121) with inrigged armrests (2183)

1. Seat depth 410mm: cannot combine T-Armrests (2191) with flexible lapbelt mount (2113)
2. No seat cushion available with 560mm seat depth (2138)
3. Seat width 530mm (2124) not possible in combination with following armrests (2192/4)
4. Seat width 560mm (2138) When in combination with Flex Tension Adjustable Backrest the seat depth is reduced by 10mm

SEAT CUSHION (sizes according to seat frame settings)

ARJ2142	Matrx Flo-tech Lite - pressure relieving seat cushion	£92
ARJ2143	Matrx Flo-tech Contour Visco 430 x430 mm	£92
ARJ2147	Matrx Libra Cushion, 90mm thick	Not available with Seat Width 530mm with 410mm seat Depth £319
ARJ2140	Without seat cushion	○

BACKREST SUPPORT

ARJ2162	Backposts 540mm (tall) must be selected for Flex3 adjustable backrest plates	◆
ARJ2163	Backposts 480mm (short)	○

SPECIAL BACKPOST MODIFICATION

SWO0001	<b>No push handles</b> (extended leadtime). Push handles will be cut down on the vertical section. If you require a specific length of push handle please note in the specials box at the end of the form.	Not available with attendant controls, dual controls, tension adjustable backrest in combination with a headrest.	£67
---------	--	---	-----

BACKREST UPHOLSTERY / CUSHION

ARJ2171	Tension Adjustable Backrest (540mm or 480mm upholstery height) Straps easily accessible and adjustable from the rear	Not available with laterals, black nylon / fabric upholstery only	◆
ARJ2179	<b>Flex3 Tension Adjustable (long = 540mm / short = 480mm), available in Black Dartex only</b>	<b>Not available with Flip up or Following armrests, Chin Remote, Laterals available though specials</b>	£113
ARJ2172	Flex3 adjustable Backrest Plates (without Backrest Cushion)		£183
ARJ2173	Flex3 adjustable Backrest Plates with Laguna Backrest cushion		£297
ARJ2174	Matrx Posture Back, height 400mm (Mesh fabric)	Not available with Laterals or 530mm (21") seat width, non Invacare Headrest but profile mount ARJ1496	£406
ARJ2175	Matrx Posture Back, height 500mm (Mesh fabric)	Not available with Laterals or 530mm (21") seat width, non Invacare Headrest but profile mount ARJ1496	£406
ARJ2176	Matrx Elite Back, height 500mm (Mesh fabric)	Not available with Laterals or 530mm (21") seat width, non Invacare Headrest but profile mount ARJ1496, seat width of 380mm (15") ARJ2121 in combination with flip up or following armrests	£406
ARJ2177	Matrx Elite TR Back, height 500mm (Mesh fabric)	Not available with Laterals or 530mm (21") seat width, non Invacare Headrest but profile mount ARJ1496	£406
ARJ2178	Matrx Elite Deep Back, height 500mm (Mesh fabric)	Not available with Laterals or non Invacare Headrest but profile mount ARJ1496	£450
ARJ2170	Without Backrest Upholstery (Straight 1" diameter Back Posts/Back Frame suitable for mounting specialist backrests)	Not available with laterals or non Invacare Headrest but profile mount ARJ1496	○

Customer Ref:

Order Number:

UPHOLSTERY COLOUR			
ARJ2151		Black fabric (2-way-stretch material w 2142, 2143)	◆
ARJ2152		Black Dartex	○
ARJ2153		LK Version (black with orange seam)	○
ARJ2154		Black Plush	○
ARMRESTS (NOTE: HIP GUIDES ONLY POSSIBLE WITH FLIP UP ARMRESTS)			
ARJ2191		<b>Single post 'T-Armrest'</b> removable and height adjustable (with stepless positioning relative to the backposts) Range of adjustment from seat plate <b>245mm to 315mm</b>	◆
SWO0002		<b>Extended height Single post 'T-Armrest'</b> : Range of adjustment from seat plate <b>295mm to 365mm</b> (extended leadtime)	£115
ARJ2183		<b>In-rigged Armrests</b> (reduces seat width by 25 mm each side)	£20
ARJ2193		<b>Flip-up Armrest</b> (height and depth adjustable) (height adjustable from <b>230 - 360mm</b> ) <b>Max user weight 150kg</b>	Not available with laterals £47
ARJ2192		<b>Following Armrest</b> (height adjustable from 230 - 300 mm) <b>Max user weight 150kg</b>	£135
ARJ2194		Following Armrest (height adjustable from 300 - 360 mm) <b>Max user weight 150kg</b>	£135
ARMPADS			
ARJ2195		Short Armpad (325 mm)	◆
ARJ2196		Long Armpad (375 mm)	Not available with seat depth 410mm (2132) ○
ARJ2197		Long Ergonomic Armpad (397 mm)	Not available with seat depth 410mm (2132) ○
ARJ2199		Soft Armpad (REA Azalea Wide Pad)	Not available with seat depth 410mm (2132) £70
HIP SUPPORTS LH (for flip up and following armrests only)			
ARJ2184		Multi adjustable hip support LH, pad size SMALL 140 x 100mm	Not available with 410mm seat depth and manual centre legrests £185
ARJ2186		Multi adjustable hip support LH, pad size MEDIUM 230 x 100mm	Not available with 410mm seat depth and manual centre legrests £185
ARJ2188		Multi adjustable hip support LH, pad size LARGE 330 x 100mm	Not available with 410mm seat depth and manual centre legrests £185
HIP SUPPORTS RH (for flip up and following armrests only)			
ARJ2185		Multi adjustable hip support RH, pad size SMALL 140 x 100mm	Not available with 410mm seat depth and manual centre legrests £185
ARJ2187		Multi adjustable hip support RH, pad size MEDIUM 230 x 100mm	Not available with 410mm seat depth and manual centre legrests £185
ARJ2189		Multi adjustable hip support RH, pad size LARGE 330 x 100mm	Not available with 410mm seat depth and manual centre legrests £185
SEAT ACCESSORIES			
HEADREST/NECKREST (Dartex)			
ARJ1350		Neckrest: height and depth adjustable (Azalea type, dartex covered)	£59
ARJ1432		Headrest: height, depth and sideways adjustable (Azalea type, dartex covered)	◆
ARJ1427		Headrest with cheek support: height and depth adjustable (Azalea type, dartex covered)	£166
ARJ1496		All bracketry as option 1432 without headrest cushion	○
ARJ1495		Push bar dia. 22mm in conjunction with telescopic seat frame/backposts, No Headrest - no mounting hardware.	○








Customer Ref:

Order Number:

LATERAL SUPPORT						
ARJ1446		Lateral Support (pair) (Dartex) (Azalea type)	£123			
ARJ1445		Lateral Support with extended arm (pair) (Dartex), width and height adjustable (Azalea type)	£123			
SEAT POSITIONING ADJUSTMENTS						
SEAT LIFTER / SEAT TILT						
		Tilt Range	SEAT HEIGHT			
			ARJ0354 450mm (17.75")	460mm (18")	ARJ0364 475 (18.75")	
ARJ1091	Manually settable	0° - 20°				◆
ARJ1092	Powered basic (max user weight 136kg)	0° - 20°				£319
ARJ1095	Powered basic (max user weight 180kg)	0° - 20°				£319
ARJ1093	Powered centre of gravity shift (max user weight 160kg)	0° - 30°				£824
ARJ1094	Seat Lifter with centre of gravity shift (300mm rise) (max user weight 160kg)	0° - 30°				£1,028
Only available with REM 24SD or G90A controls						
SEAT TILT OPTION (only if option 1093 or 1094 have been selected)						
ARJ2054		5° Fixed Anterior Tilt Assist mount (will reduce max tilt angle by 5°)	£295			
BACKREST RECLINER						
ARJ0975		Manual Settable Recliner (0° - 30°) removable with handwheels at the base of the backpost for easy folding	◆			
ARJ0977		Powered Recliner (0° - 30°) with quick release for easy folding of backrest	£245			
SEATING TYPE : 'MAX' ADJUSTABLE SEATING SYSTEM						
ARH1145		MAX Seating system	No tray available with Max seating	£5,468		
SEAT WIDTH SETTING (3 x SEAT FRAMES)						
ARH0348		Seat Width of Seat Frame 490mm (+100mm adjustment in arms)	○			
ARH0378		Seat Width of Seat Frame 550mm (+100mm adjustment in arms)	○			
ARH0380		Seat Width of Seat Frame 610mm (+100mm adjustment in arms)	Not available with Vari F legrest	○		
SEAT DEPTH SETTING (SETTABLE IN 3 STEPS)						
ARH0396		Seat Depth 560mm	○			
ARH0397		Seat Depth 595mm	○			
ARH0398		Seat Depth 630mm	○			
SEAT CUSHION						
ARH1484		MAX Seat Cushion	○			
ARH1114		Without Seat Cushion (10mm seat board only)	○			
ADAPTOR TO REDUCE THE SEAT DEPTH BY 130mm (REQUIRED IF SEAT WIDTH 490mm IS SELECTED IN COMBINATION WITH SEAT DEPTH 595mm OR 630mm)						
ARH0399		Adaptor to shorten the seat depth by 130mm	○			
BACKREST UPHOLSTERY						
ARH1106		Tension Adjustable	◆			
UPHOLSTERY COLOUR						
ARH0278		Black Dartex	◆			
ARMREST (INCLUDING REINFORCED, DEPTH ADJUSTABLE BRACKET)						
ARH0635		Fixed Armrest, 250 - 340mm tool adjustable in height, detachable	◆			
ARH0636		Fixed Armrest, 290 - 380mm tool adjustable in height, detachable	○			
ARMPADS						
ARH0647		Short armpad (325mm)	◆			
ARH0649		Long armpad (375mm)	○			
PUSH BAR						
ARH1307		Push bar	◆			
HEADREST						
ARH1432		Headrest: height, depth and sideways adjustable (Azalea type)	◆			
ARH1427		Headrest with cheek support: height and depth adjustable (Azalea type)	£166			
ARH1495		Non Invacare Headrest - no mounting hardware (with pushbar dia. 22mm in conjunction with telescopic seat frame/backposts)	○			

Customer Ref:







Order Number:

MAX SEAT ADJUSTMENTS			
<b>SEAT TILT</b>			
ARH1091		Manually settable Seat Tilt (Max user weight 180kg)	◆
ARH1095		Electric heavy duty basic tilt (0° - 20°) select either 470mm or 495mm seat to ground height. <b>(Max user weight 180kg)</b>	£319
<b>SEAT HEIGHT SETTING</b>			
ARH1024		Fixed seat height 470mm (18.5")	○
ARH1046		Fixed seat height 495mm (19.5")	○
<b>BACKREST RECLINER</b>			
ARH0958		Manually settable recliner by spindle (0° - 25°)	◆
ARH0947		Electrical recliner (0° - 25°)	£245
<b>PLEASE CHOOSE EITHER CENTRE MOUNTED OR SIDE MOUNTED LEGRESTS</b>			
<b>CENTRE MOUNTED LEGRESTS</b>			
<b>FOOTREST HANGER TYPE (INCLUDING CALF PAD AND ANGLE ADJ FOOTPLATE)</b>			
ARJ0556		<b>Manual</b> : Hanger, centre mounted, angle settable (0° - 90°): <b>knee to heel length 390 - 510mm</b>	£233
ARJ0557		<b>Powered</b> : Centre mounted, length compensating (0° - 90°): <b>knee to heel length 390 - 510mm</b>	£1,300
<b>CALF PADS</b>			
ARJ0583		Calf pads	○
<b>FOOTPLATE TYPE</b>			
ARJ0559		Angle Adjustable Footplate	○
<b>SIDE MOUNTED LEGRESTS</b>			
<b>FOOTREST HANGER TYPE; PAIRS</b>			
ARH0564		<b>Reinforced hangers</b> , fixed angle (not detachable without tools) incl Calf pads	For MAX seating only ○
ARJ0526		<b>Standard hanger 80° pair</b>	<b>Not available with angle adjustable footplate</b> ◆
ARJ0503		<b>VARI F</b> : Angle adjustable legrest hanger (swing away and detachable) (adjustable via tools 0° - 70°) : knee to heel length 290 - 460mm	○
ARJ0533		<b>VARI A</b> : Manually angle adjustable legrest hanger (swing away and detachable) (adjustable via hand wheel / knob 0° - 70°) - complete with calf pad support : knee to heel length 290 - 460mm	Not available with fixed footplates £92
ARJ0532		<b>ADM Manual</b> : Manually angle adjustable (0° - 80°) legrest hanger (swing away and detachable) (adjustable via quick release lever) - complete with calf pad support : knee to heel length 290 - 460mm	Not available with fixed footplates £163
ARJ0531		<b>ADE Powered</b> : Electrically angle adjustable (0° - 80°) legrest hanger (swing away and detachable) (powered via joystick) - complete with calf pad support : knee to heel length 290 - 460mm	Not available with fixed footplates £478
<div style="display: flex; justify-content: space-between; align-items: flex-start;"> <div style="text-align: center;"> <p><b>Standard 80° 0526</b></p>  </div> <div style="text-align: center;"> <p><b>VARI F</b></p>  </div> <div style="text-align: center;"> <p><b>VARI A</b></p>  </div> <div style="text-align: center;"> <p><b>ADM</b></p>  </div> <div style="text-align: center;"> <p><b>ADE</b></p>  </div> <div style="text-align: center;"> <p><b>MANUAL CENTRE</b></p>  </div> <div style="text-align: center;"> <p><b>POWERED CENTRE</b></p>  </div> </div>			
<b>AMPUTEE SUPPORT</b>			
ARJ0555		<b>Amputee Support LH</b>	£143
ARJ0554		<b>Amputee Support RH</b>	£143
<b>FOOTPLATE TYPE; PAIRS</b>			
ARJ0512		Fixed Footplate	○
ARJ0504		Angle Adjustable Footplate	£54
ARJ0513		Angle and depth adjustable Footplate (Adapt)	£87
<b>FOOTBOARDS for Max Seat</b>			
ARH0599		Angle adjustable flip-back, heavy-duty footboard (Recommended for SW ≥ 550mm)	Only with Reinforced hangers £245
ARH0522		Angle Adjustable Footboard (for max user weight up to 180 kg) (Recommended for SW ≥ 550 mm)	Only with Vari F hangers £165

Customer Ref:


Order Number:

## ACS2 CONTROLS

USER CONTROL + BRACKET			
ARJ1638	REM Type A (Up to 4 electric actuators)	Not in combination with lights	◆
ARJ1639	REM Type B (up to 4 electric actuators and light functionality)		£114
ARJ1601	ACS REM 24SD		£147
ARJ1614	G90A		£147
	Head-Remote, ultrasonic proximity sensor (only available as special) Please contact UK Specials Dept		POA
ARJ1207	Standard Bracket RH		○
ARJ1208	Standard Bracket LH		○
ARJ1707	Swing Away + height adjust. RH High End		£105
ARJ1708	Swing Away + height adjust. LH High End		£105
ATTENDANT REMOTE + BRACKET			
ARJ1643	REM Type A (Up to 4 electric actuators)	Not in combination with lights	£14
ARJ1644	REM Type B (up to 4 electric actuators and light functionality)		£125
ARJ1602	ACS REM 24SD		£158
ARJ1616	G90A		£158
ARJ1207	Standard Bracket RH		○
ARJ1208	Standard Bracket LH		○
TRAY & CENTRE CONTROLS			
ARJ1608	A C S Centre Control (REM 24SD)		£168
ARJ1607	A C S Tray Control c/w Acrylic Tray & Actuator Functionality - DX REM 41		£268
ARJ1207	Standard Bracket RH		○
ARJ1208	Standard Bracket LH		○
CHIN CONTROL SIMPLE			
ARJ1719	Chin Remote with standard Joystick (Always to be mounted on the Master Remote's side) (for Modulite only)	(Only compatible with Option1601 )	£1,632
ARJ1289	Swing Away Bracket RH, manual	Only available with RH user bracket	○
ARJ1288	Swing Away Bracket LH, manual	Only available with LH user bracket	○
ARJ1287	Swing Away Bracket RH, powered	Only available with RH user bracket	£1,100
ARJ1286	Swing Away Bracket LH, powered	Only available with LH user bracket	£1,100
DUAL CONTROL : CAN ONLY BE SELECTED IN ADDITION TO USER CONTROL			
ARJ1720	Compact Dual Control		£377
ARJ1207	Standard Bracket RH	Only available with LH user bracket	○
ARJ1208	Standard Bracket LH	Only available with RH user bracket	○
INVACARE IDC DUAL CONTROL + BRACKET			
ARJ1645	Invacare IDC Dual control		£542
ARJ1221	Detachable Bracket RH		○
ARJ1222	Detachable Bracket LH		○
<div style="display: flex; justify-content: space-around; text-align: center;"> <div> <p><b>REM Type A</b></p>  </div> <div> <p><b>REM 24SD</b></p>  </div> <div> <p><b>G90 A</b></p>  </div> <div> <p><b>Chin control</b></p>  </div> <div> <p><b>Dual control</b></p>  </div> <div> <p><b>IDC control</b></p>  </div> </div>			
DRIVE UNIT OPTIONS			
FRAME COLOURS (ALL COLOURED PARTS ARE DETACHABLE AND CHANGEABLE PLASTIC SHROUDS)			
ARJ0228	Wet Black		○
ARJ0212	Electric Red		○
ARJ0257	Jade Green		○
ARJ0236	Electric Blue		○
ARJ0237	Electric Silver		○
ARJ0238	Arctic White		○
ARJ0205	Transparent Purple		○
CENTRE WHEELS (DRIVING WHEELS)			
ARJ0771	3.00 by 8.00 Grey Pneumatic		○
ARJ0773	3.00 by 8.00 Grey Puncture Protected (sealed)		○
ARJ0750	3.00 by 8.00 Grey Puncture Proof (Solid Inlay)		◆

Customer Ref:

Order Number:

PERFORMANCE OPTIONS			
<b>MOTORTYPE</b>			
ARJ1282		6 km/h (4mph) (4 pole)	◆
ARJ1284		8 km/h (5mph) (4 pole)	○
BATTERY			
ARJ1260	X	73 Ah	◆
PERFORMANCE ACCESSORIES			
<b>G-TRAC</b>			
ARJ1215		 <b>G-trac</b> ensures the powerchair takes the intended direction as indicated by the user control, without the need for continual correction by the client on uneven or sloping terrain.	£605
<b>CHARGER</b>			
ARJ1343	X	Charger 8A (200 - 250V)	◆
SAFETY OPTIONS			
<b>LIGHT</b>			
ARJ1437		Light system with Indicators & LED Headlights	Only available with REM24SD or G90A joystick controls
			£325
ARJ1254		Without light system	◆
SAFETY ACCESSORIES			
<b>TRANSPORTATION KIT</b>			
ARJ1450		<b>Preparation</b> for DAHL Docking station: Docking station for car implementation, please order docking station separately (Adequate headrest is compulsory, see options below)	£430
ARJ1494		4-Point tie-down (Tested version for using wheelchair as a vehicle seat according to ISO 7176-19:2001) (Adequate headrest is compulsory, see options below)	◆
<b>HEADRESTS ADEQUATE FOR PERSON TRANSPORTATION</b>			
		INVACARE Headrest (Please select a headrest from this form)	
ARJ1495		Equivalent Headrest supporting the requirements of ISO7176-19:2001, provided independently from this order	○
<b>FIXATION BELT</b>			
ARJ1329		Lapbelt (adjustable both sides) fitted to base of backrest (MAX lapbelt when MAX seating is selected)	○
ARJ2113		Lapbelt (adjustable both sides) adjustable mounting on seat rail extrusion	Not available in combination with 410mm seat depth (2132) and T-armrest (2191)
			○
ARJ1332		H-Belt	£54
<b>REARVIEW MIRROR</b>			
ARJ1415		RH	○
ARJ1339		LH	○
ACCESSORIES			
<b>CANEHOLDER</b>			
ARJ1312		RH	£18
ARJ1363		LH	£18
<b>TRAY - mounted opposite side to controller</b>			
ARJ1404		Swing Away RH	£249
ARJ1405		Swing Away LH	£249
<b>KLICKfix BAG AND BRACKETS</b>			
ARJ1480		Handy Zip Bag, Removable	£51
ARJ1478		Quick Release Bracket RH	○
ARJ1479		Quick Release Bracket LH	○
			<b>TOTAL PRICE</b>

Customer Ref:

Order Number:

**Special Instructions - Please supply the following information, and your preferred acknowledgement method:**

**Name:**

**Email:**

**Fax:**

**Description of special request (please include relevant option code or part number):**

**Customer Ref:**

**Order Number:**