



Tel: 01656 776222 Fax: 01656 776220
Email: ordersuk@invacare.com
Online spares available at: www.invacare.co.uk

July 2016

Customer Ref:	Account code:	Built in Switzerland
Contact No:	Order no:	
Delivery address :	Order qty:	
	Order date:	

Küschall® Compact
RETAIL
Maximum user weight 130kg (20.5 Stone)

Standard Features



- Active, Lightweight, Folding Aluminium Frame
- Multiple Adjustments including Front and Rear Seat Heights and Centre of Gravity
- Swingaway Legrests, choice of 70°, 80° or 90° Angles with Flip up Plastic Footrests
- Choice of 12 Seat Widths 280mm - 500mm (11" - 20")
- Choice of 10 Seat Depths 320mm - 500mm (12" - 20")
- Fixed Backrest Tubes set at 90° with Standard Light Nylon Upholstery
- Standard Fixed Long Push Handles
- Fixed Carbon Clothesguards
- 5" Comfort Front Castors with Black Tyres
- 24" Self Propel Q/R Wheels with Grey Aluminium Handrims & Schwalbe Black & Grey Pneumatic Tyres
- Pull to Lock Brake
- Choice of 15 Frame Colours



Standard Base Price: £1,350

* Configurations may vary from the image shown

<input checked="" type="checkbox"/>	Please make your selection with a cross in the box	<input type="checkbox"/>	Standard feature	<input type="checkbox"/>	Free of Charge
Customer weight =		Kg			
Code	<input checked="" type="checkbox"/>	Description	Retail		
DDP0012		Küschall® Compact	£1,350		
DDZ0000		Küschall® Compact Special Basic Model - PLEASE COMPLETE DETAILS AT THE END OF THE FORM	P.O.A.		

SEATING CHOICE

FRONT FRAME

Important note: In certain configurations, the Push to lock brakes can conflict with the legrest hanger

DDP0211	Standard		<input checked="" type="checkbox"/>
DDP0224	Abducted, Left	(NC: 1300,1302)	<input type="checkbox"/>
DDP0205	Abducted, Right	(NC: 1300,1302)	<input type="checkbox"/>
DDP0219	Abducted Left and Right	(NC: 1300,1302)	<input type="checkbox"/>

SEAT WIDTH (SB) & SEAT DEPTH (ST)

		Seat Width in mm											
		280mm DDP0118	300mm DDP0120	320mm DDP0122	340mm DDP0124	360mm DDP0126	380mm DDP0128	400mm DDP0130	420mm DDP0132	440mm DDP0134	460mm DDP0136	480mm DDP0138	
Seat Depth in mm	320mm DDP0408	<input type="checkbox"/>	Seat widths 280mm - 420mm are fitted with a single crossbar. Max user weight 100kg/16 Stone										
	340mm DDP0410	<input type="checkbox"/>	Seat widths 440mm - 500mm are fitted with a double crossbar. Max user weight 130kg/20.5 Stone										
	360mm DDP0412	<input type="checkbox"/>	Light back upholstery, option 0721 is not available on seat widths 280mm-320mm or 500mm or on backrest heights less than 360mm										
	380mm DDP0415	<input type="checkbox"/>	One Arm Drive options 1118 & 1119 are only available on seat widths 380mm - 500mm										
	400mm DDP0417	<input type="checkbox"/>	Seat widths 280mm - 320mm are only compatible with 1 or 2 piece carbon footplates. DDP1578,1579,1566 & 1567										
	420mm DDP0420	<input type="checkbox"/>											
	440mm DDP0422	<input type="checkbox"/>											
	460mm DDP0424	<input type="checkbox"/>											
	480mm DDP0427	<input type="checkbox"/>											
	500mm DDP0429	<input type="checkbox"/>											

Matrix backrests are not compatible with seat widths 280mm - 340mm

Matrix Libra cushion is compatible with seat sizes 360-500mm wide and 400-500mm deep

Seat width 360mm on Matrix Libra & Vi cushions is only available with 400mm & 420mm seat depths

Matrix MX1 backrest is not compatible with seat widths 480mm & 500mm

Matrix Vi cushion is compatible with seat sizes 360-500mm wide and 400-500mm deep

BACK HEIGHT (RH)

Backrest Height in mm							
300mm DDP0550	○	315mm DDP0503	○	330mm DDP0506	○	345mm DDP0509	○
360mm DDP0512	○	375mm DDP0514	○	390mm DDP0516	○	405mm DDP0518	○
420mm DDP0520	○	435mm DDP0522	○	450mm DDP0524	○	465mm DDP0526	○
480mm DDP0528	○	495mm DDP0530	○	510mm DDP0532	○		

Back heights 300mm - 330mm are not compatible with folding backrest, option DDP0734

Back heights 300mm - 345mm are not compatible with light back upholstery, option DDP0721

Back height 345mm is not compatible with folding backrest, option DDU0734 in combination with backrest angles 94° & 98°

Backbrace, option DDP1491 is not compatible on back heights 300mm to 390mm

If fitting a Matrix or fixed backrest, please select option DDP9999, backposts to suit Matrix, no back height selection required

Please note this chair cannot be used as a seat during transportation if the back height is lower than 405mm

DDP9999

Backposts to suit Matrix or fixed backrests. Supplied without upholstery and backposts are 255mm high. Only compatible with fixed backrest DDP0736

○

FRONT SEAT-TO-FLOOR HEIGHT (SHv)

Front Seat to Floor in mm				
410mm DDP0843	420mm DDP0844	430mm DDP0845	440mm DDP0846	450mm DDP0847
○	○	○	○	○
460mm DDP0848	470mm DDP0849	480mm DDP0850	490mm DDP0851	500mm DDP0852
○	○	○	○	○
510mm DDP0853	520mm DDP0854	530mm DDP0855		
○	○	○		

REAR SEAT-TO-FLOOR HEIGHT (SHh)

Rear Seat to Floor in mm				
380mm DDP0943	390mm DDP0944	400mm DDP0945	410mm DDP0946	420mm DDP0947
○	○	○	○	○
430mm DDP0948	440mm DDP0949	450mm DDP0950	460mm DDP0951	470mm DDP0952
○	○	○	○	○
480mm DDP0953	490mm DDP0954	500mm DDP0955		
○	○	○		

On 70° frame rear seat height can be from 0cm upto 10cm lower than front seat height.

On 80° & 90° frames rear seat height can be from 3cm upto 10cm lower than front seat height

Mudguards DDP1405 in combination with rear wheel extension DDP0249 are not compatible with rear seat heights of 410mm or lower

Rear wheel size/Rear seat height

22" : 380-470 mm

24" : 390-490 mm

26" : 420-500 mm

KNEE-TO-HEEL LENGTH (UL)

Knee to Heel in mm									
210mm DDP1662	○	220mm DDP1663	○	230mm DDP1664	○	240mm DDP1665	○	250mm DDP1666	○
260mm DDP1667	○	270mm DDP1668	○	280mm DDP1669	○	290mm DDP1670	○	300mm DDP1671	○
310mm DDP1672	○	320mm DDP1673	○	330mm DDP1674	○	340mm DDP1675	○	350mm DDP1676	○
360mm DDP1677	○	370mm DDP1678	○	380mm DDP1679	○	390mm DDP1680	○	400mm DDP1681	○
410mm DDP1682	○	420mm DDP1683	○	430mm DDP1684	○	440mm DDP1685	○	450mm DDP1686	○
460mm DDP1687	○	470mm DDP1688	○	480mm DDP1689	○	490mm DDP1690	○	500mm DDP1691	○
510mm DDP1692	○								

On 70° legrests, Maximum knee to heel length must be front seat height minus 20mm

On 80° & 90° legrests, Maximum knee to heel length must be front seat height minus 30mm

Knee to heel lengths 210mm - 260mm are available with 1 or 2 piece carbon footplates options DDP1578, DDP1579, DDP1566 & DDP1567

Knee to heel lengths 270mm - 330mm are available with 1 or 2 piece carbon footplates options DDP1578, DDP1579, DDP1566 & DDP1567 or 2 piece plastic footplates DDP1562 & DDP1563

Knee to heel must be measured taking away the cushion thickness

BACKREST CHOICE
BACK UPHOLSTERY

DDP0700	Standard, black nylon upholstery, no tension adjustment	(NC: 0703,1491)	◆
DDP0721	Lightweight mesh back upholstery, no tension adjustment	(NC: 0703,1491)	£50
DDP0722	Velcro tension adjustable backrest with black & grey nylon cover		£100
DDP0735	Velcro tension adjustable backrest with breathable spacer black cover		£163

Order No:

Customer Ref:

 Kuschall Compact
2 of 9

 Subject to Invacare's standard Terms and Conditions.
Invacare reserve the right to alter product specification without prior notice.

PUSH HANDLES																								
DDP0702	Push-handles fixed, long	In combination with Matrx backrests, push handle height will always be approximately 475mm from the seat	◆																					
DDP0750	Without push handles	(NC: 1304)	○																					
DDP0751	Push handles fixed, mini	(NC: 1304) In combination with Matrx backrests, push handle height will always be approximately 475mm from the seat	○																					
DDP0703	Height adjustable, rear set	(NC: 0734,1304,1411) Not compatible with Matrx backrests	£185																					
DDP0733	Height adjustable, integrated	(NC: 1304,1491)	£160																					
DDP0704	Foldable push handles.	(NC: 1304,1491) Not compatible with Matrx backrests	£144																					
BACKBRACE																								
DDP1491	Backrest stabilization bar, swingaway and removable (Not compatible with Matrx backrests or with backrest heights lower than 405mm)		£120																					
BACK TUBES																								
DDP0736	Fixed backrest tubes set at 90°		◆																					
DDP0705	Lumbal tubes fixed at 90°	(NC: 0734) Not compatible with Matrx backrests	○																					
DDP0720	Angle adjustable back (Backrest angle must be selected from the 6 options below)	(NC: 0734) Not compatible with Matrx backrests	£147																					
DDP0734	Fold down back (folds forward onto the seat) Select angle choice from list opposite. Angle selected will be factory set but can be adjusted to the corresponding angle adjustment possibilities. Not compatible with Tubular Armrests, Option DDP1411 Not compatible with Matrx backrests	<table border="1"> <thead> <tr> <th colspan="3">Folding Back Angle Possibilities</th> </tr> <tr> <th>Back Angle Choice</th> <th colspan="2">Adjustment Possibilities</th> </tr> </thead> <tbody> <tr> <td>82°</td> <td>86°</td> <td>90°</td> </tr> <tr> <td>86°</td> <td>82°</td> <td>90°</td> </tr> <tr> <td>90°</td> <td>82°</td> <td>86°</td> </tr> <tr> <td>94°</td> <td>90°</td> <td>98°</td> </tr> <tr> <td>98°</td> <td>90°</td> <td>94°</td> </tr> </tbody> </table>	Folding Back Angle Possibilities			Back Angle Choice	Adjustment Possibilities		82°	86°	90°	86°	82°	90°	90°	82°	86°	94°	90°	98°	98°	90°	94°	£54
Folding Back Angle Possibilities																								
Back Angle Choice	Adjustment Possibilities																							
82°	86°	90°																						
86°	82°	90°																						
90°	82°	86°																						
94°	90°	98°																						
98°	90°	94°																						
BACKREST ANGLE Must be selected with angle adjustable backrest (0720) and fold down backrest (0734)																								
DDP0707	Angle adjustable backrest, adjusted 82°		○																					
DDP0718	Angle adjustable backrest, adjusted 86°		○																					
DDP0709	Angle adjustable backrest, adjusted 90°		○																					
DDP0719	Angle adjustable backrest, adjusted 94°		○																					
DDP0711	Angle adjustable backrest, adjusted 98°		○																					
DDP0712	Angle adjustable backrest, adjusted 102°		○																					
MATRIX BACKREST (Hardware is supplied mounted to the backposts with the backrest supplied in the box with the wheelchair, extra leadtime applies)																								
DDP0753	Matrx MX1 with Dual Mount Q/R hardware Height of backrest 230mm	Matrx backrests are not available with height adjustable push handles, option DDP0703 or foldable push handles, option DDP0704	£595																					
DDP0754	Matrx MX1 with Dual Mount Q/R hardware Height of backrest 310mm		£595																					
DDP0755	Matrx MX1 with Dual Mount Q/R hardware Height of backrest 410mm		£595																					
DDP0756	Matrx Elite Deep with Easy-set mounting hardware. Height of backrest 310mm	Compatible with seat widths 360mm - 500mm	£450																					
DDP0758	Matrx Elite Standard with Easy-set mounting hardware. Height of backrest 260mm	Compatible with seat widths 360mm - 500mm	£406																					
DDP0759	Matrx Elite Standard with Easy-set mounting hardware. Height of backrest 310mm	Compatible with seat widths 360mm - 500mm	£406																					
DDP0760	Matrx Elite Standard with Easy-set mounting hardware. Height of backrest 360mm	Compatible with seat widths 360mm - 460mm	£406																					
SEATING																								
DDP0611	Standard adjustable seat upholstery with 1 central velcro band (soft)		◆																					
DDP0614	Seat upholstery with 1 central (loop) and 2 side velcro bands (hook)		£25																					
SEATING OPTIONS																								
DDP0600	Standard seat cushion, 6cm foam with Visco top layer		£48																					
DDP0662	Matrx Libra with Black Smooth cover Very high risk pressure & positioning cushion	Matrx Libra is compatible with seat sizes 360-500mm wide & 400-500mm deep	£319																					
DDP0663	Matrx Vi with Black Meshtex cover High risk pressure & positioning cushion	Matrx Vi is compatible with seat sizes 360-500mm wide & 400-500mm deep	£335																					
DDP0606	Seat Bag (click in and removable)		£53																					
LEGRESTS																								
DDP1540	Legrests 70°		◆																					
DDP1541	Legrests 80°		○																					
DDP1542	Legrests 90°		○																					
DDP1542	Not possible with 6" or 7" front castors		○																					
FOOTREST For Seat Width 280-320mm only 1pc or 2pc Carbon footplates are possible																								
DDP1562	2 pc. plastic, set to front		◆																					
DDP1563	2 pc. plastic, set to back		○																					
DDP1564	2 pc. plastic angle adjustable set to front		£95																					
DDP1565	2 pc. plastic angle adjustable set back		£95																					
DDP1578	2 pc. carbon footplates, angle adjustable set to the front		£244																					
DDP1579	2 pc. carbon footplates, angle adjustable set to the rear		£244																					
DDP1566	1 pc. carbon Flip-up, adjustable, fixed on 1 side when folded, set to front		£122																					
DDP1567	1 pc. carbon Flip-up, adjustable, fixed on 1 side when folded, set to back		£122																					
FOOT AND LEG SUPPORT																								
DDP1534	Calf strap		○																					
DDP1530	Heel loops		○																					

SIDE GUARDS																											
CLOTHESGUARDS & MUD GUARDS																											
DDP1400	Clothes guards carbon, fixed	Clothesguards and mudguards are not compatible with armrest options (DDP1416,1421,1440,1441) clothesguards are included on armrests.	◆																								
DDP1405	Mud guards carbon, fixed (NC: 1118,1119)		£115																								
DDP1408	Mud guards carbon, removable (Adds 20mm to the total width of the chair)		£205																								
ARMRESTS (Includes clothesguards, adds 20mm to total width of the chair)																											
DDP1416	Armrest height adjustable (with tools), removable, short, non-locking		○																								
DDP1421	Armrest height adjustable (with tools), removable, long, non locking		○																								
DDP1440	Armrest with desk function, flip-back, height adjustable, removable, short (NC: 1302,0283)		£155																								
DDP1441	Armrest with desk function, flip-back, height adjustable, removable, long (NC: 1302,0283)		£185																								
DDP1411	Tubular Armrest with comfort padded arms. Non locking, height adjustable, removable and swingaway. Does not include clothesguards		£175																								
FRAME COLOURS																											
DDP0300	Glossy-Black	Select either full frame colour, option DDP0398 or rear frame and crossbar matt black option DDP0399	○																								
DDP0314	Black Matt		○																								
DDP0356	Ivory White		○																								
DDP0302	Royal Blue		○																								
DDP0303	Royal Red		○																								
DDP0350	Purple Fire		○																								
DDP0315	Silver Grey		○																								
DDP0354	Night Shimmer		○																								
DDP0357	Umber Glaze		○																								
DDP0358	Smokey Grey		○																								
DDP0359	Flashy Yellow		○																								
DDP0365	Mercury Blue (Available with option DDP0399 only)		These 4 colours are only available with rear frame and crossbar in matt black with front frames in selected colour Select option DDP0399	○																							
DDP0366	Dynamite Orange (Available with option DDP0399 only)			○																							
DDP0367	Atomic Violet (Available with option DDP0399 only)			○																							
DDP0368	Nitro Pink (Available with option DDP0399 only)	○																									
FRAME COLOUR TYPE																											
DDP0398	Full frame colour, front, rear frame and crossbar painted in selected colour		○																								
DDP0399	Rear frame and crossbar matt black, front frame in selected colour		○																								
COLOUR-POWER-SET																											
Consists of selected colour on castor forks, Starec rear wheel hub, rim & Q/R button, seat rail caps, standard rear wheel fixation, castor housing block																											
DDP0330	Colour-Power-Set techno-anthracite		£163																								
DDP0334	Colour-Power-Set b-blue		£163																								
DDP0329	Colour-Power-Set comet red		£163																								
DDP0337	Colour-Power-Set copper glaze		£163																								
DDP0336	Colour-Power-Set arctic blue		£163																								
DDP0338	Colour-Power-Set chrome	Starec wheel rim supplied in black when Chrome CPS is selected	£163																								
DDP0339	Colour-Power-Set champagne		£163																								
DDP0340	Colour-Power-Set olive		£163																								
FRONT WHEEL ATTACHMENT																											
DDP0864	High mounted (Frame will be 30mm lower)	<table border="1"> <thead> <tr> <th colspan="4">Front Castor/Front Seat Height Matrix</th> </tr> <tr> <th>Front castor size</th> <th>Castor attachment High mounted</th> <th>Castor attachment Standard mounted</th> <th>Castor attachment Low mounted</th> </tr> </thead> <tbody> <tr> <td>4"</td> <td>410 - 480mm</td> <td>440 - 510mm</td> <td>470 - 530mm</td> </tr> <tr> <td>5"</td> <td>440 - 490mm</td> <td>470 - 520mm</td> <td>500 - 530mm</td> </tr> <tr> <td>6"</td> <td>460 - 510mm</td> <td>490 - 530mm</td> <td>520 - 530mm</td> </tr> <tr> <td>7"</td> <td>480 - 520mm</td> <td>510 - 530mm</td> <td>540 - 550mm*</td> </tr> </tbody> </table>	Front Castor/Front Seat Height Matrix				Front castor size	Castor attachment High mounted	Castor attachment Standard mounted	Castor attachment Low mounted	4"	410 - 480mm	440 - 510mm	470 - 530mm	5"	440 - 490mm	470 - 520mm	500 - 530mm	6"	460 - 510mm	490 - 530mm	520 - 530mm	7"	480 - 520mm	510 - 530mm	540 - 550mm*	○
Front Castor/Front Seat Height Matrix																											
Front castor size	Castor attachment High mounted		Castor attachment Standard mounted	Castor attachment Low mounted																							
4"	410 - 480mm	440 - 510mm	470 - 530mm																								
5"	440 - 490mm	470 - 520mm	500 - 530mm																								
6"	460 - 510mm	490 - 530mm	520 - 530mm																								
7"	480 - 520mm	510 - 530mm	540 - 550mm*																								
DDP0865	Standard mounted		◆																								
DDP0866	Low mounted (Frame will be 30mm higher)	* These measurements can be achieved as a special build, contact our specials department for further details	○																								

FRONT WHEEL			
DDP0801	100mm (4") Everyday, black tyre		○
DDP0860	100mm (4") Comfort wheel with black tyre		○
DDP0803	100mm (4") Starec techno-anthracite		£46
DDP0807	100mm (4") Starec b-blue		£46
DDP0808	100mm (4") Starec comet red		£46
DDP0830	100mm (4") Starec copper glaze		£46
DDP0831	100mm (4") Starec arctic blue		£46
DDP0809	100mm (4") Starec black		£46
DDP0861	100mm (4") Starec chrome		£46
DDP0825	100mm (4") Starec champagne		£46
DDP0826	100mm (4") Starec olive		£46
DDP0810	125mm (5") Light Wheel, black tyre		○
DDP0863	125mm (5") Comfort Wheel black tyre		◆
DDP0813	125mm (5") Starec techno anthracite		£46
DDP0817	125mm (5") Starec b-blue		£46
DDP0818	125mm (5") Starec comet red		£46
DDP0833	125mm (5") Starec copper glaze		£46
DDP0832	125mm (5") Starec arctic blue		£46
DDP0819	125mm (5") Starec black		£46
DDP0862	125mm (5") Starec chrome		£46
DDP0827	125mm (5") Starec champagne		£46
DDP0828	125mm (5") Starec olive		£46
DDP0820	150mm (6") Light wheel, black tyre		○
DDP0821	150mm (6") Pneumatic Wheel		£19
DDP0822	175mm (7") Light Wheel		○
ONE ARM DRIVE			
ONE-ARM-DRIVE SYSTEM (Complete system with pneumatic light wheel 7 bar tyres & aluminium handrims and black rear wheel fixation, CPS not available)			
DDP1118	One-arm-drive Right side (Not possible with spokeguards)	(NC: 1020,1304,0909,0283,0249,0266,0267,0268,1102,1104,1106,1120)	£636
DDP1119	One-arm-drive Left side (Not possible with spokeguards)	(NC: 1020,1304,0909,0283,0249,0266,0267,0268,1102,1104,1106,1120)	£636
REAR WHEEL ADAPTION PLATE			
DDP0250	Standard fixation, 4 rear wheel positions		◆
DDP0249	Rear wheel extension, 6 rear wheel positions (Always supplied black in colour)		£48
REAR WHEEL POSITION			
DDP0263	1st position	(NC: 0249,0283,1304)	○
DDP0264	2nd position		○
DDP0265	3rd position		◆
DDP0266	4th position		○
DDP0267	5th position	5th & 6th position not for standard fixation DDP0250	○
DDP0268	6th position		○
QUAD QUICK RELEASE AXLE			
DDP0909	Quad quick release axle	(NC: 0906,1304)	£62
REAR WHEEL SIZE			
DDP0901	559mm (22")	(NC: 0906,0908,1104,1106,1120)	○
DDP0902	609mm (24")		◆
DDP0903	660mm (26")	(NC: 0283,0906,0908)	○
REAR WHEEL CAMBER (3° Not available with drum brakes and performance brakes combined)			
DDP0282	Camber 0°		◆
DDP0283	Camber 3°	(NC: 1302,1304)	○
REAR WHEEL TYPE			
DDP0915	22", 24", 26" Kuschall Starec (Silver crossed spokes, hub colourless, rim black. Hub & rim in CPS if selected)		◆
DDP0908	24" High Performance (Starec hub, spokes & rim black, hub in CPS if selected)	(NC: 1004,1006,1007,1304)	£185
DDP0906	24" Spinergy LXX wheel (Black lightweight hub, 18 spokes black, No CPS) Coloured spokes are available via MyKuschall, please indicate requirements in the space provided at the end of the form.	(NC: 1004,1006,1007,1200,1202)	£445
SWO0907	24" Spinergy SPOX wheel (Chrome hub, 18 spokes in black, no CPS) Coloured spokes are available via MyKuschall, please indicate requirements in the space provided at the end of the form.	(NC: 1004,1006,1007,1200,1202)	£340
SWO0909	24" Spinergy LX wheel (Chrome hub, 12 spokes in black no CPS). Coloured spokes are available via MyKuschall, please indicate requirements in the space provided at the end of the form.	(NC: 1004,1006,1007,1200,1202)	£529
REAR WHEEL TYRE			
DDP1014	22", 24", 26" Schwalbe Rightrun black & grey tyre		◆
DDP1002	22", 24" Profile tyre 7 bar		○
DDP1004	24" Solid tyre, profil, grey		£48
DDP1006	24" Solid tyre (KIK type), black		£63
DDP1007	22", 24" Solid, light wheel		£48
DDP1010	24", 26" Marathon Plus tyre with Smartguard, 5mm layer of protective rubber infill	(NC: 1020)	£73
DDP1011	24" Schwalbe One tyre, V guard puncture resistant		£62
ANTI-FLAT			
DDP1020	Anti Flat tyre tube protection		£15

Order No:

Customer Ref:

Kuschall Compact
5 of 9Subject to Invacare's standard Terms and Conditions.
Invacare reserve the right to alter product specification without prior notice.

HANDRIM																																																																								
DDP1113	22", 24", 26" Light Alu grey, short tab, 1 hole set at 3 cm distance		◆																																																																					
DDP1114	22", 24", 26" Light Alu grey, long tab, 2 hole set at 4 cm distance		○																																																																					
DDP1117	22", 24", 26" Light Alu black, short tab, 1 hole set at 3 cm distance		£30																																																																					
DDP1102	22", 24", 26" Supergripp, long tab, 2 hole set at 3 cm distance		£145																																																																					
DDP1104	24", 26" Steel polished, short tab, 1 hole set at 3 cm distance		£45																																																																					
DDP1106	24" Titanium, short tab, 1 hole set at 3 cm distance		£128																																																																					
DDP1120	24" Surge LT 30mm handrim, 3cm distance		£237																																																																					
SPOKE GUARDS																																																																								
DDP1200	22", 24", 26" Transparent	(NC: 1304)	£85																																																																					
DDP1202	24" Kuschall Black with white logo	(NC: 1304)	£85																																																																					
SAFETY																																																																								
BRAKES																																																																								
DDP1306	Pull brake		◆																																																																					
DDP1300	Push brake, bent	(NC: 1310)	○																																																																					
DDP1301	Push brake lever straight		○																																																																					
DDP1302	Performance Brake	(NC: 1304,1310)	£55																																																																					
DDP1303	Scissor Brake	(NC: 1310)	£55																																																																					
<table border="1"> <thead> <tr> <th colspan="7">Scissor Brake compatibilities on a 24" Rear Wheel</th> </tr> <tr> <th rowspan="2">Seat Depth</th> <th colspan="6">Rear Wheel Position</th> </tr> <tr> <th>1</th> <th>2</th> <th>3</th> <th>4</th> <th>5</th> <th>6</th> </tr> </thead> <tbody> <tr> <td>32cm-38cm</td> <td>N</td> <td>N</td> <td>N</td> <td>N</td> <td>Y</td> <td>Y</td> </tr> <tr> <td>40cm</td> <td>N</td> <td>N</td> <td>N</td> <td>N</td> <td>Y</td> <td>Y</td> </tr> <tr> <td>42cm</td> <td>N</td> <td>N</td> <td>N</td> <td>Y</td> <td>Y</td> <td>Y</td> </tr> <tr> <td>44cm</td> <td>N</td> <td>N</td> <td>Y</td> <td>Y</td> <td>Y</td> <td>Y</td> </tr> <tr> <td>46cm</td> <td>N</td> <td>Y</td> <td>Y</td> <td>Y</td> <td>Y</td> <td>Y</td> </tr> <tr> <td>48cm</td> <td>Y</td> <td>Y</td> <td>Y</td> <td>Y</td> <td>Y</td> <td>Y</td> </tr> <tr> <td>50cm</td> <td>Y</td> <td>Y</td> <td>Y</td> <td>Y</td> <td>Y</td> <td>Y</td> </tr> </tbody> </table>				Scissor Brake compatibilities on a 24" Rear Wheel							Seat Depth	Rear Wheel Position						1	2	3	4	5	6	32cm-38cm	N	N	N	N	Y	Y	40cm	N	N	N	N	Y	Y	42cm	N	N	N	Y	Y	Y	44cm	N	N	Y	Y	Y	Y	46cm	N	Y	Y	Y	Y	Y	48cm	Y	Y	Y	Y	Y	Y	50cm	Y	Y	Y	Y	Y	Y
Scissor Brake compatibilities on a 24" Rear Wheel																																																																								
Seat Depth	Rear Wheel Position																																																																							
	1	2	3	4	5	6																																																																		
32cm-38cm	N	N	N	N	Y	Y																																																																		
40cm	N	N	N	N	Y	Y																																																																		
42cm	N	N	N	Y	Y	Y																																																																		
44cm	N	N	Y	Y	Y	Y																																																																		
46cm	N	Y	Y	Y	Y	Y																																																																		
48cm	Y	Y	Y	Y	Y	Y																																																																		
50cm	Y	Y	Y	Y	Y	Y																																																																		
DRUM BRAKE (Specific complete wheel mounted with a black adapter plate) CPS if selected will not apply to this wheel option																																																																								
DDP1304	Attendant drum brake		£267																																																																					
BRAKE LEVER EXTENSION																																																																								
DDP1310	Brake extension lever		£32																																																																					
ANTI TIPPER Only 2 adaptions per side. Adaptions include Anti tips, Tipper aid and Cane holder																																																																								
DDP1702	Anti tipper right		£89																																																																					
DDP1703	Anti tipper left		£89																																																																					
DDP1704	Anti tipper both sides		£178																																																																					
SEAT BELT																																																																								
DDP1730	Seat belt	(NC: 1770)	£54																																																																					
DDP1770	Transportation kit including seat belt and tie down labels		£54																																																																					
ACCESSORIES																																																																								
DDP1729	Rear frame protector		£13																																																																					
DDP1746	Mobile phone holder		£20																																																																					
DDP1712	Kuschall watch mounted on right castor housing		£77																																																																					
DDP1713	Kuschall watch mounted on left castor housing		£77																																																																					
DDP1734	Mountain bike wheels supplied with bag, includes quick release axles		£199																																																																					
TIPPER AID																																																																								
DDP1708	Tipper-Aid, right	(NC: 1706)	£31																																																																					
DDP1709	Tipper-Aid, left	(NC: 1707)	£31																																																																					
DDP1718	Tipper-Aid, both sides	(NC: 1706,1707,1719)	£62																																																																					
CANE HOLDER																																																																								
DDP1706	Cane holder, right		£47																																																																					
DDP1707	Cane holder, left		£47																																																																					
DDP1719	Cane holder, both sides		£94																																																																					
MANUAL AIR PUMP																																																																								
DDP1715	Air pump		£10																																																																					
TOOL KIT																																																																								
DDP1717	Tool kit		◆																																																																					
PRICE																																																																								
			Total Price £																																																																					




Order No:

Customer Ref:

Kuschall Compact
6 of 9

Subject to Invacare's standard Terms and Conditions.
Invacare reserve the right to alter product specification without prior notice.

Non Conforming Items Guide	
DDP0249	Rear wheel extension
DDP0266	4th position
DDP0267	5th position
DDP0268	6th position
DDP0283	Camber 3°
DDP0703	Height adjustable rear set
DDP0712	Angle adjustable backrest, adjusted 102°
DDP0734	Foldable back (to the front)
DDP0820	6" light wheel
DDP0821	6" pneumatic wheel
DDP0822	7" Light Wheel
DDP0906	24" Spinergy LX wheels
DDP0908	24" High Performance (Starec hub)
DDP0909	Quad quick release axles
DDP1004	24" Solid tyre, profil, grey
DDP1006	24" Solid tyre (KIK type), black
DDP1007	22", 24" Solid, light wheel
DDP1020	Anti-flat tyre tube protection
DDP1102	22", 24", 26" Supergripp, dist. 3 cm
DDP1104	24", 26" Steel, dist. 3 cm
DDP1106	24" Titanium, dist. 3 cm
DDP1118	One arm drive system RH
DDP1119	One arm drive system LH
DDP1120	Surge LT handrim
DDP1200	Transparent spokeguards
DDP1202	Kuschall spokeguards
DDP1300	Push brake bent
DDP1301	Push brake lever straight
DDP1302	Performance brake
DDP1304	Attendant drum brake
DDP1310	Brake extension lever
DDP1405	Mud guards carbon fixed
DDP1411	Tubular armrests
DDP1416	Siderest height adjustable (with tools), removable, short
DDP1421	Siderest height adjustable (with tools), removable, long
DDP1440	Siderest with desk function, flip-back, height adjustable, removable, short
DDP1441	Siderest with desk function, flip-back, height adjustable, removable, long
DDP1491	Backbrace swing away and removable
DDP1530	Heel loops
DDP1562	2 pc. plastic, set to the front
DDP1563	2 pc. plastic, set back
DDP1564	2 pc. angle adjustable set to the front
DDP1565	2 pc. angle adjustable, set back
DDP1566	1 piece carbon flip up adjustable footplate set to front
DDP1567	1 piece carbon flip up adjustable footplate set to back
DDP1702	Anti tipper right
DDP1703	Anti tipper left
DDP1706	Cane holder, right
DDP1707	Cane holder, left
DDP1719	Cane holder both sides
DDP1770	Transportation Kit includes seat belt & tie down labels (available on back heights 405mm and above)

	User measurement*	Wheelchair measurement
 <p>Seat width (SB)</p> <p>The seat width determines the user's comfort and performance. It must be wide enough to ensure good weight distribution, but narrow enough to allow easy manipulation of the wheelchair.</p>	 <ul style="list-style-type: none"> Let the person sit on a flat surface. Position vertical blocks left and right of the hips and measure the distance between the blocks. Add approximately 2 cm to cater for winter clothing for the wheelchair measurement. 	 <ul style="list-style-type: none"> Remove the seat cushion. Measure the distance between the outer edges of the seat cover at the rear of the wheelchair. <p>Ultra-Light / Compact: For the seat width add 0,5 cm to the measured distance.</p>
 <p>Seat depth (ST)</p> <p>A correct seat depth measurement enables an even distribution of the user's weight across the whole seating base, thereby avoiding pressure points and skin irritation in the area behind the knee.</p>	 <ul style="list-style-type: none"> Let the person sit on a flat surface. Measure the distance from the rear most point of the buttocks/lower back to the hollow at the back of the knee. Subtract approximately 4 cm (depending on user's physique and the legrest angle) to give the wheelchair measurement. 	 <ul style="list-style-type: none"> Remove the seat cushion. Measure the distance from the front edge of the backrest tube (without upholstery) to the front edge of the seat cover.
 <p>Leg length (UL)</p> <p>The leg length measurement establishes the correct height for the footplates, which will enable part of the body weight to be taken by the feet. Footplates must have at least 2 cm clearance above the ground.</p>	 <ul style="list-style-type: none"> Let the person sit on a flat surface. Measure the distance from the heel of the shoe to the hollow at the back of the knee. Subtract any cushion thickness to give the wheelchair measurement. 	 <ul style="list-style-type: none"> Remove the seat cushion. Measure the distance from the back edge of the footplate to the front edge of the seat cover (without sagging of cover).
 <p>Front seat height (SHv)</p> <p>A correct seat height measurement will ensure the optimum driving efficiency and can give more clearance to fit underneath tables.</p>	 <ul style="list-style-type: none"> Let the person sit on a flat surface. Measure the vertical distance from the ground to the hollow at the back of the knee. This measurement should not be less than the leg length (UL) + 4 cm to allow sufficient footplate clearance above the ground. Subtract any cushion thickness to give the front seat height (SHv) wheelchair measurement. 	 <ul style="list-style-type: none"> Remove the seat cushion. Measure the vertical distance from the floor to the top of the front edge of the seat cover (without sagging of cover).
 <p>Rear seat height (SHh)</p>	<ul style="list-style-type: none"> Rear seat height is determined by the front seat height, the seat angle and the seat depth. 	 <ul style="list-style-type: none"> Remove the seat cushion. Measure the vertical distance from the floor to the top of the rear edge of the seat cover (without sagging of cover).
 <p>Backrest height (RH)</p> <p>The backrest height is dependent on the user's degree of disability. A correct backrest height should maintain a good posture whilst allowing the permitted freedom of movement.</p>	 <ul style="list-style-type: none"> Let the person sit on a flat surface. Measure the vertical distance from the seat platform to the scapulas. If more trunk support is required, measure from the seat platform up to the required level of support. If less trunk support is required reduce the backrest height in order to improve activity. 	 <ul style="list-style-type: none"> Remove the seat cushion. Remove back cover. Measure the distance from the top of the seat frame to the top of the backrest straps.

*To be done by a specially trained person

Matrix Seating & Positioning Options - THE COMPLETE SOLUTION

The Invacare Seating & Positioning range features a vast array of pressure redistributing cushions from low to very high risk, and a number of backrest styles to suit a variety of postural needs. Below is a selection of the options available on chair, further information on the complete solution can be found on individual prescription forms or alternatively visiting us online at www.invacare.co.uk

Prescription Guide

Effective seating should provide postural support and comfort, as well as help to reduce the risk of pressure ulcer development. There are a number of extrinsic factors that can influence the effectiveness of a seat cushion in relation to pressure care, our cushions have been designed taking these factors into account explained below.

Cushion Options - (For additional sizes please refer to the individual cushion prescription form)

Matrix Flo-tech Lite

Slimline contoured cushion ideally suited to people at low risk of developing pressure ulcers and for comfort purposes only. The contoured seating surface offers improved stability and postural support.



Low Risk
Overall height 3"
Max UW 127kg (20 stone)

Matrix Vi

The unique three part foam construction targets skin protection where it is needed most. The Invacare Matrix Vi Cushion is anatomically designed with dual-density, HR foam layers and includes a Visco-elastic foam insert, pre-ischial contour and trochanteric shelf.



High Risk
Overall height 3.25"
Max UW 136kg (21 stone)

Matrix Libra

Designed with the highest level of skin protection, positioning and adjustability to match the unique needs of the person, the Invacare Matrix Libra Cushion offers a complete solution to provide balance and stability as well as optimize body function.



Very High Risk
Overall height 3.25"
Max UW 136kg (21 stone)

Backrest Options - (For more information please refer to the Matrix brochure)

Matrix MX1

A lightweight carbon fibre back, designed for active people who demand peak performance and great looks.

- Lightweight carbon fibre back shell - weighs from approximately 900 grams, including hardware
- Full range of anterior / posterior angle adjustment
- 3" / 75mm of contour depth provides stable midline support
- Carbon fibre shell delivers exceptional support, durability and aesthetic appeal
- Range of three patent pending hardware options available for fitting the back to a wheelchair
- Hardware options provide optimal adjustment for the widest range of applications, ensuring a great fit for the person.

Hardware:

Quick release dual mount back cane hardware with quick release tabs.




Matrix Elite

Designed to provide precise orientation within the wheelchair for optimal postural support. The attractive, lightweight design appeals to people with active lifestyles.

Features :

- 3" Contour Depth
- 40° angle adjustment
- 4" depth adjustment
- 20° rotation
- 136kg (21 stone) Max UW





Matrix Elite Deep

Featuring integrated lateral supports to provide trunk support while allowing clearance for the shoulders and hips.

Features :

- 6" Contour Depth
- 40° angle adjustment
- 4" depth adjustment
- 20° rotation
- 136kg (21 stone) Max UW



Bodypoint - Positioning Accessories (Currently only available as spare parts - please refer to the UK Bodypoint price list)

Belts

A selection of 2 and 4 point padded hip belts are available, that facilitate aligned, stable pelvic position.

Multiple pull options and buckle styles accommodate a variety of needs



Harnesses

Designed to provide shoulder retraction, good head control and to correct shoulder rotation, Bodypoint shoulder harnesses assure comfort, durability and easy care.



Foot Supports

Used for stabilising and positioning the feet while still allowing for controlled movement

A comfortable and dynamic alternative to rigid shoe holders or foot straps



Customer Ref:

Order Number:

Subject to Invacare's standard terms conditions. Invacare reserve the right to alter product specification without prior notice.