



*A No Limits range
with personality at the heart
of its design*

Hannah Cockroft MBE
Paralympian



1 KÜSCHALL K-SERIES

The küschall K-Series represents a breakthrough in innovative wheelchair design, suitable for highly active people who demand aesthetic appeal together with exceptional strength and performance.

Available in a range of materials, including aluminium, titanium and carbon, the küschall K-Series boasts a 'No Limits' approach, offering dynamic frames that can withstand the rigours of everyday use, and suitable for even the most extreme wheelchair drivers.

Versatile and configurable to suit a wide range of personal requirements, the küschall K-Series is for anyone, of any age who demands effortless mobility combined with great style.

TECHNICAL DATA

 Seat width: 340 – 480 mm in increments of 20 mm	 Seat depth: Al 375 – 525 mm Ti C 375 – 450 mm in increments of 20 mm	 Seat height: Al front: 450 – 520 mm Ti C front: 480 – 520 mm rear: 380 – 490 mm in increments of 10 mm
 Knee-to-heel length: 290 – 480 mm in increments of 10 mm	 Backrest height: 270 – 480 mm in increments of 15 mm	 Centre of gravity: 80 mm/ 95 mm/ 110 mm/ 125 mm/ 140 mm
 Backrest angle: 74°/78°/82°/86°/90°	 Max. safe slope: 7°	 Reversing width: 920 – 1230 mm
 Total width: SW + 170 mm	 Total length: 75°: approx. 830 mm 90°: approx. 760 mm*	 Max. user weight: 120 kg
 Total weight unloaded: Al Ti approx. 8.1 kg C approx. 7.9 kg	 Transport weight w/out rear wheels: Al Ti approx. 5.5 kg C approx. 5.3 kg	 Crash Test approved Fulfills the requirements of ISO 7176-19

LEGEND

- CONFIGURATION:** Product shown has an aluminium standard frame in the frame colour ivory white, carbon camber bar, carbon mud-guards, standard tension adjustable with Velcro® black/grey backrest, standard aluminium foot-tube, active brakes, 4" Starec front wheels and high performance rear wheels with Schwalbe Marathon tyres and aluminium grey handrims.
- FRAME:** Küschall K-Series wheelchairs come in a choice of three exceptional materials for the front frame: aluminium, titanium and carbon, each with its own attributes and benefits. Please refer to the Prestige product brochure for further details.
- BACKREST:** Available either in aluminium (black) or titanium, the backrest bar has a rounded shape to facilitate access to the back pocket, which can be easily closed with the magnetic clasp. The two-tone backrest upholstery (black/grey) comes as standard on the küschall K-Series.
- CAMBER BAR FIXATION:** The standard aluminium camber bar is designed with a fixation in two tones: black and anodized silver - attention to detail that gives this wheelchair a truly stylish finish!
- COLOUR-POWER-SET:** The CPS option applies an anodised colour to specific parts of the wheelchair which creates a true personal touch. Available in a wide range of colours, check local prescription form for further details.
- TRANSPORTATION:** The backrest of the küschall K-Series can easily be folded and the rear wheels dismantled so the wheelchair can fit into the smallest of cars.

Al Aluminium / **Ti** Titanium / **C** Carbon / * not for Carbon

United Kingdom: Invacare Ltd. | Pencoed Technology Park, Pencoed, Bridgend CF35 5AQ | Phone +44 (0) 1656 776222 | Fax +44 (0) 1656 776220 | www.invacare.co.uk
Ireland: Invacare Ltd. | Unit 5, Seatown Business Campus | Seatown Road | Swords, Co. Dublin | Phone +353 (0)1 810 7084 | Fax +353 (0)1 810 7085 | www.invacare.ie

küschall is a registered trademark. Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Küschall AG reserves the right to modify existing models at any time. Any use of this publication or of parts thereof as well as any reproduction of images must have the written consent of Küschall AG. Patents pending K-Series_UK_2014_1539225 Copyright © 10/2014, Küschall AG, Switzerland - All rights reserved



2



3



4



5



6



KÜSCHALL K-SERIES