



Life Family

Find the one to suit your lifestyle



The new Life Family

Live Life to the full.

The NEW Quickie Life Family gives you more freedom than ever.

With the choice of a RIGID frame for performance or a FOLDING frame for easier transfers, it is easy to find a Life to suit you.

Quickie Life, folding
fixed front



Quickie Life, folding
swing-away



Quickie Life R, rigid
fixed front



Quickie Life R, rigid
swing-away



The Quickie Life Family >>> lightweight >>> robust >>> stylish

NEW Life R

Active Life



Rigid, durable and robust the Life R is designed to enhance your life now and for years into the future.

Available as fixed front or swing-away, the new Life R has a modular design which means the front of the frame can be changed to increase seat depth, or alter the frame type. This means it's perfect for growing teenagers and makes it easy to refurbish.



Whatever the day throws at you, the energy efficient design of the Life R means you can focus on what's important to you.



Quickie Life



Easy Life

In a business meeting this morning? Preparing for the great outdoors this afternoon? Whatever your plans your Quickie Life is ready.

If transportation and a small folding package are key, choose the Life folding. With full adjustability and configurability the Quickie Life is perfect for those who lead an active lifestyle.

But why stop there? Customise your Life by selecting from a wide range of side guards, armrests, hangers, footrests, wheels and handrims. Or define your style with a choice of 10 frame colours and anodized parts.



Combine the modern design with new frame colours and anodized parts for a stylish chair to match your life.



Quickie Life R Rigid



Swing-away



Configure your Life with a choice of frame styles, angles and lengths.



Seat height adjustable in seven positions

With a camber from 0° to 4° and four centre of gravity positions (from the very active 8.3cm to the more steady 2.3cm).



Simple and quick to adjust fork angle

With increased angle range to ensure optimal driving performance.



Angle adjustable back (optional)

Lightweight and simple to adjust. Provides 20° of angle adjustment from -12° to +8° in 4° increments.



Lightweight height adjustable single post armrest

Provides easy transfer and good arm support with a low transport weight.



Active frame (optional)

Compact frame (85°) with 2cm inset for more active users.



Hand bike axle

Increased stability to use the chair in combination with a hand bike.



Lightweight platform footrest

Angle and depth adjustable. Available in composite and aluminium.



Active swing-away hangers (optional on swing-away frame)

Offers a modern and active look. Available in 70° and 80°.

The fixed front frame is available in three angles (75°, 85° and 85° with inset) and four frame lengths. The swing-away frame is available in two angles (70° and 80°) a choice of hangers and three frame lengths.

Quickie Life Folding

Quickie Life has been designed to meet the demands of modern life. Flexible, lightweight, folding and stylish - it won't hold you back.



Swing-away



Whatever your lifestyle Quickie Life folding can be adjusted and configured to suit. Choose from a fixed front lightweight frame for a more active look, or the swing-away version for even easier folding and transportation both now available with new lightweight options.



NEW Aluminium lightweight hangers in 70° and 80°. 13% lighter than Quickie standard hangers.



NEW Lightweight height adjustable single post armrest
Provides easy transfer and good arm support with a low transport weight.



NEW Lightweight anti-tip for optimal driving performance with 480g weight reduction.



Angle adjustable back (optional)
Lightweight and simple to adjust. Provides angle adjustment from -12° to +12° in 4° increments.



Active frame
Compact frame (85°) with a 2cm inset for more active users. Available with or without a 2cm inset.



NEW folding fixation
Operated with one hand, keeps the chair securely folded.



Flip-up platform footrest
Lightweight, angle and depth adjustable.



Seat height adjustable in seven positions
With a camber from 0° to 4°, and 4 centre of gravity positions (from the very active 8.3 cm to the more stable 2.3 cm).



Technical data sheet

QUICKIE[®]
Live without limits.

	Life fixed front	Life swing-away	Life R fixed front	Life R swing-away
Overall width	Seat width +20cm (8")	SW + 20cm (8")	SW + 20cm (8")	SW + 20cm (8")
Seat width	36 – 50cm (14 – 20")	36 – 50cm (14 – 20")	34 – 48cm (13 – 19")	34 – 48cm (13 – 19")
Back height	25 – 47.5cm (10 – 19")	25 – 47.5cm (10 – 19")	25 – 45cm (10 – 18")	25 – 45cm (10 – 18")
Seat depth	36 – 50cm (14 – 20")	36 – 50cm (14 – 20")	36 – 50cm (14 – 20")	36 – 50cm (14 – 20")
Rear seat height	40 – 50cm (16 – 20")	40 – 50cm (16 – 20")	36 – 45cm (14 – 18")	36 – 45cm (14 – 18")
Front seat height	40 – 54cm (16 – 21")	40 – 54cm (16 – 21")	39 – 50cm (15 – 20")	39 – 50cm (15 – 20")
Back angle	-12° to +12°	-12° to +12°	-12° to +8°	-12° to +8°
Frame front end angle	75°, 85°, 85° with 2cm inset	80° and 70°	75°, 85°, 85° with 2cm inset	80° and 70°
Individual frame colours	10	10	10	10
Individual sling colours	1	1	1	1
Fork and axle plate anodized colours	3	3	3	3
Total chair weight, from	11.5kg	12kg	10.7kg	11.3kg
Maximum user weight	125kg (19.6 stone)	125kg (19.6 stone)	125kg (19.6 stone)	125kg (19.6 stone)



Get SOCIAL with QUICKIE

Keep up-to-date with the latest news on the Quickie Life family and other Quickie products through our Facebook page:

www.facebook.com/QuickieWheelchairsUK

SUNRISE
MEDICAL[®]

Sunrise Medical, Thorns Road, Brierley Hill, West Midlands, DY5 2LD
Tel: +44 (0) 845 605 66 88 Fax: +44 (0) 845 605 66 89
www.SunriseMedical.co.uk