

Invoice to:
 Name:
 Street:
 Town: Postal code:
 E-Mail:
 Tel: Fax:

Deliver to:
 Name:
 Street:
 Town: Postal code:
 E-Mail:
 Tel: Fax:



Order date:
 Purchase order:
 Number required:

Marked for:
 Sunrise order no:

Jay® Balance Standard Features

Fluid pad or Roho Air (single or dual valve) insert to optimize immersion and accommodate positioning needs.
 Positioning accessories available.
 Innovative dual-cover system: inner cover - water-resistant, easy clean
 Outer cover - 3 types: stretch, microclimatic and incontinent
 2 Year guarantee; 6 months on cover.

- Order
- Calculation

Max. user weight:
 150 kg (cushion width 34 to 50 cm)
 225 kg (cushion width 56 to 60cm)

USER WEIGHT kg
 Mandatory information to ensure correct configuration

Reimbursement code:

Application:

Client at high risk of skin breakdown with or without high postural needs.
 Symmetrical to significant asymmetrical posture requiring lateral, forward/rearward stability and/or positioning.
 May be unable to reposition or perform an independent pressure relief.

To order a Jay Balance cushion, follow the steps below to build your part number and then choose your accessories

fixed	Width	Depth	Insert	Outer Cove	Pos. Kit	fixed
JBALC	40	42	FL	M	P	EU
Jay Balance Cushion	STEP 1	STEP 2	STEP 3	STEP 4	STEP 5	

Example: JBALC4042FLMPEU

CUSHION WIDTH AND DEPTH - STEP 1 + STEP 2


Standard Cushion (width: 34 - 50 cm; max user weight 150 kg)
 Bariatric Cushion (width: 56 - 60cm; max user weight 225 kg)

cm	Width										
	34	36	38	40	42	44	46	48	50	56	60
34	○										
36	○	○	○								
38	○	○	○	○	○	○					
40	○	○	○	○	○	○					
42	○	○	○	○	○	○	○	○	○		
44		○	○	○	○	○	○	○	○	○	
46		○	○	○	○	○	○	○	○	○	○
48					○	○	○	○	○	○	○
50						○	○	○	○	○	○
56										○	○
60										○	○
PLA Size	A			B				c			

Base Price	
width 34-48cm	£425
width 50-60cm	£459



For delivery times - see back page
 15 day lead time - all cushions with positioning components



PELVIC LOADING AREA (PLA) INSERT* - STEP 3

- | | | |
|---|------------------|-----------------|
| <input type="radio"/> FLUID - Fluid Pad Option | Part # FL | Price £0 |
| <input type="checkbox"/> AIR Insert Option - Single Valve** | AS | £32 |
| <input type="checkbox"/> AIR Insert Pad Option - Dual Valve** | AD | £32 |

INSERT



STEP 3

(FL, AS or AD)



Fluid - FL



Air single valve - AS



Air dual valve - AD

* The insert pad sizes supplied are based on PLA size and automatically determined when you select the cushion (see matrix)

**Air Bladder Options shipped with instructions, air pump, patch kit and nozzle fittings

OUTER COVER OPTION - STEP 4

- | | | |
|-------------------------------------|-----------------|-----------------|
| <input type="radio"/> Microclimatic | Part # M | Price £0 |
| <input type="radio"/> Stretch | S | £0 |
| <input type="radio"/> Incontinent | I | £0 |

OUTER COVER



STEP 4


(M, S or I)

All cushions are delivered with a water resistant inner cover.

POSITIONING KIT - STEP 5

- | | | |
|--|-----------------|------------------|
| <input type="checkbox"/> Positioning Kit - Contains all items below | Part # P | Price £65 |
|--|-----------------|------------------|
- Lateral Thigh 2.5 cm (1") Left
 - Lateral Thigh 2.5 cm (1") Right
 - Medial Thigh 2.5 cm (1")
 - Lateral Pelvic 2.5 cm (1") Left
 - Lateral Pelvic 2.5 cm (1") Right

POS. KIT



STEP 5

('P' or leave blank)



If you want positioning components, please write "P" in the Step 5 box.

JAY BALANCE SMART PART CODE

	Width	Depth	Insert	Outer Cover	Pos. Kit	
JBALC	40	42	FL	M	P	EU
fixed	STEP 1	STEP 2	STEP 3	STEP 4	STEP 5	fixed

Part Code:

The above part code describes a Jay Balance cushion with 40 width and 42 depth, fluid insert, microclimatic outer cover and positioning components.

Additional information:

The Jay fluid insert may be underfilled or overfilled as required (same as a Jay 2 cushion).

If you require a cushion modification please write in the box below or contact your Sunrise Medical Customer Service.

Modification required:

	Unit price	Total price
JAY Balance Cushion Price		
Option price		
Total Price		

Spare Parts

Cushion base		Price
<input type="checkbox"/> BJBALCwwddEU	Base assembly	£78
ww: Width; dd: Depth		

Spare Outer Cover		
<input type="checkbox"/> CJBALCwwddSEU	Stretch outer cover	£50
<input type="checkbox"/> CJBALCwwddMEU	Microclimatic outer cover	£50
<input type="checkbox"/> CJBALCwwddIEU	Incontinent outer cover	£50
ww: Width; dd: Depth		

Spare Inner Cover		
<input type="checkbox"/> CIFLJBALCwwddEU	Inner cover for fluid cushion	£45
<input type="checkbox"/> CIFLPJBALCwwddEU	Inner cover for fluid and positioning cushion	£50
<input type="checkbox"/> CIAJBALCwwddEU	Inner cover for air cushion	£45
<input type="checkbox"/> CIAPJBALCwwddEU	Inner cover for air and positioning cushion	£50
ww: Width dd: Depth		

Solid Seat		
<input type="checkbox"/> JSIwwdd	Solid Seat insert (wood plate)	£34
ww: Width; dd: Depth		

Fluid & Air inserts *)

PLA A (for cushion widths 34-38cm)

<input type="checkbox"/> FLJBALCAEU	Fluid insert for PLA size A	£227
<input type="checkbox"/> J3CSAASEU	Air insert - single valve - for PLA size A	£302
<input type="checkbox"/> J3CSAADEU	Air insert - dual valve - for PLA size A	£302

PLA B (for cushion widths 40-48cm)

<input type="checkbox"/> FLJBALCBEU	Fluid insert for PLA size B	£227
<input type="checkbox"/> J3CSBASEU	Air insert - single valve - for PLA size B	£302
<input type="checkbox"/> J3CSBADEU	Air insert - dual valve - for PLA size B	£302

PLA C (for cushions widths 50-60cm)

<input type="checkbox"/> FLJBALCCEU	Fluid insert for PLA size C	£227
<input type="checkbox"/> J3CSCAS	Air insert - single valve - for PLA size C	£302
<input type="checkbox"/> J3CSCADEU	Air insert - dual valve - for PLA size C	£302


Jay Balance Positioning kit *)

	for width 34-40 cm	for width 42-60 cm	
Lateral Thigh 2.5cm Left	<input type="checkbox"/> JBALCLTLH-S	<input type="checkbox"/> JBALCLTLH	£13
Lateral Thigh 2.5cm Right	<input type="checkbox"/> JBALCLTRH-S	<input type="checkbox"/> JBALCLTRH	£13
Medial Thigh 2.5cm	<input type="checkbox"/> JBALCMT-S	<input type="checkbox"/> JBALCMT	£13
Lateral Pelvic 2.5 cm Left	<input type="checkbox"/> JBALCLPLH-S	<input type="checkbox"/> JBALCLPLH	£13
Lateral Pelvic 2.5 cm Right	<input type="checkbox"/> JBALCLPRH-S	<input type="checkbox"/> JBALCLPRH	£13
Positioning kit	<input type="checkbox"/> JBALCPKIT-S	<input type="checkbox"/> JBALCPKIT	£65
(Positioning kit contains 1 of each of the components above)			

*) Note: If you change the insert type of your cushion (fluid to air or vice versa) or change your cushion from positioning to non-positioning due to design differences (or vice versa), you will also have to change the inner cover type accordingly.

Jay Balance Measuring Guide

Measurements made from completed cushion with covers



Width X cm

Measurement made from outside edge to outside edge of cushion

Depth Y cm

Measurement made from front of seat notch to front topside (or bottom of contracture cut) of cushion.

Jay Balance Delivery Times

Jay Balance with Fluid Insert

		Width												
		cm	34	36	38	40	42	44	46	48	50	56	60	
Depth	34	◇												
	36	◇	○	◇										
	38	◇	◇	◇	◇	◇	◇							
	40	◇	○	○	○	○	○							
	42	◇	◇	○	○	○	○	○	◇	◇				
	44		◇	○	○	○	○	○	○	◇				
	46		◇	○	○	○	○	○	○	◇	◇			
	48					○	○	○	○	◇	◇	◇		
	50						○	○	○	◇	◇	◇		
	56									◇	◇	◇		
60										◇	◇			
PLA Size		A			B				C					

Delivery times *)

○	5 days
◇	15 days
■	not available

Jay Balance with Air Insert

		Width												
		cm	34	36	38	40	42	44	46	48	50	56	60	
Depth	34	◇												
	36	◇	◇	◇										
	38	◇	◇	◇	◇	◇	◇							
	40	◇	◇	◇	○	◇	◇							
	42	◇	◇	◇	○	○	◇	◇	◇	◇				
	44		◇	◇	○	○	○	◇	◇	◇	◇			
	46		◇	◇	○	○	○	○	◇	◇	◇	◇		
	48						◇	◇	◇	◇	◇	◇	◇	
	50							◇	◇	◇	◇	◇	◇	
	56											◇	◇	
60												◇	◇	
PLA Size		A			B				C					

Delivery times *)

○	5 days
◇	15 days
■	not available



*) applies to all non-positioning cushions; all cushions with positioning kit have 15 days delivery time.