

Invoice to:	
Name:	
Street:	
Town:	Postal code:
Country:	
Tel:	Fax:
Order date:	
Purchase order:	

Deliver to:	
Name:	
Street:	
Town:	
Country:	
Tel:	
Marked for:	

Standard Features	<ul style="list-style-type: none"> ● Angle adj. footplate or central legrest ● VR2 90Ah controller ● Parallel swing away mount ● 6KPH (4-pole) motors 	<ul style="list-style-type: none"> ● Solid drive wheels ● Solid castors ● 60Ah batteries ● Std curved headrest
<ul style="list-style-type: none"> ● JAY Comfort shallow contour backrest ● JAY Comfort seat cushion ● Matt Black Base & Frame (Jive M) 		

Fax: +44 (0)845 605 66 89

E-mail: orders@sunmed.co.uk

Spare Parts: <http://parts.sunrisemedical.eu>

JIVE R2: £6050

JIVE F: £7150

JIVE M: £7480



Max. user weight: 160kg (25 stone)

SEAT DIMENSION JAY COMFORT

Seat frame width is 41 cm. The Seat width can be adjusted between 40 and 56 cm via the armrests, based on the seat/back cushion configuration. The Seat depth can be adjusted between 41 and 55 cm via backrest movement

Seat dimensions also reflect aluminium seat plate option. Seat plate option to be selected at the Jay Comfort Cushion section (JIV020109).

Cushion width

		inches		15.5	16.5	17	18	19	20	20.5	21	22
		cm		40	42	44	46	48	50	52	54	56
Seat depth setting	48 cm depth seat cushion	16	41	○	○	○	○	○	○	○	○	○
				<i>JIV020001</i>	<i>JIV020021</i>	<i>JIV020005</i>	<i>JIV020025</i>	<i>JIV020009</i>	<i>JIV020029</i>	<i>JIV020013</i>	<i>JIV020033</i>	<i>JIV020017</i>
	55 cm depth seat cushion	18	46	○	○	○	○	○	○	○	○	○
				<i>JIV020002</i>	<i>JIV020022</i>	<i>JIV020006</i>	<i>JIV020026</i>	<i>JIV020010</i>	<i>JIV020030</i>	<i>JIV020014</i>	<i>JIV020034</i>	<i>JIV020018</i>
		20	51	○	○	○	○	○	○	○	○	○
				<i>JIV020003</i>	<i>JIV020023</i>	<i>JIV020007</i>	<i>JIV020027</i>	<i>JIV020011</i>	<i>JIV020031</i>	<i>JIV020015</i>	<i>JIV020035</i>	<i>JIV020019</i>
	21.5	55	○	○	○	○	○	○	○	○	○	
			<i>JIV020004</i>	<i>JIV020024</i>	<i>JIV020008</i>	<i>JIV020028</i>	<i>JIV020012</i>	<i>JIV020032</i>	<i>JIV020016</i>	<i>JIV020036</i>	<i>JIV020020</i>	

SEAT HEIGHT & DRIVE WHEEL POSITION

Seat height

inches	16.5	17	17.5	18	18.5	19
cm	42	43	44.5	46	47	48
RWD		○		○(3)		
<i>JIV010011</i>		<i>JIVH43</i>		<i>JIVH46</i>		
FWD		○(1)	○	○	○	○
<i>JIV010002</i>		<i>JIVH43</i>	<i>JIVH44,5</i>	<i>JIVH46</i>	<i>JIVH47</i>	<i>JIVH48</i>
MWD	○ (2)	○	○	○	○	○
<i>JIV010005</i>	<i>JIVH42</i>	<i>JIVH43</i>	<i>JIVH44,5</i>	<i>JIVH46</i>	<i>JIVH47</i>	<i>JIVH48</i>

Seat height measured floor to seat pan at 0° tilt

When powered lift/tilt options are selected achievable seat height is as follows:

- FWD: Tilt, only: 44.5cm, Lift & Tilt: 46cm
- MWD: Tilt: 43cm, Lift & Tilt: 44.5cm
- RWD: Tilt only: 43cm, Lift & Tilt: 44.5cm

JIV090022 ○ Seat height extension - Adds 3cm to seat height selected

(3) 30mm seat height kit to be selected (JIV090022)

(1) not available with power modules or 0° seat angle setting with no module

(2) not available with filler module in 0° and 3° seat angle setting (JIV030001)

SEAT TO BASE INTERFACE

JIV030001	<input type="radio"/>	Fixed tilt (filler module)	(tool adjustable seat angle)	
	<input type="radio"/>	0° - (Jive F / Jive M only)	<input type="radio"/>	3°
	<input type="radio"/>		<input type="radio"/>	6°
	<input type="radio"/>		<input type="radio"/>	9°
Power seating modules				
JIV030002	<input type="checkbox"/>	30° Power COG tilt		651
JIV030003	<input type="checkbox"/>	50° Power COG tilt	(Jive MWD Only)	694
JIV030004	<input type="checkbox"/>	-5° to 25° Powered COG tilt	(3cm added seat to floor height at 0°)	765
JIV030006	<input type="checkbox"/>	Power seat lift and tilt combined	(30° tilt and 30cm lift)	1748

JAY CUSHION SIZE**SEAT WIDTH**

Depth Coding	cm	Width Coding						
		JL40	JL42	JL44	JL46	JL48	JL50	JL56
JL__40	40	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
JL__42	42	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
JL__44	44	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
JL__46	46	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
JL__50	50	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
JL__56	56	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>

Size Configuration - COMPLETE PRODUCT**Conversion Table**

cm	40.0	42.0	44.0	46.0	48.0	50.0	56.0
inches	15.6	16.4	17.2	18.0	18.8	19.5	21.9

(1) Not available with Soft Combi P or Easy Fluid Cushion

(2) Not available with JAY Lite standard version

JAY CUSHION OPTIONS

JAY020004	<input type="checkbox"/>	Jay Lite Standard version	(Micro climat cover, Aluminium seat plate to be selected)	215
JAY020005	<input type="checkbox"/>	Jay Easy Fluid with Flat Base	(Incontinence cover, Aluminium seat plate to be selected)	308
JAY020001	<input type="checkbox"/>	Jay Easy Visco with Flat Base	(Incontinence cover, Aluminium seat plate to be selected)	132
JAY020003	<input type="checkbox"/>	Jay Soft Combi P	(Incontinence cover, Aluminium seat plate to be selected)	77

JAY COMFORT SEAT CUSHION OPTIONS

JIV020109	<input type="checkbox"/>	Aluminium seat pan - no cushion required	(Includes omit cushion option JIV020209)	-117
JIV020101	<input type="radio"/>	Comfort seat pan, without threshold	(thickness 2cm)	
JIV020102	<input type="radio"/>	Comfort seat pan, with 5cm (2") threshold	(thickness 2cm)	
JIV020201	<input type="radio"/>	Standard cushion - stretched rubber upholstery	(thickness 8cm)	
JIV020202	<input type="radio"/>	Optimum comfort - stretched rubber upholstery	(anatomic formed) (thickness 9cm)	

JAY COMFORT

JIV020101



JIV020102



JIV020201



JIV020202

JAY COMFORT BACK CUSHION OPTION

(57cm back height - Adj. up to 70cm)

JIV040001	<input type="radio"/>	Shallow contour (standard velour)	(hard back)	
JIV040002	<input type="checkbox"/>	Shallow contour (standard velour)	(vertical elastic straps)	111
JIV040003	<input type="checkbox"/>	Shallow contour (standard velour)	(horizontal tension adjustable straps)	111
JIV040004	<input type="checkbox"/>	Deep contour - active	(horizontal tension adjustable straps)	146
JIV040005	<input type="checkbox"/>	Deep contour - passive	(horizontal tension adjustable straps)	146
JIV040009	<input type="checkbox"/>	Omit back upholstery	(Chair supplied with NO back upholstery option)	-59

JAY COMFORTShallow Contour
JIV040001 /
JIV040002 /
JIV040003Deep Contour Active
JIV040004Deep Contour
Passive
JIV040005**COMFORT BACK CUSHION WIDTH**

JIV040101	<input type="radio"/>	40cm (15.5")	JIV040106	<input type="radio"/>	42cm (16.5")	JIV040102	<input type="radio"/>	44cm (17")
JIV040107	<input type="radio"/>	46cm (18")	JIV040103	<input type="radio"/>	48cm (19")	JIV040108	<input type="radio"/>	50cm (20")
JIV040104	<input type="radio"/>	52cm (20.5")	JIV040109	<input type="radio"/>	54cm (21")	JIV040105	<input type="radio"/>	56cm (22")

LATERAL SUPPORTS

JIV030601	<input type="checkbox"/>	Lateral supports short (pair)	266
JIV030602	<input type="checkbox"/>	Lateral supports long (pair)	266

THIRD PARTY BACKRESTS & INTERFACES

JIV010302	<input type="checkbox"/>	Jay J3 two point fixation kit	111
-----------	--------------------------	-------------------------------	-----

See Pages 8 and 9 to select J3 Backrest and options. Additional options can be found on the J3 prescription form

JAY J3 BACKRESTS AND OPTIONS

See Pages 8 and 9 to select J3 Backrest and options. Additional options can be found on the J3 prescription form

POSITIONING STRAPS

JIV020901	<input type="radio"/>	5cm (2") aircraft buckle	
JIV090056	<input type="radio"/>	5cm (2") aircraft buckle - 190cm (Long)	
JIV020902	<input type="checkbox"/>	5cm (2") padded aircraft buckle	52

HEADRESTS

JIV040501	<input type="radio"/>	Curved headrest	20cm x 11cm (8" x 4.5")	
JIV040502	<input type="radio"/>	Curved headrest	23cm x 11cm (9" x 4.5")	
JIV040503	<input type="checkbox"/>	Contoured (winged) headrest		87
JIV040504	<input type="checkbox"/>	Jay Comfort: standard comfort seating headrest		100
JIV040505	<input type="checkbox"/>	Jay Comfort: flexible comfort seating headrest (winged style)		122



Curved Headrest
JIV040501 /
JIV040502



Contoured Headrest
JIV040503



JAY Comfort: Std
Comfort headrest
JIV040504



JAY Comfort:
Winged headrest
JIV040505

BACKREST RECLINE OPTIONS

JIV040601	<input type="radio"/>	Manual adjustable backrest	(-5° to 12°, preset angle at 5°)	
JIV040602	<input type="checkbox"/>	Manual recline gas spring system	(-10° to 30°, require push handle option)	228
JIV040603	<input type="checkbox"/>	Power recline	(-10° to 40°) - Not available with VR2 control system	409
JIV040604	<input type="checkbox"/>	Power recline c/w mechanical shear reduction	(-10° to 40°) - Not available with VR2 control system	759
JIV090503	<input type="checkbox"/>	Height adjustable push handles		100

ARMREST TYPE

JIV050001	<input type="radio"/>	Single post armrest	(skirt guard to be selected)	
JIV050002	<input type="checkbox"/>	Reclining and flip-back armrest - High	(width and 32 - 38 cm height adjustable)	344
JIV040141	<input type="checkbox"/>	Reclining and flip-back armrest - Low	(width and 27 - 32 cm height adjustable)	344

SKIRT GUARDS

(not available with reclining armrest)

JIV050102	<input type="radio"/>	Long	36cm (14")	(plexi glass)	
JIV050103	<input type="checkbox"/>	Jay Comfort padded skirt guard	36cm (14")		52

ARMPADS - LEFT AND RIGHT HAND SIDE









	Standard arm pads		Left Right				
	Full length	36cm (14")	JIV050201	<input type="radio"/>	<input type="radio"/>	JIV050202	
Standard Full	Comfort arm pads	41cm (16")	JIV050203	<input type="radio"/>	<input type="radio"/>	JIV050204	
Standard Comfort	JAY COMFORT arm pads						
JAY Comfort Padded	Padded (vinyl) upholstery	full length 41cm (16")	JIV050251	<input type="checkbox"/>	<input type="checkbox"/>	JIV050252	52 each
	P-Shaped skai upholstery	full length 45 cm (18")	JIV050253	<input type="checkbox"/>	<input type="checkbox"/>	JIV050254	52 each
JAY Comfort P-Shaped	Hemi Armpads (Hemi Armpad & hand pad must be chosen together)						
	Hemi forearm pad	small 23cm (9")	JIV050205	<input type="checkbox"/>	<input type="checkbox"/>	JIV050206	87 each
	Hemi forearm pad	medium 28cm (11")	JIV050207	<input type="checkbox"/>	<input type="checkbox"/>	JIV050208	87 each
Hemi Armpads	Hemi forearm pad	large 33cm (13")	JIV050209	<input type="checkbox"/>	<input type="checkbox"/>	JIV050210	87 each
	Flat hand pads	small 11 x 18cm (4.5" x 7")	JIV050281	<input type="checkbox"/>	<input type="checkbox"/>	JIV050282	87 each
	Flat hand pads	large 14 x 23cm (5.5" x 9")	JIV050283	<input type="checkbox"/>	<input type="checkbox"/>	JIV050284	87 each

CENTRE MOUNT LEGRESTS

(Inc footboard)

JIV060001	<input type="radio"/>	Fixed centre mount legrest	Short	20.4 cm -31.6 cm (8" - 12.5")	
JIV060002	<input type="radio"/>	Fixed centre mount legrest	Medium	37 cm - 47 cm (13.5" - 18.5")	
JIV060003	<input type="radio"/>	Fixed centre mount legrest	Long	42 cm - 55 cm (16.5" - 22")	
JIV060004	<input type="checkbox"/>	Powered centre mount legrest		(37 - 49cm fixed mechanical adjustment + 14cm via motor:)	1080

			£
SWING AWAY LEGRESTS			<i>(width adjustable legrest mount; 8 cm (3") each side)</i>
JIV060101	<input type="radio"/> 70° s/a legrest		
JIV060102	<input type="checkbox"/> 75° s/a legrest tapered	<i>(Only with footboard, reduces inner legrest width by 2")</i>	111
JIV060104	<input type="checkbox"/> Manual elevating adult	<i>(c/w standard foam calf pads)</i>	263
JIV060105	<input type="checkbox"/> Powered elevating	<i>(c/w standard foam calf pads)</i>	496
JIV060106	<input type="checkbox"/> Powered articulating	<i>(c/w standard foam calf pads)</i>	683
JIV060107	<input type="checkbox"/> Auto connectors for powered s/a legrests		111
JIV060109	<input type="checkbox"/> Skai-vinyl calf pads		111
Extension tubes for swing-away legrests:			<i>(per pair)</i>
JIV060201	<input type="radio"/> Short	33 - 41cm (13" - 16")	
JIV060202	<input type="radio"/> Medium	39 - 46cm (15" - 18")	
JIV060203	<input type="radio"/> Long	43 - 53cm (17" - 21")	
FOOTPLATES FOR SWING-AWAY LEGRESTS			<i>(per pair)</i>
JIV060302	<input type="radio"/> Aluminium angle adjustable		
JIV060303	<input type="checkbox"/> Foot board	<i>(angle adjustable)</i>	<i>(not available on 90° legrest)</i>
JIV060401	<input type="checkbox"/> Extra legrest strap		24
VR2 CONTROLS			<i>(Not Available with 13Kph Motors - Rnet must be selected to deliver increased Power for Kerb/Hill Climbing)</i>
JIV100002	<input type="radio"/> VR2 Joystick: with actuator and light control		<i>(lights & indicators plus max. 2 actuators)</i>
Dual/Attendant control options			
JIV100201	<input type="checkbox"/> VR2 dual control		438
JIV100202	<input type="checkbox"/> VR2 JSM attendant only		111
	VR2 User Joystick JIV010002		VR2 Dual Control Joystick JIV010201
R-NET 120A CONTROLS			<i>(Required for over 2 actuators & users requiring max torque & performance)</i>
JIV100102	<input type="checkbox"/> R-net control system		<i>(Must be selected if R-Net control system required)</i>
JIV100152	<input type="radio"/> R-net LED joystick		<i>(Incl. L&I and actuator control option)</i>
JIV100153	<input type="checkbox"/> R-net LCD colour joystick		<i>(Integrated colour display with 2 jack inputs for buddy buttons)</i>
JIV100154	<input type="radio"/> R-Net Control System - Omit joystick		<i>(When joystick is not required for clients using special controls)</i>
Dual/Attendant control options			
JIV100201	<input type="checkbox"/> Dual control		<i>(Separate dual control module mounted at rear - Always RH)</i>
JIV100202	<input type="checkbox"/> R-net attendant only joystick		111
	R-net LED User Joystick JIV100152		R-net LCD User Joystick JIV100153
	R-net Dual Control Joystick JIV100201		R-net Attendant Only Joystick JIV100202
JOYSTICK MOUNTINGS			
Front control mountings			<i>(for single Joystick selection mount)</i>
JIV100251	<input type="radio"/> Swing-away		
JIV100252	<input type="checkbox"/> Centre bar control mount		168
JIV100253	<input type="checkbox"/> Swing-away tray mounted centre control - R-Net LED		<i>(not available on 22" width)</i>
JIV110143	<input type="checkbox"/> Swing-away tray mounted centre control - R-Net LCD		<i>(not available on 22" width)</i>
Joystick mounting side			
JIV100271	Left <input type="radio"/>	Right <input type="radio"/>	JIV100272
LIGHTS & INDICATORS			
JIV110201	<input type="radio"/> LED (StVZO approved)		

								£
R-NET OMNI DISPLAY & SPECIAL CONTROL INPUT MODULE								(with IR functionality)
JIV100401	<input type="checkbox"/>	R-net OMNI colour display input module					(includes Special Control Interface)	876
JIV100550	<input type="checkbox"/>	OMNI in scanning mode						876
JIV100422	<input checked="" type="checkbox"/>	R-NET OMNI Flexi mount only					(30 cm (12") long flexible height extension)	
		JIV100431	Left	<input type="checkbox"/>	<input type="checkbox"/>	Right	JIV100432	
R-NET OMNI ACCESSORIES								
JIV100503	<input type="checkbox"/>	R-net Bluetooth module					(Compatible with Windows PC and Android systems only)	344
B4M	<input type="checkbox"/>	R-net Input/Output module					(Build position selected as B4M process)	578
JIV090516	<input type="checkbox"/>	Bluetooth wireless mouse receiver						65
 <p>R-net OMNI + Flexi Mount JIV100401</p> <p>Provides Visual Display, infrared commands (Controls TV, DVD, Stereo, etc) & Special Control interface</p>		 <p>R-net Bluetooth Module JIV100503</p> <p>Allows users to control Windows based PC's, Tablets and Phones from joystick (mounted in base)</p>						
R-NET SPECIAL CONTROLS								
(R-net OMNI required from previous page; 2 buddy buttons recom.)								
JIV100551	<input type="checkbox"/>	Switch-It Versaguide Joystick as hand control						1693
JIV100552	<input type="checkbox"/>	Switch-It Micro Guide as hand control						2276
JIV100553	<input type="checkbox"/>	Switch-It Micro Pilot (0 way) joystick as hand control						2509
JIV100602	<input type="checkbox"/>	Swing-away front mount						
JIV100651		Left	<input type="checkbox"/>	Right	<input type="checkbox"/>	JIV100652		
 <p>R-Net OMNI in Scanning Mode JIV100550</p>		 <p>Proportional Versaguide Joystick JIV100551</p>		 <p>Switch-It Micro Guide Joystick JIV100552</p>		 <p>Switch-It Micro Pilot 0-Way Joystick JIV100553</p>		
R-NET HEAD RELATED INPUT DEVICES								
(if you have already chosen a standard headrest then Switch- It head array is not possible)								
JIV100701	<input type="checkbox"/>	Switch-It 3 switch head array / adult					(5 switch and other versions available via B4M) (OMNI required, 2 BB recommended)	3443
JIV100703	<input type="checkbox"/>	Switch-It Micro Pilot (0 way) joystick as chin control					(OMNI required, 2 BB recommended)	2509
JIV100704	<input type="checkbox"/>	Proportional chin control					(proportional joystick; OMNI required; 2 BB recommended)	876
JIV100705	<input type="checkbox"/>	Sip & Puff control					(tube and mouthpiece; OMNI required; 2 BB recommended)	293
JIV100706	<input type="checkbox"/>	Switch-It Micro Guide as chin control					(OMNI required, 2 BB recommended)	2276
JIV100707	<input type="checkbox"/>	Switch-It proportional head control					(OMNI required, 2 BB recommended)	3210
JIV100708	<input type="checkbox"/>	Switch-It Versaguide Joystick as chin control					(OMNI required, Switch-It flex switch kit incl.)	1810
 <p>Switch-It Head Array JIV100701</p>		 <p>Proportional Versaguide Joystick as chin control JIV100708</p>						
CHIN CONTROL SWING AWAY OPTIONS								
JIV100801	<input type="checkbox"/>	Manual s/a option						228
	<input type="checkbox"/>	Powered s/a option					(choose your setting below)	
JIV100802	<input type="checkbox"/>	Operated with buddy button					(operated via one buddy button - select colour below)	1693
JIV100803	<input type="checkbox"/>	Operated with ribbon switch						1810
		Left	<input type="checkbox"/>	Right	<input type="checkbox"/>	Swing-away direction		
JIV100851		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	JIV100852		
Buddy button switches								
JIV100901	<input type="checkbox"/>	Red x 1						87
Piko button switches								
JIV100909	<input type="checkbox"/>	Red	x 1	JIV100913	<input type="checkbox"/>	Red	x 2	87 each
JIV100910	<input type="checkbox"/>	Green	x 1	JIV100914	<input type="checkbox"/>	Green	x 2	87 each
JIV100911	<input type="checkbox"/>	Yellow	x 1	JIV100915	<input type="checkbox"/>	Yellow	x 2	87 each
JIV100912	<input type="checkbox"/>	Blue	x 1	JIV100916	<input type="checkbox"/>	Blue	x 2	87 each

● = Standard ○ = No cost option □ = Cost option			£
REAR or FRONT WHEEL DRIVE			
CASTOR WHEELS <i>(delivered always as a pair)</i> <i>(JIVE R/F only)</i>			
JIV080001	○	9" x 2.5" Solid	
JIV080002	○	9" x 2.5" Pneumatic	
JIV080003	□	10" x 3" Solid <i>(Not available with 13 kph motors and 6 kph FWD / RWD. Not recommended on RWD base)</i>	87
JIV080004	□	10" x 3" Pneumatic <i>(Not recommended on RWD base)</i>	87
MID WHEEL DRIVE			
JIV080101	●	6" solid castor wheel on front and rear	
KERB CLIMBER			
JIV080501	□	Kerb climber <i>(Jive R2 only)</i>	168
DRIVE BASE <i>(Max kerb climbing is affected by many factors incl user weight, speed selection etc. Contact Tech Serv Ctr for more info)</i>			
MOTORS/SPEED			
JIV090002	○	6kph (4mph) (4-pole) high torque motors for optimal kerb climbing	
JIV090003	○	10kph (6mph) (4-pole) high torque motors for optimal kerb climbing	
JIV090004	□	13kph (8mph) (4-pole) <i>(Not available on JIVE F)</i> <i>(Not available with VR2 - Rnet must be selected)</i>	344
DRIVE WHEELS			
JIV090101	○	14" Pneumatic "V-type" <i>(Recommended for better kerb climbing)</i>	
JIV090102	○	14" Solid "V-type"	
JIV090103	○	14" Pneumatic "knobby" <i>(Recommended for better kerb climbing)</i>	
JIV090104	○	14" Solid "Knobby"	
SUSPENSION			
JIV010901	○	Performance Spring Sports suspension, <i>(Greater traction & more responsive ride. Reduced kerb climbing to 7cm max)</i>	
JIV010909	○	Omit Performance Spring - Comfort suspension (1) <i>(Traditional suspension - greater comfort and kerb climbing - Max 10cm)</i>	
BATTERIES			
JIV090201	○	60 Ah maintenance free	
JIV160005	□	80 Ah maintenance free <i>(Recommended for 13Kph Motors where battery power is drawn faster)</i>	111
CHARGER			
JIV160100	○	220V/ 8A automatic charger (UK)	
JIV090303	□	220V/ 10A automatic charger (UK) <i>(Must be ordered with 13kph speed - charges batteries faster)</i>	168
JIVE F/ JIVE R2		JIVE MWD	
(body panel colour		(body panel colour	
JIV010101	○	White RAL 9010	JIV140004 ○ Matt Black RAL 9005
JIV010102	○	Red	JIV010101 ○ White RAL 9010
JIV010103	○	Blue	JIV010102 ○ Red
			JIV010103 ○ Blue

£

ACCESSORIES

JIV090501	<input type="checkbox"/> Crutch holder	30
JIV090502	<input type="radio"/> Rear view mirror <i>(Always required with 10Kph Chairs)</i>	
JIV090512	<input type="checkbox"/> Swing-away tray table <i>(available from 36 cm (14") to 51 cm (20") seat width)</i>	168
JIV090504	<input type="checkbox"/> 24/12V DC/ DC converter	228

JOYSTICK KNOBS

JIV090505	<input type="checkbox"/> Golf ball joystick knob	30
JIV090506	<input type="checkbox"/> T bar B38 joystick knob	30
JIV090507	<input type="checkbox"/> Soft ball joystick knob	30
JIV090508	<input type="checkbox"/> Ball joystick knob	30
JIV090509	<input type="checkbox"/> Chin cup joystick knob	30
JIV090510	<input type="checkbox"/> Stick joystick knob	30
JIV090511	<input type="checkbox"/> Mushroom joystick knob	30

TRANSPORTATION OPTIONS

JIV090539	<input checked="" type="checkbox"/> Prepared for Dahl Docking Station option <i>(Jive F and M via B4M only. Jive R2 docking preparation as standard)</i>	
-----------	--	--

For special requirements, please complete Built-4-Me prescription form



JIVE - base price	£
Option Price	£
Total Price	£



J3 Back

J3 Backrest Range

Complete backrest package for clients up to 136kg
 Universal one point fixation for rapid application
 Easily adjustable for assessment & changing needs
 Wide range of positioning and transfer options
 Backs and options fully transit tested
 1 Large & 1 small lumbar included in MT, UT & SH box
 1 small lumbar in LT box 2 year guarantee



- = 5 days from the EDC
- ◇ = 15 days from the EDC
- = upcharge option

Max. user weight: **136kg**
 USER WEIGHT kg
Mandatory information to ensure correct configuration


Applications:

For all clients requiring postural back support from mild lateral stability to extreme lateral and posterior stability.
 Full range of postural back shapes to cover clinical back support at the appropriate backrest height for the user.
 Adjustable in depth, height, angle and chair width, ensuring rapid and comprehensive client positioning and spine alignment.


CLIENT MEASUREMENT

1. Assess Lateral Stability Requirement (SC, MC, DC or PDC)
2. Assess Posterior Stability Requirement (MT, UT or SH)
3. Assess appropriate backrest height (Short, Medium, Tall)
4. Determine client trunk width across the sternum; and cross check with backpost width.
5. Assess for options and 'Velcro' positioning inserts


Item Nr. construction:
J3 SC LT S 41 (This is only an example)
J3: J3 Back
SC: Shallow Contour
LT: Lower Thoracic
S: Short (Height cm)
41: Chair Backrest Width (cm)




Shallow Contour / SC
Moderate Lateral Support



Mid Contour / MC
Medium Lateral Support



Deep Contour / DC
Full Lateral Support



Posterior Deep Contour / PDC
PDC Full Lateral support

Width

Chair backrest widths (Measured from outside of backrest to outside of backrest)

Jay 3 Backrest Width

	41	46	51		41	46	51		41	46	51		41	46	51
	41 - 46	46 - 51	51 - 56		41 - 46	46 - 51	51 - 56		41 - 46	46 - 51	51 - 56		41 - 46	46 - 51	51 - 56


Back height

Laterals, headrest & harnesses available



Client max trunk width (cm)	Short (S) - 42
Upper Thoracic (UT)	Med (M) - 46
	Tall (T) - 50

	£520			£520			£570			£570		
	38	43	48	34,5	39,5 (1)	44,5	33	38 (1)	43	33	38	43
	○	○	○	○	○	○	○	○	○	○	○	○
	○	○	○	○	○	○	○	○	○	○	○	○
	○	○	○	○	○	○	○	○	○	○	○	○



Client max trunk width (cm)	Short (S) - 53
Shoulder Height (SH)	Med (M) - 57
	Tall (T) - 61

	£520			£520			£570			£570		
	38	43	48	34,5	39,5 (1)	44,5	33	38 (1)	43	33	38	43
	○	○	○	○	○	○	○	○	○	○	○	○
	○	○	○	○	○	○	○	○	○	○	○	○
	○	○	○	○	○	○	○	○	○	○	○	○

(1) Not available with laterals

Conversion Table

cm	17.0	20.0	24.0	30.0	34.0	36.0	38.0	41.0	42.0	44.0	46.0	51.0	53.0	56.0	61.0
inches	6.5	8.0	9.5	12.0	13.5	14.0	15.0	16.0	16.5	17.5	18.0	20.0	21.0	22.0	24.0

JAY© 3 BACK ACCESSORIES

= upcharge option

ACCESSORIES		Unit Price
JAY090001	<input type="checkbox"/> Backrest depth extension	£39
JAY090002	<input type="checkbox"/> 4 Point hardware kit	NEW! £162

HEADREST (only available where headrest is indicated on the first page)	Crash tested	Unit Price
---	--------------	------------

Please choose 1 mount and 1 pad to have a complete headrest with pad

Mount		Unit Price
JAY090006	<input type="checkbox"/> Jay standard headrest (bracket & mount only)	£88
JAY090007	<input type="checkbox"/> Jay flip back headrest (bracket & mount only)	£198

Pad		Unit Price
JAY090008	<input type="checkbox"/> Jay standard pad	£41
JAY090009	<input type="checkbox"/> Jay contour pad	£41

JAY090010	<input type="checkbox"/> Universal headrest mount (required to use any other headrest not listed in this order form)	£41
-----------	--	-----

LATERAL SUPPORT (only available where laterals are indicated on the first page)	Unit Price
---	------------

Left Hand lateral		Depth					
cm		10	12	15			
Height	10	<input type="checkbox"/> JAY090011	<input type="checkbox"/> JAY090012	<input type="checkbox"/> JAY090013	←	Fix lateral JAY090017	£46
	15	<input type="checkbox"/> JAY090014	<input type="checkbox"/> JAY090015	<input type="checkbox"/> JAY090016	←	Fix lateral JAY090018	£46
		<input type="checkbox"/> JAY090019	←	Swing-away option		£87	

Right Hand lateral		Depth					
cm		10	12	15			
Height	10	<input type="checkbox"/> JAY090020	<input type="checkbox"/> JAY090021	<input type="checkbox"/> JAY090022	←	Fix lateral JAY090026	£46
	15	<input type="checkbox"/> JAY090023	<input type="checkbox"/> JAY090024	<input type="checkbox"/> JAY090025	←	Fix lateral JAY090027	£46
		<input type="checkbox"/> JAY090028	←	Swing-away option		£87	

When ordering 2 laterals Left and Right, you need to pick the sizes in both matrices above.

Total Price for laterals

If a swing-away lateral is required, first select the height and the depth, then tick the swing-away option box below the matrix

POSITIONING KIT		CM width	41 cm	46 cm	51 cm	Unit Price
Spine Align Mini Fit Kit I			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£31
			JAY090029	JAY090030	JAY090031	
Spine Align Mini Fit Kit II			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£47
			JAY090032	JAY090033	JAY090034	
JAY090035	<input type="checkbox"/> Jay 3 Back Positioning Box					£151
JAY090036	<input type="checkbox"/> Jay 3 Back Hardware Box					£474

CHEST STRAP (only available where harnesses are indicated on the first page)		Medium	Large	Extra large	Unit Price
Padded anterior strap		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£157
		JAY090037	JAY090038	JAY090039	
Centre open harness		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£157
		JAY090040	JAY090041	JAY090042	

SPARES (separate to product - separate box or delivery)		CM width	41 cm	46 cm	51 cm	Unit Price
Small lumbar <small>(1 unit included as std with product LT, MT, UT & SH)</small>			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£9.00
			JAY090043	JAY090044	JAY090045	
Large lumbar support <small>(1 unit included as std with product MT, UT & SH)</small>			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£9.00
			JAY090046	JAY090047	JAY090048	
Wedge			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	£9.00
			JAY090049	JAY090050	JAY090051	

		Unit price	Quantity	Total Price
JAY090052	<input type="checkbox"/> Curved shim	£9.00		
JAY090053	<input type="checkbox"/> Flat shim	£9.00		
JAY090054	<input type="checkbox"/> Posterior/lateral pelvic	£12.00		
JAY090055	<input type="checkbox"/> Lateral shim	£4.00		
JAY090056	<input type="checkbox"/> Small bean bag	£9.00		
JAY090057	<input type="checkbox"/> Large bean bag	£9.00		