

Invoice to:  
 Name:  
 Street:  
 Town: Postal code:  
 E-Mail:  
 Tel: Fax:  
 Order date:  
 Purchase order:

Deliver to:  
 Name:  
 Street:  
 Town: Postal code:  
 E-Mail:  
 Tel: Fax:  
 Marked for:

Fax number: +44 (0)845 605 66 89  
 Email: orders@sunmed.co.uk

USER WEIGHT  kg  
 Crash tested and passed without tie down brackets

PLEASE CHOOSE ONE OF THE FOLLOWING FRAME OPTIONS

**Helium Generation 2013**  
 Maximum user weight 125kg

Base price £2525

HUM030017

**Helium Pro**

Maximum user weight 125kg

**Fixed welded axle position**

(no rear seat height adjustment & +/- 1cm centre of gravity adjustment)

**Standard angle adjustable back**

£369

HUM030016

**Helium Pro**

Maximum user weight 100kg

**Fixed welded axle position**

(no rear seat height adjustment & +/- 1cm centre of gravity adjustment)

**Fixed welded backrest set at\***

- 4° (86°)
- 2° (88°)
- 0° (90°)       £493
- +2° (92°)      Figures in brackets are angle to ground
- +4° (94°)

\*Must choose Oval Shaped Backrest

CLIENT MEASUREMENT

Seat width in cm (measured from outside frame tube to outside frame tube)


		32	34	36	38	40	42	44	46
Seat depth in cm	34	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	36	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	38	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	40	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	42	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	44	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	46	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	48	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

- Please note the procedure of the measurement is from the back  
 - All measurements are in cm and without a cushion

High frame is 4cm higher than low frame and you can choose shorter and lighter forks to minimise the weight of your wheelchair. See next page

		Front seat height in cm																	
		Fork (Length x Width)	Frame	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	
Front Castor	3"	111 mm	Low	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>													
		111 mm	High					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>									
	4"	111 mm (1)	Low		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>											
		123 mm (2)			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>											
		Frog Leg (2)						<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>									
		111 mm (1)	High						<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>							
		123 mm (2)								<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
		Frog Leg (2)										<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					
	5"	111 mm (1)	Low				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>										
		123 mm (2)						<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>									
		Frog Leg (2)									<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>						
		111 mm (1)	High								<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					
		123 mm (2)										<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
		Frog Leg (2)											<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>
	6"	123 mm (2)	Low					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>									
		Frog Leg (2)											<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				
		123 mm (2)	High											<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			<input type="checkbox"/>


Your rear seat height can only be between 0cm and 13cm below your front seat height:  
 i.e. rear seat height CANNOT be above front seat height.  
 i.e. rear seat height CANNOT be more than 13cm below front seat height.  
 - More info, see page 2  
 # 4" & 5" soft roll and 5" & 6" soft solid are only available with a 45 mm fork



(1): Also available as one arm fork; please choose the fork type on page 4  
 (2): Compatible with Power-Add on devices

REAR WHEELS															
Rear wheel	Rear seat height in cm <span style="float: right;">(The ordered rear seat height is adjustable +/- 1cm on each chair. Excludes Helium Pro.)</span>														
	Wheels	37*	38*	39	40	41	42	43	44	45	46	47	48	49*	50*
	24"	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>
25"	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	<input type="radio"/>	
*With std seat upholstery decreasing rear seat height not possible in 37cm (24" wheel) and 38cm (25" wheel).															
*With std seat upholstery increasing rear seat height not possible in 49cm (24" wheel) and 50cm (25" wheel).														£	
SEAT UPHOLSTERY															
<b>Seat sling standard</b>															
HUM020002	<input type="radio"/>	With 1 utility bag												(1)	
HUM020003	<input type="checkbox"/>	With 2 utility bags												(1)	(only from SW 36cm) 47
HUM020004	<input type="checkbox"/>	With 1 utility bag and 1 catheter bag												(1)	(only from SW 36cm) 75
HUM020010	<input type="checkbox"/>	<b>Seat Sling Clip Style - Airflow (user weight limit 100kg)</b>												(1)	178
HUM020005	<input type="checkbox"/>	With 1 utility bag												(1)	221
HUM020006	<input type="checkbox"/>	With 2 utility bags												(1)	(only from SW 36cm Max 100Kg) 250
HUM020007	<input type="checkbox"/>	With 1 utility bag and 1 catheter bag												(1)	(only from SW 36cm Max 100Kg) 279
HUM020000	<input type="radio"/>	Without seat sling													
<b>Jay Cushions - for complete product range see Jay order forms</b>															
JAY020002	<input type="checkbox"/>	JAY Basic													55
JAY020003	<input type="checkbox"/>	JAY Soft Combi P													77
JAY020001	<input type="checkbox"/>	JAY Easy Visco, Curved Base													132
JAY020007	<input type="checkbox"/>	JAY Easy Fluid, Curved Base													308
JAY020004	<input type="checkbox"/>	JAY Lite													215
<b>Seat cushion black cover composite foam</b>															
a) Foam composite		Thickness:		HUM020101		<input type="checkbox"/> 3cm		HUM020102		<input type="checkbox"/> 5cm		(medium)		61	
b) Latex		Thickness:		HUM020103		<input type="checkbox"/> 5cm		HUM020104		<input type="checkbox"/> 8cm		(very soft)		61	
<b>Seat cushion with water resistant / breathable cover</b>															
a) Foam composite		Thickness:		HUM020105		<input type="checkbox"/> 3cm		HUM020106		<input type="checkbox"/> 5cm		(medium)		114	
b) Latex		Thickness:		HUM020107		<input type="checkbox"/> 5cm		HUM020108		<input type="checkbox"/> 8cm		(very soft)		114	
<b>Seat accessories</b>															
HUM090018	<input type="checkbox"/>	Positioning belt with lock													39
HUM020150	<input type="checkbox"/>	Spare cushion cover with zipper													41
BACKREST															
HUM030001	<input type="radio"/>	<b>Standard backrest</b>			25 - 30		30 - 35		35 - 40		40 - 45				
<b>Back height factory setting at (cm)</b>				<input type="radio"/> 25		<input type="radio"/> 30		<input type="radio"/> 35		<input type="radio"/> 40					
				<input type="radio"/> 27.5		<input type="radio"/> 33		<input type="radio"/> 37.5		<input type="radio"/> 42.5					
				<input type="radio"/> 30		<input type="radio"/> 35		<input type="radio"/> 40		<input type="radio"/> 45					
HUM030003	<input type="checkbox"/>	<b>Oval shaped backrest (Non Height adjustable and please pick bar setting below)</b>												215	
		<input type="radio"/> 25	<input type="radio"/> 27.5	<input type="radio"/> 30	<input type="radio"/> 32.5	<input type="radio"/> 35	<input type="radio"/> 37.5	<input type="radio"/> 40	<input type="radio"/> 42.5	<input type="radio"/> 45					
HUM030403	<input type="radio"/>	Stabilizer bar setting close (9.2cm)				HUM030404 <input type="radio"/> Stabilizer bar setting wide (12.4cm)				(distance away from the backrest)					
HUM030015	<input type="radio"/>	Back angle adjustable foldable												(-31° to +15°; only in black anodizing)	
HUM030012	<input type="radio"/>	Angle adjustable back (non folding)												(-31° to +15°; only in black anodizing)	
HUM030102	<input type="checkbox"/>	Double locking (angle adjustable and folding, locks in folded position)												(-31° to +15° anodised colour options on last page) 50	
HUM030103	<input type="checkbox"/>	Backrest width ONE size smaller than standard (-2cm less than the seat width)												(not available with Freestyle backrest) 142	
PUSH HANDLES															
HUM030200	<input type="radio"/>	Omit push handles													
HUM030201	<input type="radio"/>	Long												(2)	
HUM030203	<input type="checkbox"/>	Fold-down push handles												(2) (only available from 32.5cm standard back) 119	
HUM030211	<input type="checkbox"/>	Removable push handles												(only with Freestyle back) 158	
HUM030204	<input type="checkbox"/>	Height adjustable push handles												(not possible with OSB (4067) & from 32.5cm on the std backrest) 158	
HUM030205	<input type="checkbox"/>	Height adjustable push handles for Jay backs												158	
(1):	<b>no change from standard to clipstyle seat sling possible (and other way round) in the field afterwards</b>														
(1):	<b>not available with height adjustable push handles</b>														

		£
<b>BACK UPHOLSTERY</b>		
HUM030316	<input type="radio"/> EXO backrest sling, synthetic fibre, padded, adjustable, 1 accessory bag black	
Colour binding tape: <input type="radio"/> Black <input type="radio"/> Orange <input type="radio"/> Silver <input type="radio"/> Red <input type="radio"/> Blue		
HUM030317	<input type="checkbox"/> EXO PRO backrest sling, synthetic fibre, padded, adjustable, breathable, spacer fabric, 1 accessory bag, black	136
Colour binding tape: <input type="radio"/> Black <input type="radio"/> Orange <input type="radio"/> Silver <input type="radio"/> Red <input type="radio"/> Blue		
HUM030300	<input type="radio"/> Without backrest upholstery	
<b>FREESTYLE BACK, ANGLE ADJUSTABLE, FOLDABLE (not crash tested)</b>		
HUM030089	<input type="checkbox"/> <b>Freestyle back, J3 Aluminium backrest shell, shallow contour (SC)</b> <span style="float: right;"><i>(in frame colour)</i></span>	550
<input type="radio"/> Backrest shell height 31cm <span style="float: right;"><i>(backheight 325 - 425 mm)</i></span> <input type="radio"/> Backrest shell width 32cm - can be used from SW 32cm <input type="radio"/> Backrest shell width 37cm - can be used from SW 38cm		
HUM030090	<input type="checkbox"/> <b>Freestyle back, J3 Carbon backrest shell, shallow contour (SC)</b> <span style="float: right;"><i>(backheight 275 - 350 mm)</i></span>	660
<input type="radio"/> Backrest shell height 24cm <span style="float: right;"><i>(backheight 275 - 350 mm)</i></span> <input type="radio"/> Backrest shell width 32cm - can be used from SW 32cm <input type="radio"/> Backrest shell width 37cm - can be used from SW 38cm <input type="radio"/> Backrest shell width 42cm - can be used from SW 44cm <input type="radio"/> Backrest shell height 31cm <span style="float: right;"><i>(backheight 325 - 425 mm)</i></span> <input type="radio"/> Backrest shell width 32cm - can be used from SW 32cm <input type="radio"/> Backrest shell width 37cm - can be used from SW 38cm <input type="radio"/> Backrest shell width 42cm - can be used from SW 44cm		
<b>JAY BACKRESTS</b>		
<b>Shallow Contour</b>		
<b>JAY3 Back - Shallow Contour</b>		
Mid Thoracic:	<input type="checkbox"/> Medium - 34 cm (JAY030010) <input type="checkbox"/> Tall - 38 cm (JAY030033)	506
Upper Thoracic:	<input type="checkbox"/> Short - 42 cm (JAY030012) <input type="checkbox"/> Medium - 46 cm (JAY030019) <input type="checkbox"/> Tall - 50 cm (JAY030025)	572
<b>JAY3 Back - Mid Contour</b>		
Mid Thoracic:	<input type="checkbox"/> Medium - 34 cm (JAY030011) <input type="checkbox"/> Tall - 38 cm (JAY030018)	506
Upper Thoracic:	<input type="checkbox"/> Short - 42 cm (JAY030013) <input type="checkbox"/> Medium - 46 cm (JAY030020) <input type="checkbox"/> Tall - 50 cm (JAY030026)	572
<b>JAY3 Carbon - Shallow Contour</b>		
Low Thoracic:	<input type="checkbox"/> Short - 17cm (JAY060001) <input type="checkbox"/> Tall - 24 cm (JAY060002)	605
Mid Thoracic:	<input type="checkbox"/> Short - 30cm (JAY060003) <input type="checkbox"/> Tall - 32 cm (JAY060004)	605
<b>JAY3 Accessory bag</b>		
JAY090058	<input type="checkbox"/> J3 Accessory Bag	50
<b>JAY Accessories</b>		
JAY030004	<input type="checkbox"/> Jay Flip Back Headrest (J3EUHDFH & J3EUHDLSP) <span style="float: right;"><i>(only with J3 backs)</i></span>	255
JAY030006	<input type="checkbox"/> Jay Std. Headrest (J3EUHDXH & J3EUHDLSP) <span style="float: right;"><i>(only with J3 backs)</i></span>	138
<b>SIDE GUARDS &amp; ARMRESTS</b>		
HUM040101	<input type="checkbox"/> Slim style carbon side guard <span style="float: right;"><i>(1) (slim style, carbon fibre)</i></span>	187
HUM040102	<input type="checkbox"/> Straight side guard aluminium (without fender)	163
HUM040103	<input type="checkbox"/> Fender side guard - aluminium <span style="float: right;"><i>(min seat height 39cm)</i></span>	224
HUM040104	<input type="checkbox"/> Fender side guard - carbon <span style="float: right;"><i>(min seat height 39cm)</i></span>	293
HUM040107	<input type="checkbox"/> Plastic detachable side guard <span style="float: right;"><i>(Quickie style)</i></span>	152
HUM040108	<input type="checkbox"/> Tubular armrest, Aluminium, swing away (sideways), padded, removable <b>(4)</b> <span style="float: right;"><i>(not in combination with handbike axle)</i></span>	188
Single post armrest with: <span style="float: right;"><i>(Quickie style; not with tubular armrest)</i></span>		
HUM040050	<input type="checkbox"/> short armpad (25cm) <input type="checkbox"/> long armpad (33cm)	163
HUM040052	<input type="checkbox"/> Single post, tool height adjustable armrest, depth adjustable (33cm) <b>(2)</b> <span style="float: right;"><i>(not with tubular armrest)</i></span>	187
HUM040000	<input type="radio"/> Without side guard	
<b>LEGREST &amp; FOOTPLATES</b>		
HUM050029	<input type="radio"/> Lightweight tubular footrest <span style="float: right;"><i>(lower leg length 36cm - 52cm)</i></span>	
HUM050030	<input type="checkbox"/> Tubular footrest lightweight with flat cover <span style="float: right;"><i>(lower leg length 36cm - 52cm)</i></span>	68
HUM050025	<input type="checkbox"/> Platform footplate lightweight angle & depth adjustable composite <span style="float: right;"><i>(lower leg length 36cm - 52cm)</i></span>	138
HUM050024	<input type="checkbox"/> Platform footplate lightweight angle & depth adjustable aluminium <span style="float: right;"><i>(lower leg length 36cm - 52cm)</i></span>	138
HUM050026	<input type="checkbox"/> Platform footplate lightweight angle & depth adjustable carbon	282
HUM050034	<input type="checkbox"/> Platform footplate lightweight angle & depth adjustable aluminium, high mounted <span style="float: right;"><i>(LLL 22 - 39cm)</i></span>	138
HUM050125	<input type="checkbox"/> Performance footplate angle adjustable <span style="float: right;"><i>(lower leg length between 36cm - 52cm)</i></span>	206
HUM050127	<input type="checkbox"/> Foot positioning plate <span style="float: right;"><i>(only with platform footplates lightweight angle &amp; depth adjustable)</i></span>	19
HUM050128	<input type="checkbox"/> Footboard, platform, plastic, angle and depth adjustable, flip-up (to the back) <b>(5)</b> <span style="float: right;"><i>(LLL from 33 cm; not with 6" castors)</i></span>	136
HUM050022	<input type="checkbox"/> Footboard, divided, aluminium, angle-adjust., flip-up (sideways), frame colour <b>(3)</b>	136
<b>LOWER LEG LENGTH (cm)</b> <span style="float: right;"><small>Attention: LLL has to be 3cm less than front seat height. All measurements are in cm</small></span>		
<input type="radio"/> 22 <input type="radio"/> 23 <input type="radio"/> 24 <input type="radio"/> 25 <input type="radio"/> 26 <input type="radio"/> 27 <input type="radio"/> 28 <input type="radio"/> 29 <input type="radio"/> 30 <input type="radio"/> 31 <input type="radio"/> 32 <input type="radio"/> 33 <input type="radio"/> 34 <input type="radio"/> 35 <input type="radio"/> 36 <input type="radio"/> 37 <input type="radio"/> 38 <input type="radio"/> 39 <input type="radio"/> 40 <input type="radio"/> 41 <input type="radio"/> 42 <input type="radio"/> 43 <input type="radio"/> 44 <input type="radio"/> 45 <input type="radio"/> 46 <input type="radio"/> 47 <input type="radio"/> 48 <input type="radio"/> 49 <input type="radio"/> 50 <input type="radio"/> 51 <input type="radio"/> 52		
Low Frame: LLL 22-35cm only with high mounted footrest; for all other footrests from LLL 36-50cm		
High Frame: LLL 22-39cm only with high mounted footrest; for all other footrests from LLL 39-52cm		
(1):	with 24" rear wheel minimum rear seat height 41cm; with 25" rear wheel minimum rear seat height 43cm	
(2):	not with centre of gravity 9.2 & 10.4cm in combination with lightweight seat upholstery	
(3):	Lower leg length from 38cm with low frame, from 42cm with high frame, not with 4" and 6" castor	
(4):	only with slim style sideguard, Aluminium sideguard w/o fender or w/o sideguard. Not with lightweight upholstery.	
(5):	user weight limitation 75kg, lower leg length from 33cm, not with 4" castor	
(6):	not with height adjustable push handles	

			£
<b>FRAME VERSIONS (first choose frame length, then frame angle and inset)</b>			
HUM010206	<input type="radio"/>	Standard frame	
HUM010207	<input type="radio"/>	Long frame <i>(not in combination with SD 46cm; recommended with narrow castor house position)</i>	
HUM010012	<input type="radio"/>	100° Fixed front frame (UK 80°) <b>(1)</b>	
HUM010013	<input type="checkbox"/>	92° Performance fixed front frame (UK 88°) <b>(1)</b>	107
HUM010009	<input type="checkbox"/>	88° Performance compact fixed front frame (UK 92°) <i>(not available with 1cm inset option)</i>	107
HUM010100	<input type="radio"/>	Inset: Standard (0cm)	
HUM010101	<input type="radio"/>	Inset: Small (1cm) <i>(not in combination with 88° frame)</i>	
HUM010103	<input type="checkbox"/>	Inset: (3cm) only for frame angle 88° & 92° <i>(not available for seat width 32 cm)</i>	107
HUM010301	<input type="radio"/>	Castor house position wide	
HUM010302	<input type="radio"/>	Castor house position narrow <i>(only for inset 3cm, not available with 88° performance compact frame)</i>	
HUM010500	<input type="checkbox"/>	Handbike axle <b>(2)</b> <i>(always 3° more than the std axle camber)</i>	293
<b>CENTRE OF GRAVITY (COG)</b>			
Choose only one of the 8 choices		<b>(measured from front of back post in cm)</b>	
COG20	<input type="radio"/>	2.00cm	
COG32	<input type="radio"/>	3.20cm	
COG44	<input type="radio"/>	4.40cm	
COG56	<input type="radio"/>	5.60cm	
COG68	<input type="radio"/>	6.80cm	
COG80	<input type="radio"/>	8.00cm	
COG92	<input type="radio"/>	9.20cm	
COG104	<input type="radio"/>	10.40cm	
<b>FRAME ACCESSORIES</b>			
HUM090031	<input type="checkbox"/>	Suspension up to 70kg <b>(3)</b>	418
HUM090032	<input type="checkbox"/>	Suspension up to 120kg <b>(3)</b>	418
<b>Tip assist (step tube)</b>			
HUM090001	<input type="checkbox"/>	Left <i>HUM090002</i> <input type="checkbox"/>	Right
HUM090003	<input type="checkbox"/>	Pair	<i>(not available with ant-tip)</i>
HUM090009	<input type="checkbox"/>	Crutch holder	<i>(only in combination with tip assist)</i>
<b>Anti-tip tubes swing-away</b>			
HUM090004	<input type="checkbox"/>	Left <i>HUM090005</i> <input type="checkbox"/>	Right
HUM090006	<input type="checkbox"/>	Pair	<i>(not available with tip-assist)</i>
HUM090007	<input type="checkbox"/>	Quickie style anti-tip tubes pluggable	<i>(Pair)</i>
<b>Anti-tip flip-up (user operated)</b> <i>(not in combination with handbike axle)</i>			
HUM090033	<input type="checkbox"/>	Left	124
HUM090008	<input type="checkbox"/>	Right	124
HUM090035	<input type="checkbox"/>	Sport anti-tip, central mount, pluggable	<i>(For use with indoor not contact sports only; e.g. badminton/dancing)</i>
HUM090010	<input type="checkbox"/>	Transit wheels (3" wheels)	180
HUM090016	<input type="checkbox"/>	Frame protection	62
<b>PREPARATION FOR ATTITUDE ADD-ON BIKE</b>			
HUM090211	<input type="checkbox"/>	Preparation for ATTITUDE, without ATTITUDE order <i>(docking pins assembled)</i>	90
HUM090212	<input type="radio"/>	Preparation for ATTITUDE, with ATTITUDE order <i>(docking pins assembled; ATTITUDE order separately)</i>	
<b>FORKS</b>			
HUM080013	<input type="radio"/>	Carbotecture fork <i>(colours: see colour options on page 6)</i>	
HUM080008	<input type="checkbox"/>	Single arm fork <i>(only 111mm length &amp; only with 4" &amp; 5" Alu soft, Alu solid. )( colours: see colour options on page 6)</i>	156
HUM080029	<input type="radio"/>	Aluminium fork <i>(colours: see colour options on page 6. Only available in black, orange and silver)</i>	
HUM080007	<input type="checkbox"/>	Frog Leg fork <i>(only in combination with 4" &amp; 5" solid &amp; all ALU casters; only black)</i>	262
HUM090093	<input type="checkbox"/>	Fork for power add ons, 123mm (Aluminium) <i>(only in black )</i>	62
<b>CASTOR</b>			
HUM080301	<input type="radio"/>	Solid <i>(3"-6", not with One arm fork)</i>	
HUM080302	<input type="radio"/>	Soft <i>(only 5" and 6", not with One arm fork, for 45 mm fork version only)</i>	
HUM080303	<input type="checkbox"/>	Solid light up <i>(not available for 6", not with One arm fork)</i>	62
HUM080304	<input type="checkbox"/>	Solid with aluminium rim, colour rim: silver only <i>(only for 4" and 5", use with one arm fork)</i>	73
HUM080306	<input type="checkbox"/>	Soft with aluminium rim <i>(for 45mm fork version only) (only for 4" and 5", use with one arm fork, pick colours on last page)</i>	73
(1):	no change from 3" & 4" castor to 5" & 6" is possible afterwards		
(2):	Max user weight 100kg. Only with 3°, 6° & 9° camber; not in combination with carbon sideguards, lightweight seat upholstery and not in combination with lightweight compact wheel lock		
(3):	not with COG 6.8 - 10.4cm. Not with rear seat height 37 - 42cm and 24" wheels. Not with rear seat height 37 - 43cm and 25" wheels with standard seat upholstery		

						£
<b>REAR WHEELS</b>						
<i>NOTE: for users with higher user weight than 100kg, we recommend cross spoked wheels</i>						
HUM070001	<input type="radio"/>	Universal Cross spoked	(pair)		(only 24")	
HUM070002	<input type="radio"/>	Design wheel	(pair)		(only 24")	
HUM070003	<input type="radio"/>	Lightweight wheel	(pair)		(Colour anodizing selection for 24" LW wheel on the last page)	
HUM070004	<input type="checkbox"/>	Proton wheel (super lightweight)	(pair)		(user weight limit 100kg)	447
HUM070005	<input type="checkbox"/>	Spinergy wheels 18 spokes, colour hub: silver	(1) (pair)		(black spokes; other colours available via B4M)	405
HUM070006	<input type="checkbox"/>	Spinergy wheels 12 spokes, colour hub: silver	(1) (pair)		(black spokes; other colours available via B4M)	624
HUM070007	<input type="checkbox"/>	Spinergy wheels 12 spokes, colour hub: black	(pair)		(black spokes; other colours available via B4M)	624
HUM070009	<input type="checkbox"/>	Mountain bike wheel			(only with 24", mountain bike tyres & aluminium rim, black coated)	275
HUM070010	<input type="checkbox"/>	Additional mountain bike wheel (pair)			(only with 24" & mountain bike tyres)	275
HUM070000	<input type="radio"/>	Omit rear wheels				
<b>TYRES</b> (choose only one tyre type)						
HUM070101	<input type="radio"/>	Right Run			(only 24")	
HUM070102	<input type="radio"/>	Solid tyres			(only 24")	
HUM070106	<input type="checkbox"/>	Down Town Two Grip			(only 24")	83
HUM070107	<input type="checkbox"/>	Schwalbe One (slight puncture protection)				83
HUM070104	<input type="checkbox"/>	Marathon plus evolution (puncture resistant)				83
<b>HANDRIM</b>						
HUM070201	<input type="radio"/>	Aluminium silver powder coated				
HUM070202	<input type="radio"/>	Aluminium silver anodised				
HUM070204	<input type="radio"/>	Stainless steel				
HUM070203	<input type="checkbox"/>	Aluminium colour powder coated			(pick colour for handrim in colour list charged at £56.00)	56
HUM070205	<input type="checkbox"/>	Titanium			(25" only in combination with Spinergy)	171
HUM070206	<input type="checkbox"/>	Titanium clear powder coated			(pick colour for handrim in colour list; 25" only in combination with Spinergy)	171
HUM070207	<input type="checkbox"/>	Supergrip			(only 24" wheels)	206
HUM070208	<input type="checkbox"/>	Max Grepp			(only 24" wheels)	206
HUM070211	<input type="checkbox"/>	Surge <sup>®</sup> LT				164
Handrim mounting		HUM070317	<input type="radio"/>	Close 12 - 17mm		
		HUM070327	<input type="radio"/>	Wide 27mm	(Not with light weight, Proton & Spinergy wheel and alu handrim)	
		HUM070307	<input type="radio"/>	Extra close 7mm	(only with stainl. steel or Surge, only with 24" Design & Univ. wheel)	
HUM070250	<input type="checkbox"/>	Handrim cover			(only 24" wheels)	83
<b>AXLE TUBE</b>						
HUM010401	<input type="radio"/>	Standard axle tube				
HUM010402	<input type="checkbox"/>	Carbon axle tube			(only for 0° and 3° wheel chamber)	188
<b>CAMBER</b> (choose only one option)						
HUM070600	<input type="radio"/>	0°				
HUM070603	<input type="radio"/>	3°				
HUM070606	<input type="radio"/>	6°			(not available with carbon axle tube)	
HUM070609	<input type="radio"/>	9°			(not available with carbon axle tube)	
<b>SPOKE GUARDS</b>						
HUM070831	<input type="checkbox"/>	Transparent			(only available in 22", 24" & 26")	113
	<input type="checkbox"/>	Spoke guards	Number	Description		113
<p>Please consult the last page of this order form to choose your selection, or to order a copy of our new brochure please call 0845 605 66 88 or e-mail <a href="mailto:literaturerequests@sunmed.co.uk">literaturerequests@sunmed.co.uk</a>. You can also view the brochure online by visiting <a href="http://www.sunrisemedical.co.uk">www.sunrisemedical.co.uk</a> and navigating to Quickie Helium. You will find the brochure under "documentation"</p>						
<b>WHEEL LOCKS</b> (choose only one option)						
HUM060001	<input type="radio"/>	Standard wheel lock				
HUM060003	<input type="radio"/>	Compact wheel lock	(2)			
HUM060004	<input type="checkbox"/>	Lightweight compact	(2) (3) (4)			133
HUM060002	<input type="checkbox"/>	Knee lever wheel lock				52
HUM060010	<input type="checkbox"/>	Extension for knee lever wheel lock, slip-on, swing-away			(only on HUM060002)	35
HUM060006	<input type="checkbox"/>	One arm wheel lock left				140
HUM060007	<input type="checkbox"/>	One arm wheel lock right				140
<p>(1): only with handrim Aluminium silver powdercoated  (2) Not available with 6° &amp; 9° camber (please check NoGo-Matrix)  (3) Adjustment of COG ± 2 cm on a chair in the field.  (4) Not for use with add-on bikes</p>						


								£			
<b>ACCESSORIES</b>											
HUM090000	<input type="radio"/>	Tool kit									
HUM090100	<input type="radio"/>	Transportation brackets - Crash test option				<i>(not available on Helium Pro)</i>					
HUM090017	<input type="checkbox"/>	Snoll on <b>(1)</b>				<i>(Max user weight: 100kg)</i>		499			
HUM090024	<input type="checkbox"/>	High pressure pump						68			
HUM090025	<input type="checkbox"/>	Air pump						14			
HUM090030	<input type="checkbox"/>	Back pack						55			
HUM090020	<input type="checkbox"/>	Net back bag						20			
HUM090019	<input type="checkbox"/>	Caddy				<i>(not with platform footplates)</i>		55			
HUM090026	<input type="checkbox"/>	Mobile phone pouch						41			
HUM090014	<input type="checkbox"/>	Handbike re-enforcement tubes				<i>(Only with handbike axle)</i>		319			
<b>COLOURS</b>											
<i>Please visit <a href="http://www.quickiewheelchairs.eu/xenonsa">www.quickiewheelchairs.eu/xenonsa</a> to choose your colours</i>											
	Frame	Handrim	Hub/Rim (2)		Frame	Handrim	Hub/Rim (2)				
		£56	£158			£56	£158				
01 White RAL 9010	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		26 Silver Texture	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
02 Brilliant Silver	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		27 Graphite metallic	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
03 Silver Matt	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		28 Platinum metallic	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
04 Silver	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		29 Black	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
06 Sun Yellow	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		31 Matt Black	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
07 Brilliant Lime	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		32 Signal White	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
08 Matt Lime	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		33 Pink	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
14 Blue	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		34 Dark Chocolate Metallic	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
15 Brilliant Marine	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		35 Skyblue	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
16 Matt Marine	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		36 Cyan-Metallic	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
19 Brilliant Purple	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		37 Mars Red (Orange)	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
20 Matt Purple	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		39 Cappuccino	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
21 Red Ruby	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		40 Kiwi	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
22 Brilliant Red	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		41 Flamingo-Glitter	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
23 Matt Red	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		42 Matt-Anthracite	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
24 Tutti Frutti Orange	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>		43 Azure Blue	<input type="radio"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<b>24" LIGHTWEIGHT WHEELS</b>											
HUM100336	<input type="radio"/>	Orange anodised (3)		HUM100337	<input type="radio"/>	Blue anodised (3)					
HUM100338	<input type="radio"/>	Red anodised (3)									
<b>BACKREST BRACKET (Double locking back only)</b>											
HUM030035	<input type="radio"/>	Black		HUM030037	<input type="radio"/>	Silver		HUM030039	<input type="radio"/>	Blue (6)	
HUM030036	<input type="radio"/>	Orange		HUM030038	<input type="radio"/>	Red (6)					
<b>CASTOR FORKS</b>											
HUM080041	<input type="radio"/>	Black		HUM080042	<input type="radio"/>	Silver		HUM080045	<input type="radio"/>	Blue	
HUM080043	<input type="radio"/>	Orange		HUM080044	<input type="radio"/>	Red					
<b>CASTOR WHEEL (Only soft with aluminium rim)</b>											
HUM080311	<input type="radio"/>	Black		HUM080313	<input type="radio"/>	Orange		HUM080315	<input type="radio"/>	Blue	
HUM080312	<input type="radio"/>	Silver		HUM080314	<input type="radio"/>	Red					
<b>STICKERS</b>											
HUM100550	<input type="radio"/>	Black		HUM100552	<input type="radio"/>	Orange		HUM100554	<input type="radio"/>	Blue	
HUM100551	<input type="radio"/>	Silver		HUM100553	<input type="radio"/>	Red					
<b>AXLE STEM SPACER (Not available on Helium Pro)</b>											
HUM100560	<input type="radio"/>	Silver		HUM100562	<input type="radio"/>	Red					
HUM100561	<input type="radio"/>	Orange		HUM100563	<input type="radio"/>	Blue					
<b>FRAME TRANSFERS (Dots)</b>											
HUM100570	<input type="radio"/>	Black		HUM100571	<input type="radio"/>	Silver		HUM100574	<input type="radio"/>	Blue	
HUM100572	<input type="radio"/>	Orange		HUM100573	<input type="radio"/>	Red		HUM100575	<input type="radio"/>	Matt Black	
<b>(1) Not in combination with performance footrest and high mounted footrest</b>											
<b>(2) Only with universal cross spoke rear wheel</b>											
<b>(3) Anodised colours for hub and rim, black spokes. Only with light weight wheel in 24"</b>											
<b>(4) Froglegs fork and fork for power add-ons (HUM090093) only available in black. Backrest bracket only in combination with HUM030102 "double locking back" in colours, otherwise always black.</b>											
<b>(5) Only with 111mm forks (including one-arm fork) in combination with HUM080302 Castor, soft or HUM080306 Castor soft with aluminium rim</b>											
<b>(6) Only with carbotechnology forks</b>											
<b>Base Price</b>								£			
<b>Option price</b>								£			
<b>Total Price</b>								£			

**Please follow this measuring method:**

The following method is recommended on the basis of extensive practical experience. Deviations are possible in individual instances.


**Seat width (SW):**

- a) *For the user:* Allow space between sideguard and hip on each side of the chair (approx. 1 - 2cm altogether)
- b) *On the wheelchair:* From outside frame tube to outside frame tube.



**Seat depth (SD):**

- a) *For the user:* When sitting upright: leave space for 2 - 3 fingers between front edge of the seat upholstery and back of knee (approx. 3 - 6 cm). Care must be taken with adj. upholstery.
- b) *On the wheelchair:* From back tube to front of seat upholstery.



**Front seat height:**


- a) *For the user:* When sitting with 90° kneeangle: From back of knee to bottom of shoe sole.
- b) *On the wheelchair:* From front of seat tube to the ground . (min. 3cm longer than user measurement) Remember to include the cushion height.

**Rear seat height:**

- b) *On the wheelchair:* From back of seat tube to the ground. **Note:** Max. seat height difference rear to front is 13 cm

**Back height:**

- a) *For the user:* From the seat, including cushion, to be no more than 1 - 2 cm below the lower shoulder blade.
- b) *On the wheelchair:* From back of seat tube to upper edge of backrest upholstery. **Note:** The user's ability to sit upright has a major influence on the back height and it can therefore differ considerably




**Back angle:**

**Important:** The back height depends on the degree of disability. 0° is the angle 90° to the seat. Negative angles are <90° and positive angles are >90°.

**Low Seat Frame vs. High Seat Frame:**

The high frame is 4 cm higher than the low frame. So you can choose shorter and lighter forks to minimize the weight of your wheelchair.



**Technical Data Overview:**

Overall width:	Overall width with:	24" rear wheel	25" rear wheel
Overall width:	For 0° camber:	SW + 17 cm	SW + 17 cm
Max overall height	For 3° camber:	SW + 23 cm	SW + 24 cm
Max overall length	For 6° camber:	SW + 29 cm	SW + 30 cm

Please note that the matrix is based on a std configuration:  
**Std chair 40 x 40 cm; 3° camber; 5,60 cm CoG; Seat angle of 8°**  
 All changes in the configuration might change the seat heights.  
 Manufacturing tolerances of +/- 1 cm

## SEAT UPHOLSTERIES AND CUSHIONS



*HUM020002-020002*  
Seat sling, nylon, tension adjustable



*HUM020005 -020007*  
Seat sling clip style (airflow)



*HUM090018*  
Lap belt with buckle, metal

## BACKREST



*HUM030003*  
Oval Shaped Backrest



*HUM030102*  
Backrest, angle-adjustable, lock when folded



*HUM030201*  
Push handles long



*HUM030203*  
Push handles, fold-down



*HUM030204*  
Push handles, height-adjustable



*HUM030316*  
EXO backrest sling



*HUM030317*  
EXO PRO backrest sling



**FREESTYLE BACK, ANGLE ADJUSTABLE, FOLDABLE**



	with foam	w/o foam
Back height	240 mm	165 mm
	310 mm	240 mm
Back width	320 mm	280 mm
	370 mm	330 mm
	380 mm	380 mm
Contour		



**JAY BACKS**



Jay 3 Back  
Shallow Contour



Jay 3 Back  
Mid Contour



Shallow Contour

**SIDEGUARD & ARMREST**



*HUM101  
Slim style sideguard*



*HUM040104  
Carbon fibre, cold weather  
protection, fender*



Sideguard, composite,  
removable (black)



Single-post sideguard



Single-post sideguard,  
armrest (33 cm)



Tubular armrest, swing  
away, removable, padded

## FOOTBOARDS



**HUM050029**  
Tubular footrest  
lightweight



**HUM050030**  
Tubular footrest lightweight  
with flat cover



**HUM050025**  
Platform footrest LW  
composite plate



**HUM050026**  
Platform footrest LW  
carbon plate



**AGN050125**  
Performance footrest  
alu angle adj



**AGN050127**  
Foot positioning plate




**AGN050020**  
Composite footplate  
divided



**AGN050022**  
Alu footplate divided

## FRAME VERSIONS



**Frame Standard vs. Frame Long**

Difference is 20 mm



**Frame Angles**

### Frame inset table

Inset		Seat Width (mm)							
		320	340	360	380	400	420	440	460
0 mm	Frame Width Inside	270	290	310	330	350	370	390	410
	Frame Width Casterhouse	440	460	480	500	520	540	560	580
10 mm	Frame Width Inside	250	270	290	310	330	350	370	390
	Frame Width Casterhouse	420	440	460	480	500	520	540	560
30 mm	Frame Width Inside	NA	230	250	270	290	310	330	350
	Frame Width Casterhouse wide	NA	460	480	500	520	540	560	580
	Frame Width Casterhouse small	NA	400	420	440	460	480	500	520

## FRAME ACCESSORIES



Tip assist left



Anti-tip tubes SA



Anti tip tubes pluggable



Sport anti-tip



Anti tip flip up  
User operated

## FORKS



HUM080013 & HUM080306  
Castor fork, Carbotechnology & Castor solid (soft) aluminum



HUM080008  
Castor fork, aluminum, one-arm, red

## REAR WHEELS



HUM070001  
Universal wheel



HUM070002  
Design wheel



HUM070003  
Lightweight wheel



HUM070004  
Proton wheel



HUM070010  
Mountain bikewheel

## TYRES



HUM070101  
Schwalbe Right Run



HUM070106  
Schwalbe Down Town



HUM070107  
Schwalbe One



HUM070104  
Schwalbe Marathon Plus  
Evolution



HUM070102  
Tyre, puncture proof  
Evolution

## HANDRIMS



HUM070208  
Handrims Maxgrepp



HUM070211  
Handrims Surge®

## WHEEL LOCKS



HUM060001  
Standard brake



HUM060002  
Knee brake



HUM060003



HUM060004

## NoGo Matrix for Lightweight Compact Brake

Seat height values are valid with 56 mm COG, seat depth 400 mm and 8° seat angle,  
axle stem assembled with one spacer

### SD 420 - 480 mm

24"	Front Seat Height														25"	Front Seat Height															
Cog	430	440	450	460	470	480	490	500	510	520	530	540	550	560	570	Cog	430	440	450	460	470	480	490	500	510	520	530	540	550	560	570
20	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	20	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
32	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	32	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
44	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	44	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
56	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	56	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
68	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	68	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
80	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	80	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
92	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	92	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
104	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	104	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	

### SD 400 mm

24"	Front Seat Height														25"	Front Seat Height															
Cog	430	440	450	460	470	480	490	500	510	520	530	540	550	560	570	Cog	430	440	450	460	470	480	490	500	510	520	530	540	550	560	570
20	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	20	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
32	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	32	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
44	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	44	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
56	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	56	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
68	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	68	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
80	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	80	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
92	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	92	X	X	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
104	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	104	X	X	X	X	X	X	X	✓	✓	✓	✓	✓	✓	✓	

### SD 340 - 380 mm

24"	Front Seat Height														25"	Front Seat Height															
Cog	430	440	450	460	470	480	490	500	510	520	530	540	550	560	570	Cog	430	440	450	460	470	480	490	500	510	520	530	540	550	560	570
20	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	20	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
32	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	32	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
44	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	44	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
56	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	56	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
68	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	68	X	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	
80	X	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	80	X	X	X	X	X	X	✓	✓	✓	✓	✓	✓	✓	✓	
92	X	X	X	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	92	X	X	X	X	X	X	X	X	✓	✓	✓	✓	✓	✓	
104	X	X	X	X	X	X	X	✓	✓	✓	✓	✓	✓	✓	✓	104	X	X	X	X	X	X	X	X	X	X	✓	✓	✓	✓	

Possible

NoGo


## ACCESSORIES



**AXLE STEM COLOUR**



**FRAME GRAPHIC "FLOATING DOTS" COLOUR**



**24 " LIGHT WEIGHT WHEEL COLOURS**



*HUN100336  
Lightweight wheel, anodized orange*



*HUM100337  
Lightweight wheel, anodized blue*



*HUM100338  
Lightweight wheel, anodized red*





HUM070823 - Princess



HUM070808 - Dancing Bear



HUM070818 - Jungle Fever



HUM070816 - Happy Lions



HUM070832 - Water World



HUM070825 - Randolf



HUM070819 - Monkeys



HUM070821 - Pirates



HUM070811 - Feather



HUM070820 - Oriental



HUM070824 - Racing



HUM070809 - Dirt Race



HUM070829 - Superhero



HUM070814 - Graffiti



HUM070806 - Chrome Wheel



HUM070817 - Hot Wheels



HUM070802 - Basketball



HUM0708036 - Basketball 2



HUM070822 - Pool



HUM070812 - Flying Water



HUM070828 - Space



HUM070827 - Silk



HUM070801 - Abstract



HUM070807 - Creative



HUM070830 - Tiger



HUM070833 - Zebra



HUM070805 - Cheese



HUM070804 - Bibbo



HUM070813 - Galactica



HUM070826 - Rotation



HUM070810 - Fantasy



HUM070815 - Graphic