

Invoice to:	
Name:	
Street:	
Town:	Postal code:
Country:	
Tel:	Fax:
Order date:	
Purchase order:	

Deliver to:	
Name:	
Street:	
Town:	Postal code:
Country:	
Tel:	Fax:
Marked for:	

Base Price £4,395



Email: orders@sunmed.co.uk

Fax number: +44 (0)845 605 66 89

Maximum user weight limit 140kg (22 Stone)

Standard Features (Salsa R²) <ul style="list-style-type: none"> ● Rear wheel drive base ● 12.5" drive wheels ● Drive wheel suspension ● 60Ah batteries ● 60cm overall width (with 12" wheels) 	<ul style="list-style-type: none"> ● Red, Blue or White Shroud ● Motor/wheel disengagement ● Transportation approved ● Width, height & depth adjustable seat ● Flip back, angle, height & depth adj. arms 	Seating Standard Options <ul style="list-style-type: none"> ● Tension adj. backrest sling ● 70° swing-away legrests ● Angle & depth adjustable footplates ● Adjustable lap strap ● Height and angle adjustable backrest
--	--	---

STANDARD SEAT: FRAME / SEAT SIZE		
	Width Setting	Depth setting
SW41	<input type="radio"/> 41cm (16")	SD41
SW43	<input type="radio"/> 43cm (17")	SD43
SW46	<input type="radio"/> 46cm (18")	SD46
SW48	<input type="radio"/> 48cm (19")	SD48
SW51	<input type="radio"/> 51cm (20")	SD51
STR020114 : Adjustable Seating System Seat frame is fully adjustable from 41cm to 51cm (16" - 20") in width and depth. Independently adjustable in 2.5cm (1") increments		
		£

STANDARD SEAT: CUSHION		
STR020109	● Firm seat board	(standard with standard seat)
STR020111	<input type="radio"/> 5cm (2") standard foam cushion 41cm (16") wide	(used for seat width 41 - 46cm)
STR020112	<input type="radio"/> 5cm (2") standard foam cushion 48cm (19") wide	(used for seat width 48 - 51cm)
STR020110	<input type="radio"/> Omit seat cushion	

JAY SEAT DIMENSION																																													
<table border="1"> <thead> <tr> <th colspan="6">Width</th> </tr> <tr> <th></th> <th>41</th> <th>43</th> <th>46</th> <th>48</th> <th>51</th> </tr> </thead> <tbody> <tr> <th rowspan="5">Depth</th> <td>41</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>43</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>46</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>48</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> <tr> <td>51</td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> <td><input type="checkbox"/></td> </tr> </tbody> </table>	Width							41	43	46	48	51	Depth	41	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	43	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	46	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	48	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	51	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	JAY020006 <input type="checkbox"/> JAY Easy Visco, Flat Base JAY020005 <input type="checkbox"/> JAY Easy Fluid, Flat Base JAY020010/12 <input type="checkbox"/> JAY Balance Air (single valve, microclimatic cover) JAY020011/13 <input type="checkbox"/> JAY Balance Fluid (microclimatic cover) JAY020014/16 <input type="checkbox"/> JAY Balance Air, Positioning (single valve, microclimatic cover) JAY020015/17 <input type="checkbox"/> JAY Balance Fluid, Positioning (microclimatic cover) JAY020019 <input type="checkbox"/> JAY J2 JAY020020 <input type="checkbox"/> JAY J2 Deep Contour	132 308 503 468 568 533 492 525
	Width																																												
		41	43	46	48	51																																							
	Depth	41	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																						
		43	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																						
		46	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																						
48		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																							
51		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>																																							

JAY COMFORT SEAT		
		
	 JCL020011	 JCL020012
	 JCL020003	 JCL020007
JCL020011	<input type="checkbox"/> Comfort seat pan, without threshold	(thickness 2cm) 51
JCL020012	<input type="checkbox"/> Comfort seat pan, with 5cm (2") threshold	(thickness 2cm) 84
JCL020003	<input type="checkbox"/> Standard cushion - stretched rubber upholstery	(thickness 8cm) 155
JCL020007	<input type="checkbox"/> Optimum comfort - stretched rubber upholstery	(anatomically formed - thickness 9cm) 199

● = Standard ○ = No cost option □ = Cost option

£

SEATING OPTION - SELECT SEAT POSITIONING OPTION WITH ASSOCIATED SEAT HEIGHT

SEAT POSITIONING OPTION	12.5" Drive Wheels		14" Drive Wheels		Option Price
	Low Seat Height	Tall Seat Height	Low Seat Height	Tall Seat Height	
	<i>(No Extension Kit)</i>	<i>(Inc Extension Kit - SAL090022)</i>	<i>(No Extension Kit)</i>	<i>(Inc Extension Kit - SAL090022)</i>	
SAL130001 Manual adjustable single pivot tilt (3° preset)	○ 42cm (16.5")	N/A	○ 43cm (17")	N/A	
SAL130002 Manually adjustable C.O.G Tilt (3° preset)	○ 44cm (17.3")	○ 47cm (18.5")	○ 45cm (17.7")	○ 48cm (19")	
SAL130003 Powered single pivot tilt 3° to 20°	□ 42cm (16.5")	N/A	□ 43cm (17")	N/A	425
SAL130004 Powered C.O.G Tilt: 0° to 30°	□ 46cm (18")	□ 49cm (19.5")	□ 47cm (18.5")	□ 50cm (20")	651
SAL130005 Powered negative C.O.G tilt: -5° to 25°	□ 46cm (18")	□ 49cm (19.5")	□ 47cm (18.5")	□ 50cm (20")	785
SAL130006 Powered lift & tilt: 30cm lift, 30° tilt	□ 47cm (18.5")	□ 50cm (20")	□ 48cm (19")	□ 51cm (20.5")	1748

BACKREST RECLINE

STR030018	○ Mechanical fixed recline back	<i>(-4° to 16°)</i>	
STR030025	□ Manual adjustable recline back	<i>(Attendant operation only - gas strut, -5° to 30°, 41 to 46 cm or 48 to 51 cm depth adjustment range)</i>	211
STR030026	□ Powered recline back	<i>(-5° - 30°, 41 to 46 cm or 48 to 51 cm depth adjustment range)</i>	421

STANDARD BACKREST

STR030300	○ Omit Backrest Upholstery	<i>(Chair supplied with NO back upholstery option)</i>	
STR030301	○ Tension adj. sling backrest (48cm /19" height)	<i>(2.5cm up and down adjustment, with std foam cushion only)</i>	
Backrest Height Setting	BH46* ○ 46cm (18")	BH485 ○ 48.5cm (19")	BH51 ○ 51cm (20")
* Not available with manual or power reclining back			
SB43886	□ Solid covered backrest	<i>(Available via B4me department)</i>	POA
Backrest Height Setting	○ 50cm (19.5")	○ 57cm (22.5")	

JAY SHAPE BACKREST - MAX USER WGT 135KG

(Premium tension adjustable backrest upholstery - Not available on SW 51cm)

JAY080001	□ JAY SHAPE Backrest	<i>(Not available in 51cm Seat Width)</i>	193
JAY080002	□ JAY SHAPE Backrest - Oversized (Supplied 1" larger than seat width to allow greater immersion)	<i>(N/A in 48 or 51cm Seat Width)</i>	193
Back Height Setting	BH46 ○ 46cm (18")	BH485 ○ 48.5cm (19")	BH51 ○ 51cm (20")
	BH535 ○ 53.5cm (21")	BH56 ○ 56cm (22")	BH585 ○ 58.5cm (23")
JAY080010	□ JAY SHAPE Lateral Thoracic supports (Pair) - Support 15 x 15 cm (Pair)		193
JAY080011	□ JAY SHAPE Lateral Thoracic supports (Pair) - Support 10 x 12 cm (Pair)		193

JAY COMFORT BACKRESTS

JCL030004	□ Shallow contour (standard velour) - 50cm back height	<i>(hard back)</i>	153	
JCL030008	□ Shallow contour (standard velour) - 57cm back height	<i>(vertical elastic straps)</i>	222	
JCL030006	□ Shallow contour (standard velour) - 57cm back height	<i>(horizontal tension adjustable straps)</i>	222	
JCL030012	□ Deep contour - passive - 57cm back height	<i>(horizontal tension adjustable straps)</i>	<i>(4" Integrated lateral support)</i>	257
JCL030010	□ Deep contour - active - 57cm back height	<i>(horizontal tension adjustable straps)</i>	<i>(4" integrated lumbar lateral support)</i>	257

LATERAL SUPPORTS

STR030601	□ Lateral supports short (pair)	<i>(Recommended for TA backrests)</i>	273
STR030602	□ Lateral supports long (pair)	<i>(Recommended for JAY Comfort backrests)</i>	273

JAY3 BACKRESTS

J3 Back - Shallow Contour

Upper Thoracic:	□ JAY030012 Short - 42 cm	□ JAY030019 Medium - 46 cm	□ JAY030025 Tall - 50 cm	572
Shoulder Height:	□ JAY030015 Short - 53 cm	□ JAY030022 Medium - 57 cm	□ JAY030028 Tall - 61 cm	572

J3 Back - Mid Contour

Upper Thoracic:	□ JAY030013 Short - 42 cm	□ JAY030020 Medium - 46 cm	□ JAY030026 Tall - 50 cm	572
Shoulder Height:	□ JAY030016 Short - 53 cm	□ JAY030023 Medium - 57 cm	□ JAY030029 Tall - 61 cm	572

J3 Back - Deep Contour

Upper Thoracic:	□ JAY030014 Short - 42 cm	□ JAY030021 Medium - 46 cm	□ JAY030027 Tall - 50 cm	627
Shoulder Height:	□ JAY030017 Short - 53 cm	□ JAY030024 Medium - 57 cm	□ JAY030030 Tall - 61 cm	627

J3 Back - Posterior Deep Contour

Upper Thoracic:	□ JAY030031 Medium - 46 cm	627
Shoulder Height:	□ JAY030032 Short - 53 cm	627

JAY3 HEADRESTS

JAY030006	□ Jay Std. Headrest (J3EUHDXH & J3EUHDLSP)	<i>(only with J3 backs)</i>	129
JAY030004	□ Jay Flip Back Headrest (J3EUHDFH & J3EUHDLSP)	<i>(only with J3 backs)</i>	239

● = Standard ○ = No cost option □ = Cost option			£
STANDARD HEADRESTS			
STR030416	<input type="checkbox"/> Comfort headrest (Winged style with adj. side pads)	(Mounted via push handles)	174
STR030417	<input type="checkbox"/> Curved headrest (23 x 11cm / 9" x 4.5" pad)	(Mounted via push handles)	174
JCL040001	<input type="checkbox"/> Jay Comfort: standard comfort seating headrest	(Mounted via Comfort backrest)	196
JCL040002	<input type="checkbox"/> Jay Comfort: flexible comfort seating headrest (winged style)	(Mounted via Comfort backrest)	221
ARMRESTS			
STR040122	<input type="checkbox"/> Flip-up armrest (height, angle adjustable)	(Height Angle and depth adjustable via outside of armrest)	
STR040127	<input type="checkbox"/> Reclining flip-back armrest	(Armrest stays level when backrest Reclines)	171
STR040126	<input type="checkbox"/> Single post armrest	(2" inward width adjustment via sideguard and armpad / 4" depth adjustment via sideguard)	216
ARMPADS			
STR040112	<input type="checkbox"/> Standard Armpad - 385 x 60mm pad		
STR040113	<input type="checkbox"/> Contoured Armpads 41cm (16")		57
STR040114	<input type="checkbox"/> Jay Comfort Padded Armpad 41cm (16")		113
STR040118	<input type="checkbox"/> Jay Comfort P-Shaped Padded Armpad 45 cm (18")	(Only available with Single Post Armrest)	158
STR040119	<input type="checkbox"/> Hemi Armpad (Medium forearm pad & Small hand Pad) 28cm (11")	Pad 11 x 18cm (4.5" x 7")	408
POSITIONING ITEMS			
STR030604	<input type="checkbox"/> Lateral Thigh Guide/Support	(4" Depth adjustment / 2" inward width adjustment) (Not available with Single post Armrest)	158
STR030605	<input type="checkbox"/> Pommel	(Height / depth adj. / Quick release removable) (Not available with centre mount legrests)	193
LEGRESTS & FOOTPLATES			
STR050040	<input type="checkbox"/> Standard hanger receivers	(Comes standard, with 70 degree s/a hangers, not with Uni style)	
STR050041	<input type="checkbox"/> Offset hanger receivers	(20mm outwards offset each side, comes std with all ELR and ALR hangers, not with Uni style)	
STR050015	<input type="checkbox"/> 70° swing-away legrests		
XL0025	<input type="checkbox"/> Uni style hangers and footplates	(Available via B4M - (LLL 32-47cm (12.5"-18.5"). Not available with 0° tilt options)	
STR050016	<input type="checkbox"/> Manual elevating legrest		271
STR050017	<input type="checkbox"/> Powered elevating legrest		498
STR050018	<input type="checkbox"/> Powered articulating legrest		710
STR050036	<input checked="" type="checkbox"/> Extension Tube - (Lower leg length: 34cm - 48cm / 13" - 19")	(measured from seat pan to footplate)	
OPTIONS FOR MANUAL AND POWER ELEVATING LEGRESTS			
STR050038	<input type="checkbox"/> Comfort calf pads	(Not available with centre mount or 70° legrests)	102
STR050039	<input type="checkbox"/> Auto connector hanger receivers		121
FOOTPLATES (FOR SWING AWAY LEGRESTS)			
STR050022	<input type="checkbox"/> Aluminium angle and depth adjustable		
STR050023	<input type="checkbox"/> Footboard (flip-up)	(Only available with 70° swing-away legrests)	121
VR2 CONTROLS (70A ON 6KPH, 90A ON 10KPH) (Max 2 powered seat options - for more than 2 powered options choose R-net)			
SAL030050/1	<input type="checkbox"/> VR2 remote	(Remote type selected according to options chosen)	
STR110020	<input type="checkbox"/> VR2 attendant only	(For attendant only VR2 controls option above must be selected also)	96
STR110021	<input type="checkbox"/> VR2 dual control (mounted RH)	(For dual control VR2 remote must be selected also)	452
R-NET CONTROLS (FOR SPECIAL CONTROLS OPTIONS SEE PAGE 5)			
STR110015/6	<input type="checkbox"/> R-net LED remote	(Remote type selected according to options chosen)	466
STR110017	<input type="checkbox"/> R-net colour LCD remote	(Integrated colour display with 2 stereo jack inputs for buddy buttons)	606
STR110020	<input type="checkbox"/> R-net attendant only (mounted RH)	(For attendant only R-net LED/LCD remote must be selected also)	96
STR110021	<input type="checkbox"/> R-net dual control (mounted RH)	(For dual control R-net LED/LCD remote must be selected also)	452
SAL090062	<input type="checkbox"/> Swing-away tray mounted centre control	(R-Net LED as standard, VR2 via B4M)	436
STR110009	<input type="checkbox"/> Gyro Control		220
SAL090068	<input type="checkbox"/> R-net Bluetooth control module	(Compatible with Windows PC and Android systems only)	633
CONTROL MOUNTING OPTIONS			
Parallel swing-away controller mount STR110052 RH ○ ○ LH STR110053			
SAL090062	<input type="checkbox"/> Swing-away tray mounted centre control R-NET LED	(Mounted RH Flip Away)	436
SAL110143	<input type="checkbox"/> Swing-away tray mounted centre control R-NET LCD	(Mounted RH Flip Away)	1606
B4M	<input type="checkbox"/> Centre Bar flip-up controller mount	(B4M Option / (Mounted RH Flip away)	POA

● = Standard

○ = No cost option

□ = Cost option

£

LIGHTS & INDICATORS			£
SAL110201	<input type="checkbox"/>	LED lights & indicators (STVZO APPROVED)	395
SAL020001	<input type="checkbox"/>	LED lights & indicators - Standard	395
KERB CLIMBER			
SAL090038	<input type="checkbox"/>	Kerb climber	163
WHEELS & CASTORS			
STANDARD CASTORS			
SAL080307	<input type="checkbox"/>	9" solid castors	
SAL080308	<input type="checkbox"/>	10" solid castors	
SAL080309	<input type="checkbox"/>	10" pneumatic castors	
STANDARD REAR WHEELS			
SAL070029	<input type="checkbox"/>	12.5" solid drive wheels <i>(60cm / 24" overall chair width - Use only with 9" castors)</i>	
SAL070030	<input type="checkbox"/>	12.5" pneumatic drive wheels <i>(60cm / 24" overall chair width - Use only with 9" castors)</i>	
SAL070031	<input type="checkbox"/>	14" solid V-type drive wheels <i>(61cm / 24.5" overall chair width)</i>	
SAL070032	<input type="checkbox"/>	14" pneumatic V-type drive wheels <i>(61cm / 24.5" overall chair width)</i>	
SAL070033	<input type="checkbox"/>	14" solid knobby drive wheels <i>(62cm / 25" overall chair width)</i>	85
SAL070034	<input type="checkbox"/>	14" pneumatic knobby drive wheels <i>(62cm / 25" overall chair width)</i>	85
MOTOR SPEED			
SAL100602	<input type="checkbox"/>	3.2kph (2mph indoor setting) <i>(6kph motor with reduced speed via programming)</i>	
SAL100601	<input type="checkbox"/>	6kph (4mph) 2 Pole Motors	
SAL100603	<input type="checkbox"/>	10kph (6mph) 2 Pole Motors <i>(Lights and rear review mirror required by law)</i>	271
SAL100604	<input type="checkbox"/>	10kph (6mph) 4-pole High Torque Motors <i>(Lights and rear review mirror required by law)</i>	334
COLOUR			
SAL140001	<input type="checkbox"/>	White	
SAL140003	<input type="checkbox"/>	Blue	
SAL140002	<input type="checkbox"/>	Red	
BATTERIES			
SAL160001	<input type="checkbox"/>	60Ah gel batteries	
		SAL160002 <input type="checkbox"/> 70Ah gel batteries	113
CHARGER			
SAL160100	<input checked="" type="checkbox"/>	8A Charger	
ACCESSORIES			
SAL090000	<input checked="" type="checkbox"/>	Transportation tie down points	
SAL090005	<input checked="" type="checkbox"/>	Adjustable lap strap	
SAL090040	<input type="checkbox"/>	Swing-away tray table 41 - 48cm	174
SAL090041	<input type="checkbox"/>	Swing-away tray table 46 - 51cm	174
SAL090069	<input type="checkbox"/>	Bluetooth USB mouse receiver <i>(if there is no Bluetooth function on the computer)</i>	67
SAL090067	<input type="checkbox"/>	Rear view mirror	39
SAL090003	<input type="checkbox"/>	Crutch holder	61
TRANSPORTATION OPTIONS			
SAL090039	<input checked="" type="checkbox"/>	Prepared for Dahl docking station option	
QUICKIE SALSA R² BASE PRICE			£4,395
Option Price			£
VAT			£
Total Price			£

- = Standard
- = No cost option
- = Cost option



		£				
R-NET SPECIAL CONTROLS <i>(R-NET OMNI required; 2 buddy buttons recom.)</i>						
STR110101	<input type="checkbox"/> R-Net OMNI colour display & special control input module	876				
STR110102	<input type="checkbox"/> OMNI in scanning mode <i>(c/w 10 assignable functions through DB9 interface)</i>	876				
Front control mountings						
STR110110	● R-NET OMNI Flexi mount only					
	<table style="width: 100%; border: none;"> <tr> <td style="width: 50%; text-align: center;">Left</td> <td style="width: 50%; text-align: center;">Right</td> </tr> <tr> <td style="text-align: center;">STR110121 ○</td> <td style="text-align: center;">○ STR110122</td> </tr> </table>	Left	Right	STR110121 ○	○ STR110122	
Left	Right					
STR110121 ○	○ STR110122					
SPECIAL JOYSTICKS AS HAND CONTROL <i>(R-NET OMNI required; 2 buddy buttons recommended.)</i>						
STR110130	<input type="checkbox"/> Switch-It Versaguide proportional joystick as hand control	1693				
STR110131	<input type="checkbox"/> Switch-It Micro Guide as hand control	2276				
STR110132	<input type="checkbox"/> Switch-It Micro Pilot (0 way) joystick as hand control	2509				
Mounting						
STR110150	<input type="checkbox"/> Fixed front mount					
STR110151	<input type="checkbox"/> Swing-away front mount	59				
	<table style="width: 100%; border: none;"> <tr> <td style="width: 50%; text-align: center;">Left</td> <td style="width: 50%; text-align: center;">Right</td> </tr> <tr> <td style="text-align: center;">STR110161 ○</td> <td style="text-align: center;">○ STR110162</td> </tr> </table>	Left	Right	STR110161 ○	○ STR110162	
Left	Right					
STR110161 ○	○ STR110162					
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>STR110101</p> </div> <div style="text-align: center;">  <p>STR110131</p> </div> <div style="text-align: center;">  <p>STR110130</p> </div> <div style="text-align: center;">  <p>STR110132</p> </div> </div>						
SPECIAL HEAD RELATED INPUT DEVICES						
<i>(if you have already chosen a standard headrest then switch-It head array is not possible)</i>						
STR110133	<input type="checkbox"/> Switch-It 3 switch head array / adult <i>(5 switch and other versions available via B4M) (OMNI required, 2 BB recommended)</i>	3443				
STR110134	<input type="checkbox"/> Switch-It Micro Pilot (0 way) joystick as chin control <i>(OMNI required, 2 BB recommended)</i>	2509				
STR110135	<input type="checkbox"/> Proportional chin control <i>(proportional joystick; OMNI required; 2 BB recommended)</i>	876				
STR110136	<input type="checkbox"/> Sip & Puff control <i>(tube and mouthpiece; OMNI required; 2 BB recommended)</i>	293				
STR110137	<input type="checkbox"/> Switch-It Micro Guide as chin control <i>(OMNI required, 2 BB recommended)</i>	2276				
STR110138	<input type="checkbox"/> Switch-It proportional head control <i>(OMNI required, 2 BB recommended)</i>	3210				
STR110139	<input type="checkbox"/> Switch-It Versaguide proportional joystick as chin control <i>(OMNI required, C/W Switch-It flex switch kit)</i>	1810				
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>STR110139</p> </div> <div style="text-align: center;">  <p>STR110133</p> </div> </div>						
CHIN CONTROL / SIP AND PUFF SWING-AWAY MOUNTING OPTIONS						
STR110152	<input type="checkbox"/> Manual s/a option	228				
	○ Powered s/a option (select type below) <i>(not available with forus attendant control; s/a operation through separate switch)</i>					
	<table style="width: 100%; border: none;"> <tr> <td style="width: 50%; text-align: center;">Left</td> <td style="width: 50%; text-align: center;">Right</td> </tr> <tr> <td style="text-align: center;">STR110161 ○</td> <td style="text-align: center;">○ STR110162</td> </tr> </table>	Left	Right	STR110161 ○	○ STR110162	
Left	Right					
STR110161 ○	○ STR110162					
STR110153	<input type="checkbox"/> Operated with Piko button <i>(operated via one Piko button - select colour below)</i>	1693				
STR110154	<input type="checkbox"/> Operated with ribbon switch	1810				
<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>Powered swing-away chin option</p> </div> <div style="text-align: center;">  <p>Piko buttons</p> </div> </div>						
PIKO BUTTONS (Ø 3 cm)		per switch				
STR110191	<input type="checkbox"/> Red x 1	85				
STR110193	<input type="checkbox"/> Green x 1	85				
STR110195	<input type="checkbox"/> Yellow x 1	85				
STR110197	<input type="checkbox"/> Blue x 1	85				
STR110192	<input type="checkbox"/> Red x 2	85				
STR110194	<input type="checkbox"/> Green x 2	85				
STR110196	<input type="checkbox"/> Yellow x 2	85				
STR110198	<input type="checkbox"/> Blue x 2	85				

JAY COMFORT SEAT OPTIONS & ACCESSORIES

JAY COMFORT

Seat pan, with 5cm (2") threshold



Upholstered wooden seat board complete with 5cm (2") front and rear threshold. Ideal for supporting air filled cushions.

Standard comfort cushion



This cushion has a 12mm (1/2") beech wood base with a slightly contoured foam cushion, featuring a 2cm (3/4") seat angle designed to aid client in front transfers.

Optimum comfort cushion



The Optimum comfort cushion provides additional contouring and support for increased user positioning. A built-inommel, thigh guides and a seat well for optimum comfort and positioning.

Shallow contour backrest Hard back, vertical strap or tension adj. straps



Anatomically shaped with a shallow contour and padded for extra comfort. Available as a solid wooden backrest, vertical strap option or with tension adjustable straps for posture management. Solid back height 50cm, vertical strap back and tension adjustable height 57cm.

Deep contour backrest - active



This back features tension adjustable straps for posture management and a 4" deep contour providing integrated lateral support around the lower back, tapering away at the shoulders so as not to restrict active users. - 57cm back height

Deep contour backrest - passive



Similar to the Deep Contour active backrest, this features tension adjustable straps and a 4" deep contour continuing up the whole of the backrest, ideal for clients with restricted movement. - 57cm back height