

Invoice to:
 Name:
 Street:
 Town: Postal code:
 Country:
 Tel: Fax:
 Order date:
 Purchase order:

Deliver to:
 Name:
 Street:
 Town:
 Country:
 Tel:
 Marked for:

Standard Features		
● Power double tilt module	● Chest fixation (2-point strap)	● Swing-away control mount
● 0° - 22° Posterior Power tilt	● Single post armrest	● LED lights and indicators
● 0° - 80° Anterior Power tilt	● Standard full length armpads	● R-Net complex seating / lighting
● Powered lift 20cm	● R-Net 120A Controller	● 6kph (4mph) 4 pole motors
● Powered centre mount legrest (incl. Knee fixation)	● R-Net LCD Joystick	● Pneumatic tyres with V tread
● 0° to 85° powered recline with 100/110/120 mm anti-shear mechanism	● Matt black base & frame	● 60Ah batteries
		● 8A Charger

Fax: +44 (0)845 605 66 89

E-mail: orders@sunmed.co.uk

Spare Parts: <http://parts.sunrisemedical.eu>

JIVE UP: £17,105

Max. user weight: 120kg (19 stone)

SEAT DIMENSION JAY COMFORT

Seat frame width is 41 cm. The Seat width can be adjusted between 40 and 52 cm via the armrests, based on the seat/back cushion configuration. The Seat depth can be adjusted between 41 and 55 cm via backrest movement

Seat dimensions also reflect aluminium seat plate option. Seat plate option to be selected at the Jay Comfort Cushion section (JIV020109).

Cushion width

cm	40	42	44	46	48	50		
Seat depth	40	○ JUP020001	○ JUP020008	○ JUP020015	○ JUP020022	○ JUP020029	○ JUP020036	Seat cushion depth 48 cm
	42	○ JUP020002	○ JUP020009	○ JUP020016	○ JUP020023	○ JUP020030	○ JUP020037	
	44	○ JUP020003	○ JUP020010	○ JUP020017	○ JUP020024	○ JUP020031	○ JUP020038	
	46	○ JUP020004	○ JUP020011	○ JUP020018	○ JUP020025	○ JUP020032	○ JUP020039	
	48	○ JUP020005	○ JUP020012	○ JUP020019	○ JUP020026	○ JUP020033	○ JUP020040	
	50	○ JUP020006	○ JUP020013	○ JUP020020	○ JUP020027	○ JUP020034	○ JUP020041	Seat cushion depth 55 cm
	52	○ JUP020007	○ JUP020014	○ JUP020021	○ JUP020028	○ JUP020035	○ JUP020042	

SEAT HEIGHT

Seat height

cm	46	48
Jive M	○ SH46 (1)	○ SH48

Seat height measured floor to seat pan at 0° tilt and without seat cushion

(1) Not available with 52cm seat depth setting

SEAT TO BASE INTERFACE

JUP130010	Power stand-up modules	
	● Power double tilt module	(counts as 2 actuators)
	● 0° - 22° Posterior Power tilt	
	● 0° - 80° Anterior Power tilt	
	● Powered lift 20cm	
	● 0° to 85° powered recline with 100/110/120 mm anti-shear mechanism	
	● Powered centre mount legrest (incl. Knee fixation)	

JAY CUSHION SIZE**SEAT WIDTH****Size Configuration - COMPLETE PRODUCT**

		Width Coding JL 40 --- JL 42 --- JL 44 --- JL 46 --- JL 48 --- JL 50 ---					
Depth Coding	cm	40	42	44	46	48	50
JL __ 40	40	○	○	○	○(1)	○(1)	■
JL __ 42	42	○	○	○	○	○	■
JL __ 44	44	○	○	○	○	○	■
JL __ 46	46	○	○	○	○	○	○
JL __ 50	50	■	■	■	○	○	○

Conversion Table

cm	40.0	42.0	44.0	46.0	48.0	50.0	56.0
inches	15.6	16.4	17.2	18.0	18.8	19.5	21.9

(1) Not available with Soft Combi P or Easy Fluid Cushion

JAY CUSHION OPTIONS

JAY020004	<input type="checkbox"/>	Jay Lite Standard version	(Micro climatic cover, Aluminium seat plate to be selected)	215
JAY020005	<input type="checkbox"/>	Jay Easy Fluid with Flat Base	(Incontinence cover, Aluminium seat plate to be selected)	308
JAY020001	<input type="checkbox"/>	Jay Easy Visco with Flat Base	(Incontinence cover, Aluminium seat plate to be selected)	132
JAY020003	<input type="checkbox"/>	Jay Soft Combi P	(Incontinence cover, Aluminium seat plate to be selected)	77
JUP020109	<input type="radio"/>	Seat Plate Aluminium	(supplied in the size of the cushion chosen)	

JAY COMFORT SEAT CUSHION OPTIONS

JCL020003	<input type="checkbox"/>	Standard cushion - stretched rubber upholstery	(thickness 8cm)	155
JCL020007	<input type="checkbox"/>	Optimum comfort - stretched rubber upholstery	(anatomically formed) (thickness 9cm)	199
JUP020110	<input type="radio"/>	Omit cushion		

JAY COMFORT

JCL020003



JCL020007

JAY COMFORT BACK CUSHION OPTION

(57cm back height - Adj. up to 70cm)









JCL030004	<input type="checkbox"/>	Shallow contour (standard velour)	(57 cm back height)	(hard back)	153
JCL030008	<input type="checkbox"/>	Shallow contour (standard velour)	(57 cm back height)	(vertical elastic straps)	222
JCL030006	<input type="checkbox"/>	Shallow contour (standard velour)	(57 cm back height)	(horizontal tension adjustable straps)	222
JCL030010	<input type="checkbox"/>	Deep contour - active	(57 cm back height)	(horizontal tension adjustable straps)	257
JCL030012	<input type="checkbox"/>	Deep contour - passive	(57 cm back height)	(horizontal tension adjustable straps)	257
JUP030209	<input type="radio"/>	Omit back upholstery	(Chair supplied with NO back upholstery option)		
JUP030210	<input type="checkbox"/>	Jay Comfort mounting hardware only			39






JAY COMFORTShallow Contour
JCL030004 /
JCL030008 /
JCL030006Deep Contour Active
JCL030010Deep Contour
Passive
JCL030012**COMFORT BACK CUSHION WIDTH**


BW40	<input type="radio"/>	40cm (15.5")	BW42	<input type="radio"/>	42cm (16.5")	BW44	<input type="radio"/>	44cm (17")
BW46	<input type="radio"/>	46cm (18")	BW48	<input type="radio"/>	48cm (19")	BW50	<input type="radio"/>	50cm (20")

POSITIONING STRAPS

JUP090047	<input type="radio"/>	5 cm (2") aircraft buckle				
JUP090048	<input type="checkbox"/>	5 cm (2") padded aircraft buckle				105
JUP090052	●	Chest fixation (2-point strap)	JUP090054	M	L	JUP090055
JUP090053	<input type="checkbox"/>	2nd Chest fixation (2-point strap)	JUP090056	M	L	JUP090057
JUP090097	<input type="checkbox"/>	Adjustable and lift off chest bar with upholstered cover				165
JUP090098	<input type="checkbox"/>	Adjustable and lift off chest bar with chest pad (20cm x 10cm)				165

										£	
HEADRESTS											
JUP030426	<input type="radio"/>	Curved headrest	20cm x 11cm (8" x 4.5")								
JUP030417	<input type="radio"/>	Curved headrest	23cm x 11cm (9" x 4.5")								
JUP030427	<input type="checkbox"/>	Contoured (winged) headrest								105	
JUP030428	<input type="checkbox"/>	Jay Comfort: standard comfort seating headrest								138	
JUP030429	<input type="checkbox"/>	Jay Comfort: flexible comfort seating headrest (winged style)								165	
JUP030415	<input type="radio"/>	Omit headrest									
	Curved Headrest JUP030426 / JUP030417		Contoured Headrest JUP030427		JAY Comfort: Std Comfort headrest JUP030428		JAY Comfort: Winged headrest JUP030429				
BACKREST RECLINE OPTIONS											
JUP030204	<input type="checkbox"/>	Height adjustable push handles								97	
JUP030601	<input type="checkbox"/>	Height adjustable swing-away laterals (short)								259	
ARMREST TYPE											
JUP040122	<input type="radio"/>	Flip Back armrest								<i>(width and height adjustable, pair)</i>	
JUP040138	<input type="checkbox"/>	Reclining and flip-back armrest								<i>(width and height adjustable, pair)</i> 334	
SKIRT GUARDS <i>(not available with reclining armrest)</i>											
JUP040139	<input type="checkbox"/>	Jay Comfort padded skirt guard	36cm (14")							105	
LATERAL SEAT SUPPORTS											
JUP020340	<input type="checkbox"/>	Adjustable seat side rails								88	
Adjustable Side pads <i>(requires JUP020340 selection)</i>								Left	Right		
<input type="checkbox"/>	Knee laterals	JUP020341	<input type="radio"/>	<input type="radio"/>	JUP020342			105			
<input type="checkbox"/>	Upper leg laterals	JUP020343	<input type="radio"/>	<input type="radio"/>	JUP020344			105			
Adjustable and swing away chest lateral supports											
JUP030601	<input type="checkbox"/>	Height adjustable Swing away laterals, short								259	
ARMPADS - LEFT AND RIGHT HAND SIDE											
	Standard armpads		Left Right								
	Full length	42 cm (16,5")	JUP040201	<input type="radio"/>	<input type="radio"/>	JUP040202					
	JAY COMFORT armpads										
	Padded (vinyl) upholstery	full length 40 cm (16")	JUP040203	<input type="checkbox"/>	<input type="checkbox"/>	JUP040204	51 each				
	P-Shaped skai upholstery	full length 45 cm (18")	JUP040205	<input type="checkbox"/>	<input type="checkbox"/>	JUP040206	51 each				
	Hemi Armpads		<i>(Hemi Armpad & hand pad must be chosen together)</i>								
	Hemi forearm pad	small 23cm (9")	JUP040207	<input type="checkbox"/>	<input type="checkbox"/>	JUP040208	85 each				
	Hemi forearm pad	medium 28cm (11")	JUP040209	<input type="checkbox"/>	<input type="checkbox"/>	JUP040210	85 each				
	Flat hand pads	small 11 x 18cm (4.5" x 7")	JUP040211	<input type="checkbox"/>	<input type="checkbox"/>	JUP040212	85 each				
	Flat hand pads	large 14 x 23cm (5.5" x 9")	JUP040213	<input type="checkbox"/>	<input type="checkbox"/>	JUP040214	85 each				
R-NET 120A CONTROLS											
JUP110004	<input checked="" type="radio"/>	R-net 120A control system								<i>(no on board actuator control output)</i>	
JUP110017	<input checked="" type="radio"/>	R-net LCD colour joystick								<i>(Integrated colour display with 2 jack inputs for buddy buttons)</i>	
Dual control options											
JUP110021	<input type="checkbox"/>	Dual control								<i>(Separate dual control module mounted at rear - Always RH)</i> 425	
JOYSTICK MOUNTINGS											
Front control mountings <i>(for single Joystick selection mount)</i>											
JUP110051	<input type="radio"/>	Swing-away									
JUP110143	<input type="checkbox"/>	Swing-away tray mounted centre control - R-Net LCD								<i>(not available on 22" width)</i> 1559	
Joystick mounting side											
JUP110053	Left <input type="radio"/>	<input type="radio"/>	Right	JUP110052							

				£	
5-Switch direct actuator box (CTRL+5) (c/w assignable Button functionality)					
JUP110176	<input type="checkbox"/>	5-Switch direct actuator box with paddle switches	(direct control of 6 memory seat positions)	550	
JUP110175	<input type="checkbox"/>	5-Switch direct actuator box with buttons	(direct control of 6 memory seat positions)	660	
LIGHTS & INDICATORS					
JIV110201	<input type="radio"/>	LED (StVZO approved)			
R-NET OMNI DISPLAY & SPECIAL CONTROL INPUT MODULE (with IR functionality)					
JUP110101	<input type="checkbox"/>	R-net OMNI colour display input module	(includes Special Control Interface)	850	
JUP110102	<input type="checkbox"/>	OMNI in scanning mode	(c/w 10 assignable functions through DB9 interface)	850	
JUP110110	<input type="radio"/>	R-NET OMNI Flexi mount only	(30 cm (12") long flexible height extension)		
		JUP110121 Left <input type="radio"/>	<input type="radio"/> Right JUP110122		
R-NET OPTIONAL CONTROLS					
JUP110033	●	R-NET CxSM (seating module with memory function)			
JUP110031	●	R-NET ISM-L (lighting module)			
JUP090068	<input type="checkbox"/>	R-net Bluetooth module	(Compatible with Windows PC and Android systems only)	334	
JUP090077	<input type="checkbox"/>	R-net 2nd Bluetooth module	(Compatible with Windows PC and Android systems only)	334	
JUP090078	<input type="checkbox"/>	R-net Input/Output module	(7 Switched controllable output commands, only available via B4Me)	385	
JUP090069	<input type="checkbox"/>	Bluetooth wireless mouse receiver		63	
R-NET SPECIAL CONTROLS (R-net OMNI required from previous page; 2 buddy buttons recom.)					
JUP110130	<input type="checkbox"/>	Switch-It Versaguide Joystick as hand control		1643	
JUP110131	<input type="checkbox"/>	Switch-It Micro Guide as hand control		2210	
JUP110132	<input type="checkbox"/>	Switch-It Micro Pilot (0 way) joystick as hand control		2437	
JUP110150	<input type="radio"/>	Fixed front mount			
JUP110151	<input type="radio"/>	Swing-away front mount			
JUP110162	<input type="radio"/>	Left <input type="radio"/>	Right <input type="radio"/> JUP110161		
	R-Net OMNI in Scanning Mode JIV100550		Proportional Versaguide Joystick JIV100551		Switch-It Micro Guide Joystick JIV100552
					Switch-It Micro Pilot 0-Way Joystick JIV100553
R-NET HEAD RELATED INPUT DEVICES					
(if you have already chosen a standard headrest then Switch- It head array is not possible)					
JUP110133	<input type="checkbox"/>	Switch-It 3 switch head array / adult	(5 switch and other versions available via B4M) (OMNI required, 2 BB recommended)	3343	
JUP110138	<input type="checkbox"/>	Switch-It proportional head control	(OMNI required, 2 BB recommended)	3116	
	Switch-It Head Array JUP110133				
Buddy button switches					
JUP110181	<input type="checkbox"/>	Red x 1		85	
Piko button switches					
JUP110191	<input type="checkbox"/>	Red x 1	JUP110192 <input type="checkbox"/> Red x 2	85	
JUP110193	<input type="checkbox"/>	Green x 1	JUP110194 <input type="checkbox"/> Green x 2	85	
JUP110195	<input type="checkbox"/>	Yellow x 1	JUP110196 <input type="checkbox"/> Yellow x 2	85	
JUP110197	<input type="checkbox"/>	Blue x 1	JUP110198 <input type="checkbox"/> Blue x 2	85	
SPEED					
JIV100710	<input type="radio"/>	6 kph (4mph) (4 - pole)			
JIV100711	<input type="radio"/>	10 kph (6mph) (4 - pole)			
JIV100712	<input type="checkbox"/>	13 kph (8mph) (4 - pole)	(lower torque than 10kph motors)	344	

● = Standard ○ = No cost option □ = Cost option			£
DRIVE WHEELS - JIVE M always with 6" solid castor wheel on front and rear			
JIV070032	○	14" Pneumatic "V-type"	
JIV070031	○	14" Solid "V-type"	
JIV070034	○	14" Pneumatic "knobby"	
JIV070033	○	14" Solid "Knobby"	
BATTERIES			
JIV160001	○	60 Ah maintenance free	
JIV160005	□	80 Ah maintenance free <i>(Recommended for 13Kph Motors where battery power is drawn faster)</i>	111
CHARGER			
JIV160100	○	220V/ 8A automatic charger (UK)	
JIV160102	□	220V/ 10A automatic charger (UK) <i>(Must be ordered with 13kph speed - charges batteries faster)</i>	168
SHROUD COLOUR - castor arms are always black			
JIV140001	○	White RAL 9010	
JIV140002	○	Red	
JIV140003	○	Blue	
JIV140004	○	Matt Black	
ACCESSORIES			
JUP090067	○	Rear view mirror <i>(Always required with 10Kph Chairs)</i>	
JUP090040	□	Swing-away tray table	164
JUP090061	□	24/12V DC/ DC converter	221
JOYSTICK KNOBS			
JUP090505	□	Golf ball joystick knob	29
JUP090506	□	T bar B38 joystick knob	29
JUP090507	□	Soft ball joystick knob	29
JUP090508	□	Ball joystick knob	29
JUP090509	□	Chin cup joystick knob	29
JUP090510	□	Stick joystick knob	29
JUP090511	□	Mushroom joystick knob	29
For special requirements, please complete Built-4-Me prescription form			
			
JIVE - base price			£17,105
Option Price			£
Total Price			£



J3 Back

J3 Backrest Range

Complete backrest package for clients up to 136kg
 Universal one point fixation for rapid application
 Easily adjustable for assessment & changing needs
 Wide range of positioning and transfer options
 Backs and options fully transit tested
 1 Large & 1 small lumbar included in MT, UT & SH box
 1 small lumbar in LT box 2 year guarantee



- = 5 days from the EDC
- ◇ = 15 days from the EDC
- = upcharge option

Max. user weight: **136kg**
USER WEIGHT kg
 Mandatory information to ensure correct configuration

Applications:

For all clients requiring postural back support from mild lateral stability to extreme lateral and posterior stability.
 Full range of postural back shapes to cover clinical back support at the appropriate backrest height for the user.
 Adjustable in depth, height, angle and chair width, ensuring rapid and comprehensive client positioning and spine alignment.

CLIENT MEASUREMENT

1. Assess Lateral Stability Requirement (SC, MC, DC or PDC)
2. Assess Posterior Stability Requirement (MT; UT or SH)
3. Assess appropriate backrest height (Short, Medium, Tall)
4. Determine client trunk width across the sternum; and cross check with backpost width.
5. Assess for options and 'Velcro' positioning inserts

Item Nr. construction:
J3 SC LT S 41 (This is only an example)
J3: J3 Back
SC: Shallow Contour
LT: Lower Thoracic
S: Short (Height cm)
41: Chair Backrest Width (cm)



Shallow Contour / SC Moderate Lateral Support	Mid Contour / MC Medium Lateral Support	Deep Contour / DC Full Lateral Support	Posterior Deep Contour / PDC PDC Full Lateral support
---	---	--	---


Width
 Chair backrest widths (Measured from outside of backrest to outside of backrest)

Jay 3 Backrest Width											
	41	46	51		41	46	51		41	46	51
	41 - 46	46 - 51	51 - 56		41 - 46	46 - 51	51 - 56		41 - 46	46 - 51	51 - 56

Back height



Backrest Height		£520			£520			£570			£570		
Upper Thoracic (UT)	(S) - 42	38	43	48	34.5	39.5 (1)	44.5	33	38 (1)	43	33	38	43
	(M) - 46	○	○	○	○	○	○	○	○	○	○	○	○
	(T) - 50	○	○	○	○	○	○	○	○	○	○	○	○



Backrest Height		£520			£520			£570			£570		
Shoulder Height (SH)	(S) - 53	38	43	48	34.5	39.5 (1)	44.5	33	38 (1)	43	33	38	43
	(M) - 57	○	○	○	○	○	○	○	○	○	○	○	○
	(T) - 61	○	○	○	○	○	○	○	○	○	○	○	○

(1) Not available with laterals

Conversion Table

cm	17.0	20.0	24.0	30.0	34.0	36.0	38.0	41.0	42.0	44.0	46.0	51.0	53.0	56.0	61.0
inches	6.5	8.0	9.5	12.0	13.5	14.0	15.0	16.0	16.5	17.5	18.0	20.0	21.0	22.0	24.0

 <p>JAY090056 & JAY090057</p>	 <p>JAY090006</p>	 <p>JAY090007</p>
 <p>JAY090008</p>	 <p>JAY090009</p>	 <p>JAY090010</p>
 <p>Fixed laterals</p>	 <p>Swing-away laterals</p>	 <p>JAY090029 - JAY090031</p>
 <p>JAY090032 - JAY090034</p>	 <p>JAY090035</p>	 <p>JAY090036</p>
 <p>JAY090037 - JAY090039</p>	 <p>JAY090040 - JAY090042</p>	 <p>JAY090043 - JAY090045</p>
 <p>JAY090046 - JAY090048</p>	 <p>JAY090049 - JAY090051</p>	 <p>JAY090052</p>
 <p>JAY090053</p>	 <p>JAY090054</p>	 <p>JAY090055</p>



JUP090047 & JUP090048



JCL020003




JCL020007




JCL030004




JCL030008



JCL030006



JCL030012



JCL030010



JUP030426 & JUP030417



JUP030427



JUP030428



JUP030429



JUP040201 & JUP040202



JUP040203 & JUP040204



JUP040205 & JUP040206



JUP040211 - JUP040214



JUP040207 - JUP040210



JUP110176



JUP110175



JUP090068



JUP110102



JUP110130



JUP110131



JUP110132



JUP110151



JUP110133



JUP110138



JUP110191 - JUP110198



JUP110181



V-Type




Knobby



JIV160101



JIV160103



JUP090067



JUP090061



JUP090505



JUP090506



JUP090507



JUP090508



JUP090509




JUP090510



JUP090511



JUP090040



JUP090069