

Invoice to:	
Name:	
Street:	
Town:	Postal code:
Country:	
Tel:	Fax:
Order date:	
Purchase order:	



Deliver to:	
Name:	
Street:	
Town:	Postal code:
Country:	
Tel:	Fax:
Marked for:	

Fax: +44 (0)845 605 66 89

Base Price £11,550
E-mail: orders@sunmed.co.uk

Max. user weight: 160 Kg

Standard Features	<ul style="list-style-type: none"> ● C-Me Drive Technology ● R-Net 120Ah Controller
<ul style="list-style-type: none"> ● Jive M2 Matt Black Chassis ● SpiderTrac 2 Suspension ● 60Ah Batteries ● 6 or 10Kph 4 Pole Motors 	<ul style="list-style-type: none"> ● Sedeo Ergo Biometric Seating System ● 180° Anti Shear Backrest ● Powered centre mount legrest ● Headrest
	<ul style="list-style-type: none"> ● R-Net Advanced Color LCD Joystick ● Capsulated Swing away joystick mount ● 6 Programmable Smart seat Memory positions ● Lapbelt

SEAT DIMENSION

Effective seat width ⁽¹⁾

cm	40	42	44	46	48	50	52	54	56
40 (4)	○	○	○	○					
42 (4)	○	○	○	○					
44 (5)	○	○	○	○					
Seat depth	46	○	○	○	○	○	○		
	48	○	○	○	○	○	○		
	50			○	○	○	○	○	○
	52					○	○	○	○
	54						○	○	○
	56							○	○

(2)

The seat width can be adjusted between 40 and 56 cm via the seat plate and the armrests.

The seat depth can be adjusted between 40 and 56 cm via backrest movement

Note: At 40cm seat depth Backrest fully reclines does not lie flat. Stops at 160°
SEAT HEIGHT

Seat height (cm)				
Seat to base interface	41	43.5	46	48.5
Filler module	○	○	○	○
50° Power tilt		○ (3)	○ (3)	○
-5° to 45° Powered tilt option				○
Power seat lift and tilt combined			○	○

Seat height measured floor to seat pan at 0° tilt

(1) Physical cushion width = always 4 cm smaller than effective seat width

(2) Not available with 50° tilt module

(3) Not available with 50, 52 and 54cm seat depth

(4) Maximum recline angle 160°

(5) Maximum recline angle 170°

SEAT TO BASE INTERFACE

JMV130001	<input type="radio"/>	Filler module			(tool adjustable seat angle)	
		JMV130100 <input type="radio"/> 0°	JMV130106 <input type="radio"/> 6°			
		JMV130103 <input type="radio"/> 3°	JMV130109 <input type="radio"/> 9°			

Power seating modules (with positioning feedback)

JMV130007	<input type="checkbox"/>	50° Power tilt (1)	(0° to 50°, with centre of gravity shift)	715
JMV130005	<input type="checkbox"/>	-5° to 45° Powered tilt option	(3 cm added seat to floor height at 0°)	765
JMV130006	<input type="checkbox"/>	Power seat lift and tilt combined	(50° tilt and 30 cm lift)	1760

SEDEO ERGO SEAT CUSHION OPTIONS

JMV020501	<input type="radio"/>	Cushion - Air flow upholstery (1)	(Standard option if nothing selected)	
JMV020502	<input type="checkbox"/>	Cushion - Air flow upholstery with incontinence protection (1)		105
JMV020509	<input type="radio"/>	Centre, side cushion and hip support	(Standard option if nothing selected)	
JMV020510	<input type="checkbox"/>	Omit center cushion - Select if JAY cushion required WITH Sedeo Ergo Side supports		-55
JMV020512	<input type="checkbox"/>	Omit center, side cushion and hip support		-110
JMV020513	<input type="checkbox"/>	JAY Cushion - select if JAY Cushion required		-110

		Left	Right			
<input type="checkbox"/>	Lateral s/a knee support	JMV020341	<input type="radio"/>	<input type="radio"/>	JMV020342	160 each
<input type="checkbox"/>	Lateral Thigh Support (2)	JMV020343	<input type="radio"/>	<input type="radio"/>	JMV020344	160 each

SEDEO ERGO BACK CUSHION OPTION

Back cushion height	BH51 <input type="radio"/> 51cm	BH56 <input type="radio"/> 56cm	BH61 <input type="radio"/> 61cm	
JMV030051	<input type="radio"/>	Back Upholstery - AIR FLOW		
JMV030209	<input type="radio"/>	Omit back upholstery		

SEDEO ERGO BACK CUSHION WIDTH

BW40	<input type="radio"/>	40 cm	BW46	<input type="radio"/>	46 cm	BW52	<input type="radio"/>	52 cm	
BW42	<input type="radio"/>	42 cm	BW48	<input type="radio"/>	48 cm	BW54	<input type="radio"/>	54 cm	
BW44	<input type="radio"/>	44 cm	BW50	<input type="radio"/>	50 cm	BW56	<input type="radio"/>	56 cm	

POSITIONING STRAPS*(anterior pelvic supports)*

JMV090047	<input type="radio"/>	5 cm aircraft buckle (130 cm length)	
JMV090056	<input type="radio"/>	5 cm aircraft buckle (190 cm length)	
JMV090048	<input type="checkbox"/>	5 cm padded aircraft buckle	55

SEDEO ERGO HEADRESTS

JMV030406	<input type="radio"/>	Sedeo Ergo contoured headrest with air flow upholstery (20 x 11 cm)	
JMV030407	<input type="radio"/>	Sedeo Ergo contoured headrest with air flow upholstery (23 x 11 cm)	
JMV030408	<input type="checkbox"/>	Sedeo Ergo flexible comfort headrest with adjustable side wings	105
JMV030415	<input type="checkbox"/>	Omit Headrest	-55

BACKREST RECLINE OPTIONS

JMV030034	<input checked="" type="radio"/>	85° to 175° Power recline c/w mechanical shear reduction (3)	
JMV030201	<input type="checkbox"/>	Push Handles	138

(1) please use the cushion composition**(2) Only available with JMV020512; Omit center, side cushion and hip support****(3) Max. 160° at 40 and 42cm seat depth, Max. 170° with 44cm seat depth**

JAY CUSHIONS

For Jay cushion must select omit centre, side cushion and hip support- JMV020512

		Seat Width						
		40	42	44	46	48	50	56
Seat Depth	40	○	○	○	○ ⁽¹⁾	○ ⁽¹⁾		
	42	○	○	○	○	○		
	44	○	○	○	○	○		
	46	○	○	○	○	○	○	○ ⁽³⁾
	48		○ ⁽²⁾	○ ⁽²⁾	○ ⁽²⁾	○ ⁽²⁾	○ ⁽²⁾	○ ^(2/3)
	50			○ ⁽²⁾	○	○	○	○ ⁽³⁾
56							○ ⁽³⁾	

(1) Not available with Soft Combi P, Easy Fluid & Easy Visco

(2) Only available for Jay Balance cushions

(3) Not available with Jay Lite cushions

JAY020003	<input type="checkbox"/> Jay Soft Combi P	(Incontinence cover)	77
JAY020004	<input type="checkbox"/> JAY Lite Standard version	(Microclimatic cover)	215
JAY020001	<input type="checkbox"/> Jay Easy Visco with Flat Base	(Incontinence cover)	132
JAY020005	<input type="checkbox"/> Jay Easy Fluid with Flat Base	(Incontinence cover)	308
JAY020010	<input type="checkbox"/> Jay Balance Roho (Air Single Valve)	(Microclimatic cover)	495
JAY020011	<input type="checkbox"/> Jay Balance Fluid	(Microclimatic cover)	468

JAY J3 BACKS - measurements in cm



Shallow Contour (SC)			
41	46	51	



Mid Contour (MC)			
41	46	51	



Deep Contour (DC)			
41	46	51	



Posterior Deep Contour (PDC)			
41	46	51	

Jay J3 Back Size

Backrest Height	
Mid Thoracic (MT)	Short (S) - 30
	Med. (M) - 34
	Tall (T) - 38

572			
○	○	○	
○	○	○	
○	○	○	

572			
○	○	○	
○	○	○	
○	○	○	

Backrest Height	
Upper Thoracic (UT)	Short (S) - 42
	Med. (M) - 46
	Tall (T) - 50

572			
○	○	○	
○	○	○	
○	○	○	

572			
○	○	○	
○	○	○	
○	○	○	

627			
○	○	○	
○	○	○	
○	○	○	

627			
○	○	○	

Backrest Height	
Shoulder Height (SH)	Short (S) - 53
	Med. (M) - 57
	Tall (T) - 61

572			
○	○	○	
○	○	○	
	○	○	

572			
○	○	○	
○	○	○	
	○	○	

627			
○	○	○	
○	○	○	
	○	○	

627			
○	○	○	

JAY J3 Back Item Nr. construction: J3 SC MT S 41 (This is only an example)

J3: J3 Back

SC: Shallow Contour

MT: Mid Thoracic

S: Short (Height cm)


41: Chair Backrest Width (cm)

JAY J3 BACK MOUNTING HARDWARE				£																																																						
				Unit Price																																																						
JAY090065	<input type="checkbox"/>	Jay 3 Back Shell Adaptor																																																								
JAY J3 BACK HEADREST (only available where headrest is indicated on the first page)				Crash tested																																																						
				Unit Price																																																						
Please choose 1 mount and 1 pad to have a complete headrest with pad																																																										
Mount																																																										
JAY090006	<input type="checkbox"/>	Jay Standard Headrest	(Bracket & Mount only)	88																																																						
JAY090007	<input type="checkbox"/>	Jay Flip Back Headrest	(Bracket & Mount only)	198																																																						
Pad																																																										
JAY090008	<input type="checkbox"/>	Jay Standard Pad		41																																																						
JAY090009	<input type="checkbox"/>	Jay Contour Pad		41																																																						
JAY090010	<input type="checkbox"/>	Universal Headrest Mount	(required to use any other headrest not listed in this order form)	41																																																						
JAY J3 BACK LATERAL SUPPORT (only available where laterals are indicated on the first page)				Unit Price																																																						
<p>Left Hand lateral</p> <table border="1"> <thead> <tr> <th rowspan="2">Height</th> <th colspan="3">Depth</th> <th rowspan="2"></th> <th rowspan="2"></th> </tr> <tr> <th>10</th> <th>12</th> <th>15</th> </tr> </thead> <tbody> <tr> <td>10</td> <td><input type="checkbox"/> JAY090011</td> <td><input type="checkbox"/> JAY090012</td> <td><input type="checkbox"/> JAY090013</td> <td>←</td> <td>Fix lateral JAY090017</td> </tr> <tr> <td>15</td> <td><input type="checkbox"/> JAY090014</td> <td><input type="checkbox"/> JAY090015</td> <td><input type="checkbox"/> JAY090016</td> <td>←</td> <td>Fix lateral JAY090018</td> </tr> <tr> <td></td> <td><input type="checkbox"/> JAY090019</td> <td></td> <td></td> <td>←</td> <td>Swing Away option</td> </tr> </tbody> </table> <p>Right Hand lateral</p> <table border="1"> <thead> <tr> <th rowspan="2">Height</th> <th colspan="3">Depth</th> <th rowspan="2"></th> <th rowspan="2"></th> </tr> <tr> <th>10</th> <th>12</th> <th>15</th> </tr> </thead> <tbody> <tr> <td>10</td> <td><input type="checkbox"/> JAY090020</td> <td><input type="checkbox"/> JAY090021</td> <td><input type="checkbox"/> JAY090022</td> <td>←</td> <td>Fix lateral JAY090026</td> </tr> <tr> <td>15</td> <td><input type="checkbox"/> JAY090023</td> <td><input type="checkbox"/> JAY090024</td> <td><input type="checkbox"/> JAY090025</td> <td>←</td> <td>Fix lateral JAY090027</td> </tr> <tr> <td></td> <td><input type="checkbox"/> JAY090028</td> <td></td> <td></td> <td>←</td> <td>Swing Away option</td> </tr> </tbody> </table>					Height	Depth					10	12	15	10	<input type="checkbox"/> JAY090011	<input type="checkbox"/> JAY090012	<input type="checkbox"/> JAY090013	←	Fix lateral JAY090017	15	<input type="checkbox"/> JAY090014	<input type="checkbox"/> JAY090015	<input type="checkbox"/> JAY090016	←	Fix lateral JAY090018		<input type="checkbox"/> JAY090019			←	Swing Away option	Height	Depth					10	12	15	10	<input type="checkbox"/> JAY090020	<input type="checkbox"/> JAY090021	<input type="checkbox"/> JAY090022	←	Fix lateral JAY090026	15	<input type="checkbox"/> JAY090023	<input type="checkbox"/> JAY090024	<input type="checkbox"/> JAY090025	←	Fix lateral JAY090027		<input type="checkbox"/> JAY090028			←	Swing Away option
Height	Depth																																																									
	10	12	15																																																							
10	<input type="checkbox"/> JAY090011	<input type="checkbox"/> JAY090012	<input type="checkbox"/> JAY090013	←	Fix lateral JAY090017																																																					
15	<input type="checkbox"/> JAY090014	<input type="checkbox"/> JAY090015	<input type="checkbox"/> JAY090016	←	Fix lateral JAY090018																																																					
	<input type="checkbox"/> JAY090019			←	Swing Away option																																																					
Height	Depth																																																									
	10	12	15																																																							
10	<input type="checkbox"/> JAY090020	<input type="checkbox"/> JAY090021	<input type="checkbox"/> JAY090022	←	Fix lateral JAY090026																																																					
15	<input type="checkbox"/> JAY090023	<input type="checkbox"/> JAY090024	<input type="checkbox"/> JAY090025	←	Fix lateral JAY090027																																																					
	<input type="checkbox"/> JAY090028			←	Swing Away option																																																					
					46																																																					
					46																																																					
					87																																																					
					46																																																					
					46																																																					
					87																																																					
When ordering 2 laterals Left and Right, you need to pick the sizes in both matrices above.				Total Price																																																						
If a swing-away lateral is required, first select the height and the depth, then tick the Swing Away option box below the matrix																																																										
JAY J3 BACK POSITIONING KIT				CM width	Unit Price																																																					
				41 cm	46 cm																																																					
				51 cm																																																						
Spine Align Mini Fit Kit I				<input type="checkbox"/>	31																																																					
Spine Align Mini Fit Kit II				<input type="checkbox"/>	47																																																					
JAY090035	<input type="checkbox"/>	Jay 3 Back Positioning Box		151																																																						
JAY090036	<input type="checkbox"/>	Jay 3 Back Hardware Box		474																																																						
JAY J3 BACK CHEST STRAP				Medium	Large																																																					
				Extra large	Unit Price																																																					
Padded Anterior Strap				<input type="checkbox"/>	<input type="checkbox"/>																																																					
				JAY090037	JAY090038																																																					
				JAY090039	157																																																					
Centre open harness (4 point strap, not available with Jay Comfort)				<input type="checkbox"/>	<input type="checkbox"/>																																																					
				JAY090040	JAY090041																																																					
				JAY090042	157																																																					
JAY J3 BACK SPARES (separate to product - separate box or delivery)				CM width	Unit Price																																																					
				41 cm	46 cm																																																					
				51 cm																																																						
Small Lumbar (1 unit included as std with product LT, MT, UT & SH)				<input type="checkbox"/>	9																																																					
Large Lumbar Support (1 unit included as std with product MT, UT & SH)				<input type="checkbox"/>	9																																																					
Wedge				<input type="checkbox"/>	9																																																					
				Unit price	Quantity																																																					
JAY090052	<input type="checkbox"/>	Curved Shim		9																																																						
JAY090053	<input type="checkbox"/>	Flat Shim		9																																																						
JAY090054	<input type="checkbox"/>	Posterior/Lateral Pelvic		12																																																						
JAY090055	<input type="checkbox"/>	Lateral Shim		4																																																						
JAY090056	<input type="checkbox"/>	Small Bean Bag		9																																																						
JAY090057	<input type="checkbox"/>	Large Bean Bag		9																																																						


										£	
ARMREST TYPE										(pair)	
JMV040123	●	Reclining and flip up armrest							(reclining and flip-back)		
JMV040142	□	Fixed Sideguard with Upholstery (1)							(width and height adjustable via seat plate)		165
ARMPADS - LEFT & RIGHT HAND SIDE											
Standard armpads (Microclimatic cover to match Ergo Upholstery)										Left	Right
Comfort armpads		41 cm	JMV040217	○	○	JMV040218					
Comfort armpads											
Skai (vinyl) upholstery		full length	40 cm	JMV040203	□	□	JMV040204		52 each		
P-Shaped skai upholstery		full length	45 cm	JMV040205	□	□	JMV040206		52 each		
Hemi Arm pads											
Hemi forearmpad		small	23 cm (2)	JMV040207	□	□	JMV040208		87 each		
Hemi forearmpad		medium	28 cm (2)	JMV040209	□	□	JMV040210		87 each		
Hemi forearmpad		large	33 cm (2)	JMV040219	□	□	JMV040220		87 each		
Hemi hand pads (required when forearm pad selected)											
Flat hand pads		small	11 x 18 cm	JMV040211	□	□	JMV040212		87 each		
Flat hand pads		large	14 x 23 cm	JMV040213	□	□	JMV040214		87 each		
LEGRESTS											
JMV050111	●	Powered centre mount legrest							Lower leg length 34 cm - 52 cm		
JMV050038	○	Adjustable calf pads with vent air cover									
R-NET CONTROLS / ELECTRONICS											
JMV110004	●	R-NET Heavy Duty motor control (120 A)							(no on board actuator control output)		
JMV110009	□	R-NET Traction Control System (GYRO Control System)								275	
JMV110013	●	Advanced Colour LCD JSM							(with 2 jack inputs, actuator and L&I control)		
JMV110021	□	Dual control (Must also select push handles)							(Separate Dual control module mounted at rear on RH mounted)		438
JMV110020	□	R-Net JSM control ATTENDANT ONLY							(Select JSM above to be used RH mounted)		105
Dual/Attendant control options mounting position				Left	Right						
				○	○	JMV110058					
R-NET CONTROL MOUNTINGS											
Front control mountings										(for single Joystick selection mount)	
JMV110056	●	Capsulated Swing away									
JMV110143	□	Centre tray mounted LCD remote control - Supplied with R-Net Std LCD Joystick (NOT Advanced Joystick)								1606	
				Left	Right						
				○	○	JMV110053					
Actuator Switch Boxes (Comes Standard with USB charger Port)											
JMV110175	□	5-Switch Box with push buttons with free assignable functions								550	
JMV110176	□	5-Switch Box with Toggle Switches with free assignable functions								660	
5-Switch Box mounting position				Left	Right						
				○	○	JMV110078					
LIGHTS & INDICATORS											
JMV110201	●	LED with StVZO approval									
(1) Only available with JMV020512; Omit centre, side cushion and hip support											
(2) Must also select flat hand pad											

								£
R-NET OPTIONAL CONTROLS / ELECTRONICS								
JMV110031	●	R-NET ISM-L (lighting module)						
JMV110033	●	R-NET CxSM (8-channel seating module) <i>(includes memory seat programming)</i>						
JMV090068	□	R-NET BLUETOOTH MODULE (Controls 1 x Android and Windows device)						215
JMV090141	□	2nd R-NET BLUETOOTH MODULE (Controls second Android and Windows device)						215
JMV090142	□	R-NET BLUETOOTH MODULE IOS (Controls Apple Device) <i>(Not available with 2nd Android Bluetooth module)</i>						215
OMNI DISPLAY & SC INPUT MODULE <i>with IR functionality</i>								
JMV110101	□	R-NET OMNI colour display input module <i>(includes Special Control Interface)</i>						876
JMV110110	○	Flex front control mountings <i>(If JSM selected always mounted opposite to remote pod)</i> <i>(30 cm long flexible height extension)</i>						
			Left	Right				
		JMV110121	○	○	JMV110122			
R-NET SPECIAL HAND CONTROLS <i>R-NET OMNI required; 2 buddy buttons recommended</i>								
JMV110102	□	OMNI in scanning mode <i>(c/w 10 assignable functions through DB9 interface)</i>						876
JMV110130	□	Switch-It Compact Joystick (Versa Guide) as hand control						1693
JMV110131	□	Switch-It mini joystick as hand control <i>(only fixed mounted)</i>						2276
JMV110132	□	Switch-It 0-way joystick as hand control <i>(only fixed mounted)</i>						2509
Mounting options for hand controls								
JMV110150	○	Fixed, front mount		Left	Right			
			JMV110162	○	○	JMV110161		
JMV110151	○	Swing away, front mount		Left	Right			
			JMV110162	○	○	JMV110161		
R-NET HEAD RELATED SPECIAL INPUT DEVICES								
JMV110133	□	Switch-It proximity head array / adult <i>(OMNI required, 2 BB recommended)</i>						3443
JMV110146	□	Dual Pro Switch-It proportional proximity head control <i>(OMNI required, 2 BB recommended)</i>						4345
JMV110134	□	Switch-It 0-way joystick as chin control <i>(OMNI required, 2 BB recommended)</i>						2509
JMV110135	□	Proportional chin control <i>(proportional joystick; OMNI required; 2 BB recommended)</i>						876
JMV110136	□	Sip & Puff control <i>(tube and mouthpiece; OMNI required; 2 BB recommended)</i>						293
JMV110137	□	Switch-It mini joystick as chin control <i>(OMNI required, 2 BB recommended)</i>						2276
JMV110138	□	Switch-It proportional head control <i>(OMNI required, 2 BB recommended)</i>						3210
JMV110139	□	Switch-It Compact Joystick (Versa Guide) as chin control <i>(OMNI required, switch kit incl.)</i>						1810
Mounting options for head and chin controls								
JMV110152	□	Manual s/a option		Left	Right		228	
			JMV110162	○	○	JMV110161		
JMV110153	□	Powered s/a option operated with PIKO button <i>(operated via one Piko button; only RH mounted)</i>						1693
PIKO BUTTON SWITCHES <i>per switch</i>								
JMV110191	□	Red	x 1	JMV110192	□	Red	x 2	87 each
JMV110193	□	Green	x 1	JMV110194	□	Green	x 2	87 each
JMV110195	□	Yellow	x 1	JMV110196	□	Yellow	x 2	87 each
JMV110197	□	Blue	x 1	JMV110198	□	Blue	x 2	87 each


MID WHEEL DRIVE POSITION		
JMV080310	● JIVE M2 always with 6" solid castor wheel on front and rear	
SPEED		
JMV180005	<input type="radio"/> 6 kph (4 - pole)	
JMV180007	<input type="radio"/> 10 kph (4 - pole)	
JMV180008	<input type="checkbox"/> 13 kph (4 - pole) <i>(High Speed Low torque - Not recommended for obstacle climbing)</i>	358
DRIVE WHEELS - Tyres and rims always black with coloured inserts		
JMV070051	<input type="radio"/> 14" Pneumatic "V-type" <i>(without coloured insert at the wheel)</i>	
JMV070050	<input type="radio"/> 14" Solid "V-type"	
JMV070053	<input type="radio"/> 14" Pneumatic "knobby"	
JMV070052	<input type="radio"/> 14" Solid "Knobby"	
SUSPENSION		
JMV010901	<input type="radio"/> Performance Spring Sports suspension, <i>(Greater traction & more responsive ride. Reduced kerb climbing to 7cm max)</i>	
BATTERIES		
JMV160001	<input type="radio"/> 60 Ah maintenance free <i>(pair)</i>	
JMV160005	<input type="checkbox"/> 80 Ah maintenance free <i>(pair)</i>	111
CHARGER		
JMV160100	<input type="radio"/> 220V/ 8A automatic charger (UK)	
JMV160103	<input type="checkbox"/> 220V/ 10A automatic charger (UK) <i>(required for 13 kph speed option)</i>	168
SEAT SHROUD INSERT COLOURS - frame and chassis body is always matt black		
JMV140001	<input type="radio"/> White <i>(Colour RAL Code: 9010)</i>	
JMV140005	<input type="radio"/> Apple Green <i>(Colour RAL Code: 110 70 77)</i>	
JMV140003	<input type="radio"/> Blue <i>(Colour RAL Code: 260 40 45)</i>	
JMV140004	<input type="radio"/> Black (20% Gloss) <i>(Colour RAL Code: 9005)</i>	
JMV140006	<input type="radio"/> Orange <i>(Colour RAL Code: 060 60 80)</i>	
JMV140002	<input type="radio"/> Red <i>(Colour RAL Code: 040 40 67)</i>	
ACCESSORIES		
JMV090067	<input type="radio"/> Rear view mirror	
JMV090061	<input type="checkbox"/> 24/12V DC/ DC converter	228
JMV090505	<input type="checkbox"/> Golf ball joystick knob	33
JMV090506	<input type="checkbox"/> T-bar B38 joystick knob	33
JMV090507	<input type="checkbox"/> Soft ball joystick knob	33
JMV090508	<input type="checkbox"/> Ball joystick knob	33
JMV090509	<input type="checkbox"/> Chin cup joystick knob	33
JMV090510	<input type="checkbox"/> Stick joystick knob	33
JMV090511	<input type="checkbox"/> Mushroom joystick knob	33
JMV090512	<input type="checkbox"/> Swing-away tray table	171
JMV090069	<input type="checkbox"/> Bluetooth USB mouse receiver	105
Quickie Jive M2 SEDEO ERGO		£11,550
Option Price		£
Total Price		£




General




Example: eff. seat width 460 mm




Example: Min. seat width 400 mm



General



Example: eff. seat width 460 mm



Example: Min. seat width 400 mm



JAY090006



JAY090007



JAY090008



JAY090009



Fixed laterals



Swing-away laterals



JAY090029 - JAY090031



JAY090032 - JAY090034



JAY090035



JAY090036



JAY090037 - JAY090039



JAY090040 - JAY090042



JAY090043 - JAY090045



JAY090046 - JAY090048



JAY090049 - JAY090051



JAY090052



JAY090053



JAY090054



JAY090055

 <p>JMV130007</p>	 <p>JMV130006</p>	 <p>JMV030051</p>	 <p>JMV030406 & JMV030407</p>
 <p>JMV020509</p>	 <p>JMV020510</p>	 <p>JMV020341 & JMV020342</p>	 <p>JMV020512</p>
 <p>JMV030034</p>	 <p>JMV050111</p>		
 <p>JMV040123</p>	 <p>JMV040217 & JMV040218</p>	 <p>JMV040203 & JMV040204</p>	 <p>JMV040205 & JMV040206</p>
 <p>JUP040207 - .210 & .219 - .220</p>	 <p>JUP040211 / .212 / .213 / .214</p>	 <p>JMV110013</p>	 <p>JMV110056</p>
 <p>JMV110175</p>	 <p>JMV110176</p>	 <p>JMV110101</p>	



JMV110201



JMV110130



JMV110191 - JMV110198



JMV110133



JMV110146



JMV110132



JMV110138




JMV110139



Knobby



V-Type



JMV090067



JMV090061



JMV090505



JMV090506



JMV090507



JMV090508



JMV090509



JMV090510



JMV090511



Dealer Assessment & handover Set up Notes
NOT FOR ORDERS USE



BIOMETRIC SMART SEATING - MEMORY FUNCTIONS

Memory Function	Default Setting	Customer Requirement
1	Seated	
2	Lie Down	
3	Pressure Relief	
4	Relax	
5	Blank	
6	Blank	

--

PRESSURE RELIEF COACHING FUNCTION

Interval reminder	
--------------------------	--

Duration for pressure relief	
-------------------------------------	--

--

BIOMETRIC REPOSITIONING

Factory Setting: ON	Customer Requirement	<input type="checkbox"/> YES	<input type="checkbox"/> NO
----------------------------	-----------------------------	-------------------------------------	------------------------------------

ASSIGNABLE BUTTONS - EXMAPLE OF FUNCTIONS TO ASSIGN

Seating	Memory Functions	Latched Functions (Moves automatically)	Connectivity
Tilt Forward	Lie Flat	Tilt - Latched	InfraRED - TV ON
Tilt Backward	Transfer Position	Lie Down - Latched	InfraRED - TV OFF
Tilt - Toggle (rotates between forward/back)	Custom Memory Position	Custom Memory Position - Latched	Bluetooth

DEALER NOTES FOR ASSESSMENT TO PROGRAM AT HANDOVER

TOP LEFT BUTTON Factory Setting (1st Function)	
HAZARD ON	
Re-assign 1st Function (short press)	Re-assign 2nd Function (Long press)

TOP RIGHT BUTTON Factory Setting (1st Function)	
LIGHTS ON	
Re-assign 1st Function (short press)	Re-assign 2nd Function (Long press)

BOTTOM LEFT BUTTON Factory Setting (1st Function)	
LEFT INDICATOR	
Re-assign 1st Function (short press)	Re-assign 2nd Function (Long press)

BOTTOM RIGHT BUTTON Factory Setting (1st Function)	
RIGHT INDICATOR	
Re-assign 1st Function (short press)	Re-assign 2nd Function (Long press)

MODE BUTTON Factory Setting (1st Function)	
MODE	
Re-assign 1st Function (short press)	Re-assign 2nd Function (Long press)
DO NOT CHANGE	

TOGGLE UP Factory Setting (1st Function)	
SPEED UP	
Re-assign 1st Function (short press)	Re-assign 2nd Function (Long press)

PROFILE BUTTON Factory Setting (1st Function)	
PROFILE	
Re-assign 1st Function (short press)	Re-assign 2nd Function (Long press)
DO NOT CHANGE	

TOGGLE DOWN Factory Setting (1st Function)	
SPEED DOWN	
Re-assign 1st Function (short press)	Re-assign 2nd Function (Long press)



BUTTON 1	
PUSH FORWARD	PUSH BACK
LIE DOWN	SEATED
REASSIGN FORWARD	REASSIGN BACK



BUTTON 2	
PUSH FORWARD	PUSH BACK
LIFT UP	LIFT DOWN
REASSIGN FORWARD	REASSIGN BACK

BUTTON 3	
PUSH FORWARD	PUSH BACK
TILT BACK	TILT FORWARD
REASSIGN FORWARD	REASSIGN BACK



BUTTON 4	
PUSH FORWARD	PUSH BACK
RECLINE BACK	RECLINE UP
REASSIGN FORWARD	REASSIGN BACK

BUTTON 5	
PUSH FORWARD	PUSH BACK
LEGREST EXTEND UP	LEGREST BACK
REASSIGN FORWARD	REASSIGN BACK



Dealer Assessment & handover Set up Notes
NOT FOR ORDERS USE



JIVE M2 SEDEO ERGO CONNECTIVITY

MODULE	DEVICE TO BE CONNECTED & SET UP
BLUETOOTH - ANDROID 1	
BLUETOOTH - ANDROID 2	
BLUETOOTH - APPLE	
INFRA RED	TV:
	TV2:
	TV3:
	DVD PLAYER:
	SKY BOX:
	STEREO:
	ENVIRONMENTAL PRODUCTS - GEWAH CODES (DOOR OPENING, BLIND OPENERS, INFRA RED CENTRAL HEATING CONTROL, DOOR INTERCOMMS) (Client must have pre existing Infra red devices installed in the home - Rnet advanced can then control these devices)