

Triride®

by Gianni Conte

UNA NUOVA LIBERTÀ'
A NEW FREEDOM



Max Power
720W

Potenza

Power

★★★★☆

Leggerezza

Lightness

★★★★★

Compattezza

Compactness

★★★★☆

Comfort Guida

Driving comfort

★★★★☆



mod. **BASE**



Dispositivo medico iscritto all'Albo con il numero: 922795 / TUV Certification No. 15-PEL-0010022-TIC

CARATTERISTICHE TECNICHE



Potenza Massima 720W

Elettronica by Triride

Retromarcia

Cruise Control

Sollevamento Assistito

Freno Elettronico

Peso: 8,4 kg

Cerchio in alluminio da 12" o 14"

Display LCD

Batteria Litio-ioni cilindrica da 36V

Freno a disco

Pneumatico: City o Cross

Telaio in acciaio inox

Velocità a norma di legge regolabile con
acceleratore

5 Livelli di velocità

TECHNICAL FEATURES



Max Power 720W

Electronics by Triride

Reverse

Cruise Control

Assisted Load

Electronic Brake

Weight: 8,4 Kg

12" or 14" driving wheel with strengthened
aluminium rim

LCD Display

36V Lithium-ion battery

Powerful disc brake

Tire selection: City or Cross

Stainless steel frame

Speed compliant with legal standards

5 Speed levels