

Triride®

by Gianni Conte

UNA NUOVA LIBERTA'
A NEW FREEDOM



Max
Power
540W

Potenza

Power

★★★★☆

Leggerezza

Lightness

★★★★★

Compattezza

Compactness

★★★★★

Comfort Guida

Driving comfort

★★★★☆

mod. KIDS

TUV Certification No. 15-PEL-0010022-TIC



CARATTERISTICHE TECNICHE



Potenza Massima 540W limitata e controllata
con Elettronica by Triride

Retromarcia

Cruise Control

Sollevamento Assistito

Freno Elettronico

Peso: 7,5 kg

Cerchio in alluminio da 12"

Display LCD (opzionale)

Batteria Litio-ioni cilindrica da 36V

Freno a disco

Pneumatico: City o Cross

Telaio in acciaio inox

Velocità a norma di legge regolabile con
acceleratore

5 Livelli di velocità (con display opzionale)

TECHNICAL FEATURES



Max Power 540W limited and controlled by
Triride Electronics

Reverse

Cruise Control

Assisted Load

Electronic Brake

Weight: 7,5 Kg

12" driving wheel with strengthened
aluminium rim

LCD Display (optional)

36V Lithium-ion battery

Powerful disc brake

Tire selection: City or Cross

Stainless steel frame

Speed compliant with legal standards

5 Speed levels (with optional display)