

# Triride®

by Gianni Conte

UNA NUOVA LIBERTA'  
A NEW FREEDOM



Max Power  
1500W

**Potenza**  
Power  
★★★★★+

**Leggerezza**  
Lightness  
★★★★☆

**Compattezza**  
Compactness  
★★★★☆

**Comfort Guida**  
Driving comfort  
★★★★★



## mod. SPECIAL L14



Dispositivo medico iscritto all'Albo con il numero: 1560510 / TUV Certification No. 15-PEL-0010022-TIC

### CARATTERISTICHE TECNICHE



Potenza Massima 36V 1200W / 48V 1500W

Electronica by Triride

Retromarcia

Cruise Control

Sollevamento Assistito

Freno Elettronico

Peso: 12 kg

Cerchio in alluminio da 14"

Display LCD

Batteria Litio-ioni da 36V o 48V

Freno a disco

Pneumatico: City o Cross

Telaio in acciaio inox

Velocità a norma di legge regolabile con  
acceleratore

5 Livelli di velocità

### TECHNICAL FEATURES



Max Power 36V 1200W / 48V 1500W

Electronics by Triride

Reverse

Cruise Control

Assisted Load

Electronic Brake

Weight: 12 Kg

14" driving wheel with strengthened  
aluminium rim

LCD Display

36V or 48V Lithium-ion battery

Powerful disc brake

Tire selection: City or Cross

Stainless steel frame

Speed compliant with legal standards

5 Speed levels