

Triride®

by Gianni Conte

UNA NUOVA LIBERTA'
A NEW FREEDOM



Max Power
1000W

Potenza

Power

★★★★☆

Leggerezza

Lightness

★★★★☆

Compattezza

Compactness

★★★★☆

Comfort Guida

Driving comfort

★★★★☆



mod. SPECIAL LIGHT



Dispositivo medico iscritto all'Albo con il numero: 1560509 / TUV Certification No. 15-PEL-0010022-TIC

CARATTERISTICHE TECNICHE



Potenza Massima 36V 750W / 48V 1000W

Elettronica by Triride

Retromarcia

Cruise Control

Sollevamento Assistito

Freno Elettronico

Peso: 8,9 kg

Cerchio in alluminio da 12" o 14"

Display LCD

Batteria Litio-ioni da 36V o 48V

Freno a disco

Pneumatico: City o Cross

Telaio in acciaio inox

Velocità a norma di legge regolabile con
acceleratore

5 Livelli di velocità

TECHNICAL FEATURES



Max Power 36V 750W / 48V 1000W

Electronics by Triride

Reverse

Cruise Control

Assisted Load

Electronic Brake

Weight: 8,9 Kg

12" or 14" driving wheel with strengthened
aluminium rim

LCD Display

36V or 48V Lithium-ion battery

Powerful disc brake

Tire selection: City or Cross

Stainless steel frame

Speed compliant with legal standards

5 Speed levels