

# Triride®

by Gianni Conte

UNA NUOVA LIBERTA'  
A NEW FREEDOM



Max  
Power  
250W

Leggerezza  
Lightness



Compattezza  
Compactness



mod. **TRIBIKE**  Hybrid



Dispositivo medico iscritto all'Albo con il numero: 1776590 / TUV Certification No. 15-PEL-0010022-TIC

## CARATTERISTICHE TECNICHE



Potenza Motore 250W

Elettronica by Triride

Sensore pedalata assistita PAS

Posizione pedalata completamente regolabile

Assistenza elettrica partenza/salita

Velocità con supporto motore max 25 Km/h

Peso: 12 Kg (senza batteria)

Cerchio in alluminio da 20"

Display LCD con 5 livelli di assistenza

Assistenza programmabile e personalizzabile

Batteria Litio-ioni cilindrica da 36V

Doppio freno (disk brake e V-brake)

Cambio Shimano a 7 rapporti

Telaio in acciaio con colore a scelta

Pneumatico: City

Verniciatura a caldo ad alta resistenza

## TECHNICAL FEATURES



Max Power 250W

Electronics by Triride

PAS assisted riding sensor

Fully adjustable pedal position

Electrical assistance start/climb

Speed with motor support max 25Km/h

Weight: 12 Kg (without battery)

Aluminium 20" driving wheel

LCD Display with 5 levels of assistance

Programmable and customizable assistance

36V Lithium-ion battery

Double brake (disk brake e V-brake)

7-speed Shimano gearbox

Steel frame with customizable color

Tire selection: City

High resistance painting