


zippie®
Explore your world!



Youngster 3



Fit for Adventures



Youngster³

The lightest folding
wheelchair with growth
adjustability!

Zippie Youngster 3:

The new lightweight Youngster 3 offers integrated seat width and seat depth growth. The new frame design gives your wheelchair a stylish and unique look.

On top of that the Youngster 3 offers new, innovative and lightweight options using the latest technology.





Patented 2cm seat width growth via cross brace for growth (without new parts).



New lightweight fork and frame for maximum adjustability, improved looks and better manoeuvrability.

- Lightweight package:

Lightweight fork, lightweight seat upholstery, lightweight anti-tips, Safari wheel lock



4cm seat depth growth via seat tubes for growth (without new parts).

Explore your world! Have fun, enjoy playing with your friends and be active

Options

The new Youngster 3 offers a wide range of options for more comfort, freedom, fun and safety.



New selection of back upholsteries with a wide range of different versions and improved adjustability.



Lightweight seat upholstery for easy handling and transport.

Thirty two frame colours!

Select a colour for each component.




Modern side guard design in combination with the new Safari wheel lock for improved safety.

Spoke guards

There are 32 fantastic spoke guard designs you can choose from. See our spoke guard catalogue for the whole range.



Technical specifications



Frame colour :	32 colours available	Weight of heaviest part:	4.4kg
Upholstery colour:	Black / Vented (black & silver) / Alcantara (grey & brown)	Product weight:	from 8.5kg
Max user weight:	85kg (13.5 stone)	Front wheel diameter:	10, 12, 15cm (4", 5", 6")
Growth: Seat width:	+2cm (0.7") without parts	Rear wheel diameter:	51, 55, 61cm (20", 22", 24")
Seat depth:	+4cm (1.6") without parts		

Frame Colours

1 White RAL 9010	2 Brilliant Silver	3 Silver Matt	4 Silver	6 Sun Yellow	7 Lime	8 Matt Lime	9 Candy Teal
10 Matt Candy Teal	14 Blue	15 Brilliant Marine	16 Matt Marine	19 Brilliant Purple	20 Matt Purple	21 Red Ruby	22 Brilliant Red
23 Matt Red	24 Tutti Frutti Orange	27 Graphite Metallic	28 Platinum Metallic	29 Black	31 Matt Black	32 Signal white	33 Pink
34 Dark Chocolate Metallic	35 Sky Blue	36 Cyan Metallic	37 Mars Red	39 Cappuccino	40 Kiwi	41 Flamingo Glitter	42 Matt Anthracite

It is best to select colours using the metal colour plates at your dealers and to use the order form to finalise your desired colour scheme



PB_Youngster3_EN_GB_REV3

Find out more today!

Visit our website at www.SunriseMedical.co.uk for more information or to request a demonstration.

Sunrise Medical Limited

Thorns Road,
Brierley Hill,
West Midlands
DY5 2LD
Tel.: +44 (0) 845 605 6688
Fax: +44 (0) 845 605 6689
www.SunriseMedical.co.uk