

JAY Zip Back

Clinical seating, designed for kids



The JAY Zip backrest was designed through the use of anthropometric data so it is specially sized to meet a paediatric client's unique needs. It is an ultra lightweight, anti-microbial backrest that features X-static® and Aquaguard™ cover technologies as well as adjustable, one-step release hardware to accommodate a wide range of mounting locations and user requirements.

Versatile Hardware

- Easy one-step release
- Fits on a wide range of tubing: 1.9cm (0.7") with inserts; 2.2 - 2.3cm (0.8 - 0.9") with inserts or 2.5cm (1") without insert
- One tool adjustments (and the tool is provided!)

Convenient & Anti-microbial

Outer Cover

- Machine washable and quick to dry
- Silver impregnated X-static® stretch fabric naturally inhibits microbial growth
- Extra layer of microclimatic spacer fabric for comfort
- Safe non-toxic material

Inner Cover

- Easy to wipe off
- Aquaguard™ zipper resists moisture
- Anti-wicking seam thread protects foam base



One step release

Backrest width:	20, 25, 30, 35cm (8, 10, 12, 14")	Back angle:	+20° to -20°
Backrest height:	15, 20, 25, 30cm (6, 8, 10, 12")	Width adjustability:	+5cm (2")
Transport weight:	1.1kg	Depth adjustability:	7.5 cm (3")
Max. user weight:	75kg	Accessories: Headrest, lateral support, spine align positioning components, chest strap, vanity flap and spare cover	



Sunrise Medical Limited

Thorns Road
 Brierley Hill
 West Midlands
 DY5 2LD
 Tel.: +44 (0) 845 605 66 88
 Fax: +44 (0) 845 605 66 89
 www.SunriseMedical.co.uk