



Children's wheelchair

The smallest wheelchair for kids.

The Permobil **Koala Miniflex** is the ideal wheelchair for children, from toddlers to young teenagers. It is equipped with the same powerful motors and front-wheel drive as Permobil's wheelchairs for adults, but adapted to the needs of children. The Miniflex seating system is a flexible seating system for growing children and children with special seating needs. Available in three sizes and with a range of body supports. Options include installation of an orthotic seat.



Koala
MINIFLEX

power functions





Bumper

For safe and comfortable driving available as an option.



Seat

Equipped as standard with combined electrical tilt and height adjustment.



Rear Wheels

Swivel rear wheels for responsive steering.



Footplate

Two sizes.

Colours



apple green




lemonade yellow

Upholstery



leatherette



black fabric

Specifications

| | |
|------------------------------------|--------------------------|
| Total length | 965 mm |
| Total width | 585 mm |
| Minimum transport length | 710 mm |
| Seat height (fixed) | 340 mm |
| Seat height with powered seat lift | 410 - 660 mm |
| Total weight incl. batteries | 94 kg |
| Maximum user weight | 35 kg |
| Seat width | 250 - 300 - 350 mm |
| Seat depth | 270 - 320 - 370 mm |
| Backrest width | 270 - 320 - 370 mm |
| Backrest height | 325 - 375 - 425 - 475 mm |
| Distance between armrests | 250 - 350 mm |
| Backrest angle - adjustable | 80° - 145° |
| Leg support angle - adjustable | 90° - 150° |

| | |
|-----------------------------|------------|
| Battery capacity | 2 x 40 Ah |
| Range (ISO 7176-4) | 25 - 30 km |
| Maximum speed | 6.5 km/hr |
| Turning radius (ISO 7176-5) | 570 mm |
| Obstacle height | 60 mm |
| Suspension | no |
| Electronics | R-net |
| Powered seat lift | yes |
| Powered seat tilt | yes |
| Powered leg rest adjustment | no |
| Powered backrest adjustment | no |
| Guarantee | 2 years |
| EN12182/EN12184 CE | yes |